

4. What bioanalytical methods are used to assess concentrations?

See the response to section IV.F.1 above.

**APPEARS THIS WAY
ON ORIGINAL**

**APPEARS THIS WAY
ON ORIGINAL**

V. LABELING RECOMMENDATIONS

Please refer to Appendix A, Annotated Label.

APPEARS THIS WAY
ON ORIGINAL

39 pages redacted from this section of
the approval package consisted of draft labeling

Appendix B. Individual Study Reviews

An investigation of the potential for daptomycin to inhibit cytochrome P450 1A2, 2A6, 2C9, 2C19, 2D6, 2E1, and 3A4 in cryopreserved human hepatocytes (ADME Report #12)

OBJECTIVE:

The purpose of the study was to assess the potential of daptomycin to inhibit hepatic cytochrome P450 mediated metabolism via CYP1A2, CYP2A6, CYP2C9, CYP2C19, CYP2D6, CYP2E1, CYP 3A4.

FORMULATION:

Daptomycin (Lot No. X800384)

STUDY DESIGN:

The *in vitro* evaluation of daptomycin as an inhibitor of human cytochrome P450 (CYP P450) isozymes was performed using daptomycin concentrations of 2.5, 10, and 40 µg/mL. Human hepatocytes were pooled from five adult male (all Caucasian) and five adult female (4 Caucasian, 1 Hispanic) subjects. Hepatocytes were incubated in protein-free Krebs-Henseleit buffer (KHB) containing 0.1 M CaCl₂.

Isolated hepatocytes were diluted with suspension media (Dulbecco's modified Eagle medium) to determine viability using Trypan Blue. Cells were then centrifuged and diluted with incubation medium (KHB) to prepare a cell suspension (2×10^6 cells/mL).

The CYP form-specific activities were evaluated using the following probe substrates: 50 µM phenacetin (CYP1A2), 50 µM coumarin (CYP2A6), 75 µM tolbutamide (CYP2C9), 50 µM S-mephenytoin (CYP2C19), 8 µM dextromethorphan (CYP2D6), 50 µM chlorzoxazone (CYP2E1), and 50 µM testosterone (CYP3A4). A negative control consisted of incubation media alone (KHB buffer). Positive control inhibitors consisted of 10 µM furafylline (CYP1A2), 50 µM diethyldithiocarbamate (CYP2A6), 1 µM sulfaphenazole (CYP2C9), 10 µM omeprazole (CYP2C19), 1 µM quinidine (CYP2D6), 100 µM 4-methylpyrazole (CYP2E1), and 1 µM ketoconazole (CYP3A4).

The rate of enzyme activity was assessed by the rate of formation of acetaminophen (CYP1A2), 7-hydroxycoumarin (CYP2A6), 4-hydroxytolbutamide (CYP2C9), 4-hydroxymephenytoin (CYP2C19), dextrophan (CYP2D6), 6-hydroxychlorzoxazone (CYP2E1), and 6β-hydroxytestosterone (CYP3A4). The percentage of the activity remaining was calculated as the ratio of enzyme activity in the presence of daptomycin relative to the enzyme activity in the presence of negative control. The impact of the positive control was calculated using the vehicle control (incubation medium + 0.1% DMSO) instead of the negative control.

RESULTS:

The viability of hepatocytes in suspension media using Trypan blue was 73%.

The positive controls yielded greater than 50% inhibition of enzyme activity and adequately inhibited the activity of the isoforms.

The rate of enzyme activity (pmol/million cells/min) in the presence of negative control, positive control, and daptomycin is shown in Table 1. The percent of the activity remaining for each P450 isozyme is shown in Table 2.

Daptomycin, at concentrations of 2.5, 10, and 40 µg/mL, did not appreciably inhibit of the activity of CYP1A2, CYP2A6, CYP2C9, CYP2D6, and CYP3A4 in human hepatocytes. The activity of CYP2E1

was reduced in the presence of daptomycin (67.5% to 73.2%), although the inhibition was less than 50% of the positive control (4-methylpyrazole).

Table 1. Effects of daptomycin on *in vitro* metabolism rates of major human CYP P450 specific activities

CYP P450	Probe Substrate	Enzyme activity (pmol/million cells/min)				
		(-) Control	(+) Control	2.5 µg/mL	10 µg/mL	40 µg/mL
1A2	Phenacetin	13.5 ± 0.3	1.02 ± 0.06	14.1 ± 0.2	14.4 ± 0.3	14.1 ± 0.2
2A6	Coumarin	3.64 ± 0.17	0.49 ± 0.04	2.94 ± 0.13	3.43 ± 0.30	3.12 ± 0.10
2C9	Tolbutamide	10.2 ± 0.3	1.98 ± 0.04	10.4 ± 0.3	10.4 ± 0.2	10.9 ± 0.3
2C19	S-Mephenytoin	8.69 ± 0.43	2.02 ± 0.15	8.05 ± 0.41	8.41 ± 0.15	8.53 ± 0.31
2D6	Dextromethorphan	10.0 ± 0.3	1.15 ± 0.16	11.7 ± 0.4	12.4 ± 0.2	12.4 ± 0.3
2E1	Chlorzoxazone	6.61 ± 0.25	0.49 ± 0.04	4.84 ± 0.15	4.60 ± 0.23	4.46 ± 0.52
3A4	Testosterone	20.1 ± 3.2	0.77 ± 0.15	20.7 ± 3.5	20.6 ± 2.6	32.6 ± 13.6

Table 2. Percentage of activity remaining of major human CYP P450 isoenzymes in the presence of daptomycin

CYP P450	Probe Substrate	Activity remaining (%)			
		Positive Control	2.5 µg/mL	10 µg/mL	40 µg/mL
1A2	Phenacetin	8.2%	104.4%	106.7%	104.4%
2A6	Coumarin	13.7%	80.8%	94.2%	85.7%
2C9	Tolbutamide	21.1%	102.0%	102.0%	106.9%
2C19	S-Mephenytoin	25.3%	92.6%	96.8%	98.2%
2D6	Dextromethorphan	10.2%	117.0%	124.0%	124.0%
2E1	Chlorzoxazone	17.2%	73.2%	69.6%	67.5%
3A4	Testosterone	5.1%	103.0%	102.5%	162.2%

Following the administration of daptomycin IV 4 mg/kg q24h, anticipated peak plasma concentrations of daptomycin are approximately 50 µg/ml (total) and 5 µg/mL (unbound). Thus, the maximum concentration of daptomycin assessed in the study exceeded the anticipated peak plasma concentrations of daptomycin by approximately 8-fold.

CONCLUSIONS:

Based on the *in vitro* results, daptomycin IV 4 mg/kg is unlikely to inhibit the metabolism of drugs dependent on P450 isoforms CYP1A2, CYP2A6, CYP2C9, CYP2C19, CYP2D6, CYP2E1, or CYP3A4.

COMMENTS:

Although the sponsor states that the activities of all substrates were analyzed using proprietary methods, assay validation data were not provided with the study report. The sponsor is encouraged to provide the full assay validation data with the study report in the future.

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.