

Dose escalation and pharmacokinetic study of a humanized anti-HER2 monoclonal antibody in patients with HER2/*neu*-overexpressing metastatic breast cancer

Y Tokuda¹, T Watanabe², Y Omuro², M Ando², N Katsumata², A Okumura¹, M Ohta¹, H Fujii³, Y Sasaki³, T Niwa⁴ and T Tajima¹

¹Department of Surgery, Tokai University School of Medicine, Bohseidai, Isehara, Kanagawa 259-1193, Japan; ²National Cancer Center Hospital, Tokyo 104-0045, Japan; ³National Cancer Center Hospital East, Kashiwa, Chiba 277-0882, Japan; ⁴Mitsubishi Chemical Corporation, Yokohama Research Center, Kanagawa 227-8502, Japan

Summary We conducted a phase I pharmacokinetic dose escalation study of a recombinant humanized anti-p185^{HER2} monoclonal antibody (MKC-454) in 18 patients with metastatic breast cancer refractory to chemotherapy. Three or six patients at each dose level received 1, 2, 4 and 8 mg kg⁻¹ of MKC-454 as 90-min intravenous infusions. The first dose was followed in 3 weeks by nine weekly doses. Target trough serum concentration has been set at 10 µg ml⁻¹ based on in vitro observations. The mean value of minimum trough serum concentrations at each dose level were 3.58 ± 0.63, 6.53 ± 5.26, 40.2 ± 7.12 and 87.9 ± 23.5 µg ml⁻¹ respectively. At 2 mg kg⁻¹, although minimum trough serum concentrations were lower than the target trough concentration with a wide range of variation, trough concentrations increased and exceeded the target concentration, as administrations were repeated weekly. Finally 2 mg kg⁻¹ was considered to be sufficient to achieve the target trough concentration by the weekly dosing regimen. One patient receiving 1 mg kg⁻¹ had grade 3 fever, one at the 1 mg kg⁻¹ level had severe fatigue defined as grade 3, and one at 8 mg kg⁻¹ had severe bone pain of grade 3. No antibodies against MKC-454 were detected in any patients. Objective tumour responses were observed in two patients; one receiving 4 mg kg⁻¹ had a partial response in lung metastases and the other receiving 8 mg kg⁻¹ had a complete response in soft tissue metastases. These results indicate that MKC-454 is well tolerated and effective in patients with refractory metastatic breast cancers overexpressing the HER2 proto-oncogene. Further evaluation of this agent with 2–4 mg kg⁻¹ weekly intravenous infusion is warranted. © 1999 Cancer Research Campaign

Keywords: HER2/*neu*; humanized monoclonal antibody; pharmacokinetics; phase I study

The *c-erbB-2/HER2* proto-oncogene encodes a receptor-type tyrosine kinase (Yarden and Ullrich, 1988) corresponding to, but distinct from, the epidermal growth factor receptor (Coussens et al, 1985; Yamamoto et al, 1986). The HER2 product consists of three domains, extracellular, transmembrane and intracellular domains which has tyrosine kinase activity making the HER2 product autophosphorylated. Appreciable amplification and/or overexpression of this gene has been demonstrated in a variety of adenocarcinomas including breast cancer (King et al, 1985; Slamon et al, 1989). However, the expression of this gene in normal adult tissues is weak (De Potter et al, 1989; Press et al, 1990). Therefore, the HER2 product is thought to be a useful target for antibody therapy of cancers overexpressing the HER2 gene.

Several studies of murine monoclonal antibodies (Mabs) directed against the HER2 product have already been reported to have in vitro and in vivo anti-tumour effects (Drebin et al, 1985; Hudziak et al, 1989; Hancock et al, 1991; Stancovski et al, 1991; Tagliabue et al, 1991; Harwerth et al, 1992; Kasprzyk et al, 1992). However, human anti-mouse antibody response during therapy

would be a major limitation in the clinical application of such murine mabs (Schroff et al, 1985; Shawler et al, 1985). Carter and his colleagues constructed a humanized antibody containing only the antigen-binding loops from a murine mab against the extracellular domain of the HER2 gene product and human variable region framework residues plus IgG1 constant domains (Carter et al, 1992). A clinical study of the humanized antibody revealed some objective responses without antibodies against it (Baselga et al, 1996).

We report here the results of a phase I dose escalation study of the humanized mab (MKC-454) against the HER2 product in patients with metastatic breast cancer refractory to conventional chemotherapeutic agents and positive for this oncogene product.

PATIENTS AND METHODS

Patient eligibility

Patients with metastatic breast cancer refractory to conventional chemo or endocrine therapies and overexpressing the HER2 product were eligible for this study. Patients were considered refractory to endocrine therapy if hormone receptors were negative or tumour growth was noted during treatment with tamoxifen, and were considered refractory to chemotherapy if disease progression was noted after treatment with doxorubicin-containing regimens. All patients had measurable disease,

Received 20 November 1998

Revised 14 May 1999

Accepted 7 June 1999

Correspondence to: Y Tokuda

Table 1 Patient demographics

Parameter	No.
No. of patients	18
Age median (years) (range)	51 (32–64)
Performance score (ECOG)	
0	10
1	7
2	1
Receptor status	
Oestrogen receptor	
Positive	3
Negative	11
Unknown	4
Progesterone receptor	
Positive	3
Negative	10
Unknown	5
Histological grade	
II	3
III	12
Unknown	3
Prior treatment	
Chemotherapy	
Adjuvant chemotherapy	12
Neoadjuvant chemotherapy	1
Metastatic disease (no. of regimens)	
1	7
2	7
>2	4
Median (range)	2 (1–4)
Hormonal therapy	
Adjuvant therapy	11
Metastatic disease	6

satisfactory haematologic, hepatic, renal and pulmonary function with an Eastern Cooperative Oncology Group (ECOG) performance status grade 3 or less. Patients who had chemotherapy within 4 weeks (6 weeks for mitomycin or nitrosurea) or hormonal therapy within 2 weeks before study entry were ineligible. Written informed consent was mandatory before study entry.

The expression of HER2 was determined by immunohistochemical staining of the paraffin-embedded thin sections of the primary or metastatic tumours using rabbit polyclonal antibodies (DAKO, Glostrup, Denmark and NICHIREI, Tokyo, Japan) or a murine monoclonal antibody (Lab Corp, North Carolina, USA) against the human HER2 product. Tumours were considered to overexpress HER2 if at least 10% of tumour cells had positive membrane staining.

Treatment

A recombinant humanized anti-HER2 mAb (trastuzumab, MKC-454) was constructed by molecular engineering from a murine mAb recognizing the extracellular domain of the HER2 product (Carter *et al.*, 1992). The drug was administered intravenously at a constant infusion rate for 90 min. The starting dose level of MKC-454 was 1 mg kg⁻¹, and subsequent dose escalations were to 2, 4 and 8 mg kg⁻¹. The first dose was followed in 3 weeks by nine weekly doses. A minimum of three patients assessable for toxicity were treated at each dose level. If one of the first three patients entered at any dose level experienced a dose-limiting toxicity, an additional three patients were entered at that dose level. Toxicities

were graded according to the JCOG toxicity criteria (Tobinai *et al.*, 1993). Dose-limiting toxicity was defined as any grade 3 or 4 toxicity.

Pharmacokinetics, determination of antibodies against the humanized antibody and circulating shed antigen

Serum concentrations of MKC-454 were determined in a receptor binding assay that detects binding with the extracellular domain of the HER2 product. The recombinant extracellular domain of the HER2 product provided from Genentech, Inc. was coated on 96-well microtitre plates. MKC-454 bound to the coated antigen was detected by horseradish peroxidase-labelled goat anti-human IgG Fc. The quantitative limit for MKC-454 in serum samples was 156 ng ml⁻¹. The presence of antibodies against MKC-454 was determined with a bridging-type titre enzyme-linked immunosorbent assay (ELISA). Serum concentrations of circulating extracellular domain of the HER2 product, circulating shed antigen were measured using an ELISA. A pair of polyclonal antibodies against extracellular domain of the HER2 product were provided from Genentech, Inc. The quantitative limit for circulating shed antigen was 1.54 ng ml⁻¹. Pharmacokinetic parameters of MKC-454 after first administration were determined for each patient. Maximum serum concentration (C_{max}) was the observed value.

Half-life ($t_{1/2}$), area under concentration-time curve ($AUC_{0-\infty}$), distribution volume (V) and serum clearance (CL) were estimated by non-compartmental model of WinNonlin Standard Japanese Edition Version 1.1 (Scientific Consulting Inc.). As for the trough serum levels during the weekly treatment period, minimum and maximum values were obtained for each patient.

Evaluation of response

All responses were reviewed centrally by an independent extramural evaluation committee. A complete response (CR) was defined as the complete resolution of all measurable tumours for a duration of at least 4 weeks. Partial response (PR) was defined as a $\geq 50\%$ reduction in the sum of the products of the two longest perpendicular diameters of all measurable tumours for a duration of at least 4 weeks, a minimal response (MR) as a $\geq 25\%$ but less than 50% reduction in the sum of the products of the two longest perpendicular diameters of all tumours, and progressive disease (PD) as a $\geq 25\%$ increase in any measurable lesions or the appearance of any new lesion.

RESULTS

Patient characteristics are listed in Table 1. All patients had received chemotherapy for metastatic disease and 11 patients had received two or more regimens. All of the metastatic diseases had become refractory to anthracycline-containing regimens. Although hormonal therapy for metastatic disease was also given to a third of the patients, 11 primary tumours were negative for oestrogen receptor. Twelve patients had pathological grade III primary tumours.

Serum concentration-time profiles after first administration are illustrated in Figure 1 and estimated pharmacokinetic parameters for each patient are summarized in Table 2. Mean C_{max} increased in good proportion to dose. Mean distribution volume (V) ranging from 52.1 to 74.3 ml kg, was independent of dose, and approximated the serum space. Half-life ($t_{1/2}$) increased from 2.7 days at

Table 2 Single dose pharmacokinetic parameters of MKC-454 in each patient

Dose (mg kg ⁻¹)	Patient no.	C _{max} (ng ml ⁻¹)	t _{1/2} (day)	AUC _{0-∞} (μg day ml ⁻¹)	V (ml kg ⁻¹)	CL (ml day kg ⁻¹)	
1	1-102	17 506	2.6	67	56.0	15.0	
	1-103 ^a	14 994	2.3	47	45.2	14.2	
	1-104	21 664	2.4	70	49.6	14.2	
	1-105	22 230	3.2	94	48.0	10.6	
	1-106	18 590	2.8	85	45.9	11.8	
	3-107	15 770	2.3	53	61.1	18.9	
	Mean ^b	19 152	2.7	74	52.1	14.1	
	s.d. ^b	2 750	0.4	16	6.3	3.2	
	2	3-108	52 370	1.6	158	33.4	12.7
		2-109	42 500	2.3	198	49.1	10.1
1-110		35 270	5.3	190	75.4	10.5	
Mean		43 413	3.1	182	52.6	11.1	
s.d.		8 537	2.0	21	21.2	1.4	
4		3-111	62 340	10.2	549	85.9	7.3
	3-112	62 720	7.7	611	77.0	6.5	
	1-113	92 320	8.4	746	60.0	5.4	
	Mean	72 460	8.8	635	74.3	6.4	
	s.d.	17 200	1.3	101	13.2	1.0	
	8	3-114	167 600	6.4	928	68.8	8.6
1-115		193 500	15.4	2101	77.1	3.8	
1-116		197 050	10.3	1954	56.9	4.1	
2-117		149 850	10.2	1551	72.4	5.2	
3-118		176 700	11.3	1810	67.5	4.4	
3-119		133 250	9.0	1111	79.5	7.2	
Mean		169 658	10.4	1576	70.4	5.6	
s.d.		24 862	3.0	471	8.1	1.9	

^a0.67 mg/kg of dose was used for calculation of pharmacokinetic parameters because infusion was discontinued at 60 min.

^bPatient no. 1-103 was deleted from calculation of mean and s.d. due to different dose.




Figure 1 Serum concentration-time profiles of MKC-454 after first administration. The values of serum concentration at a dose level of 1 mg kg⁻¹ (open circle), 2 mg kg⁻¹ (closed circle), 4 mg kg⁻¹ (open triangle), or 8 mg kg⁻¹ (closed triangle) represent mean ± s.d.

the lowest dose (1 mg kg⁻¹) to 10.4 days at the highest dose (8 mg kg⁻¹). Mean serum clearance (CL) decreased dose-dependently from 14.1 to 5.6 ml day⁻¹ kg⁻¹. Trough serum concentrations during repeated administration of MKC-454 are individually summarized in Figure 2. Maximum values of trough serum concentrations at each dose level were 9.19 ± 2.26, 19.1 ± 15.1, 102 ± 47.9, and 248 ± 64.4 μg ml⁻¹ respectively. The mean value




Figure 2 Trough serum concentrations of MKC-454. Left and right ends of each horizontal bar represent minimum and maximum trough serum concentration of MKC-454 during weekly administration

of minimum trough serum concentrations at each dose level were 3.58 ± 0.63, 6.53 ± 5.26, 40.2 ± 7.12 and 87.9 ± 23.5 μg ml⁻¹ respectively. Figure 3 shows the serum concentrations of circulating extracellular domain of the HER2 product, shed antigen. One patient (Patient No. 3-108) at 2 mg kg⁻¹ with high values of circulating shed antigen showed reduced values of the trough serum concentrations compared with those of the other two

Table 3 Toxicity of MKC-454

	Dose 1 mg kg ⁻¹ n = 6 Grade			Dose 2 mg kg ⁻¹ n = 3 Grade			Dose 4 mg kg ⁻¹ n = 3 Grade			Dose 8 mg kg ⁻¹ n = 6 Grade		
	I	II	III	I	II	III	I	II	III	I	II	III
Fever	–	2	1	–	1	–	1	1	–	–	1	–
Nausea/vomiting	–	1	–	–	–	–	1	–	–	–	1	–
Gastrointestinal disorder	–	–	1	–	–	–	–	–	–	–	–	–
Liver												
ASAT	4	–	–	–	–	–	–	1	–	1	–	–
ALAT	2	–	–	–	–	–	–	1	–	–	–	–
Tachycardia	–	–	–	–	–	–	–	1	–	1	–	–
Rigors	1	–	–	1	–	–	–	–	–	–	–	–
Feeling of warmth	1	–	–	1	–	–	–	–	–	–	–	–
Bone pain	–	–	–	–	–	–	–	–	–	–	–	1
Tumour pain	–	–	–	–	–	–	–	–	–	–	1	–
Cough	–	–	–	1	–	–	–	–	–	–	–	–
Oedema	–	–	–	1	–	–	–	–	–	–	–	–
Malaise	–	–	–	1	–	–	–	–	–	–	–	–




Figure 3 Serum concentrations of circulating extracellular domain of HER2 product (HER2 ECD). Open circles represent concentrations of HER2 ECD prior to the commencement of the treatment. Left and right ends of each horizontal bar represent minimum and maximum concentrations of circulating HER2 ECD during weekly administration

patients at the same dose level. Antibodies against MKC-454 were not detected in any patient.

From a total of 134 administrations of MKC-454, drug-related toxicity is shown in Table 3. A 38°C or higher fever was observed in six patients. One patient at 1 mg kg⁻¹ experienced grade 3 fever up to 40.7°C while receiving the first administration. Despite the temporary event, further treatment with MKC-454 was not given. Although one patient had no bone metastasis, the patient at the dose of 8 mg kg⁻¹ experienced severe generalized bone pain judged as grade 3. Although the pain subsided within 8 h, administration of MKC-454 was discontinued.

After the clinical trial had been completed, treatment responses of all the 18 patients were assessed by the committee confirming toxicity and efficacy. Data are summarized in Table 4. Although no response was observed at the dose level of 1 mg kg⁻¹, one patient with cervical and mediastinal lymph node metastases had a MR at the dose level of 2 mg kg⁻¹. One patient with multiple lung metastases had a PR at 4 mg kg⁻¹. Due to her good response to




Figure 4 Metastatic tumours in the lung before treatment with MKC-454 (A) and 7 months later showing marked response (B)

date, she has received continuous weekly administration of MKC-454 4 mg kg⁻¹ for 18 months (Figure 4). A patient at 8 mg kg⁻¹ with skin and bilateral supraclavicular lymph node metastases had a CR for 3 months (Figure 5). However, skin metastasis relapsed




Figure 5 Skin metastasis before treatment with MKC-454 (A) showing strong positivity for HER2 product (B) and 4 months later showing pathologically complete response (C). Relapsing tumour cells with positive HER2 product (D)

Table 4 Response to MKC-454

Dose mg kg ⁻¹	NE	PD	NC	MR	PR	CR	Total
1	1	4	1	–	–	–	6
2	–	2	–	1	–	–	3
4	–	2	–	–	1	–	3
8	1	1	1	2	–	1	6
Total	2	9	2	3	1	1	18

NE, not evaluable; PD, progressive disease; NC, no change; MR, minimal response; PR, partial response; CR, complete response.

during the course of the treatment. It is of some interest that the relapsed tumour cells were still strongly overexpressing the HER2 product. Her circulating shed antigen had been found low through the entire treatment course. Therefore it was considered that the circulating shed antigen had not affected her MKC-454 pharmacokinetics.

DISCUSSION

The HER2 product is thought to be a unique and useful target for antibody therapy of cancers overexpressing the HER gene. Several

series of murine mAbs directed against the extracellular domain of the HER2 gene product have been developed and examined to reveal their anti-tumour effects, mainly in vitro. However, human anti-mouse antibody response during the therapy would be a major limitation on the clinical application of such murine mAbs. Preclinical studies suggested that a humanized antibody would have the same or greater anti-tumour potency in clinical trials than its murine counterpart (Tokuda et al, 1996). Therefore, a humanized antibody against the HER2 gene product is expected to be useful in the clinic.

Baselga and his colleagues reported the first clinical evidence of the anti-tumour activity of a humanized antibody against the HER2 product (Baselga et al, 1996). They set 10 µg ml⁻¹ as the target trough serum concentration based on in vitro data (Carter et al, 1992). Our preclinical data also indicated 10 µg ml⁻¹ as a minimally effective concentration for anti-proliferation effects and antibody-dependent cell-mediated cytotoxicity (ADCC) (Tokuda et al, 1996). The dosage and schedule of the antibody administration in this clinical trial was based on unpublished phase I trials in the United States. In this clinical trial, non-linear pharmacokinetics of MCK-454 was demonstrated as previously described for other monoclonal antibodies (Koizumi et al, 1986; Eger et al, 1987). According to our data, 4 mg kg⁻¹ is sufficient to achieve the target

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.