


D I S P A T C H
 A MONORAIL PICCOLINO® PUBLICATION

Introducing the Schneider
MONORAIL®-GEX™
 Guidewire Exchange Catheter

In conjunction with the MONORAIL PICCOLINO® PTCA catheter, the MONORAIL® -GEX™ catheter will allow you to exchange or re-form the distal tip of the coronary guidewire without compromising coronary artery position. The peel-away slit through lumen will facilitate safe and easy exchange from the MONORAIL® -GEX™ to the MONORAIL PICCOLINO® Series P34.



ORDERING INFORMATION

Model Number	Shaft Diameter (Fr/mm)	Maximum Guidewire (in)	Usable Length (cm)
K53-002	3.6/1.20	.014	135

Schneider (USA) Inc. • Pfizer Hospital Products Group • 5905 Nathan Lane, Minneapolis, MN 55442 • Tel 612-650-5600 • Fax 612-650-5771
 Customer Service 800-822-4822 • 24 hour availability • 365 days a year
 Schneider (Europe) AG • Pfizer Hospital Products Group • Schärenmoosstrasse 115, CH-8052 Zurich, Switzerland • Tel 01-5012350 Fax 01-5012434
 Schneider Canada • 245 West Beaver Creek Road, Richmond Hill, Ontario L4B1L1, Canada • Tel 416-764-6367 • Fax 416-734-7101 • Customer Service 800-268-5625
 MONORAIL and MONORAIL PICCOLINO are registered trademarks of Schneider (Europe) AG. U.S. Patent Number 4,782,123.
 MONORAIL-GEX is a trademark of Schneider (USA) Inc.

© Schneider (USA) Inc. 10/90 460112-1090