


Expert opinion on investigational
drugs.
v. 8, no. 8 (Aug. 1999)
General Collection
W1 EX52M
Received: 12-30-1999



Ashley Publications

Expert Opinion on Investigational Drugs

Vol. 8 No. 8

Reviews & Updates

Triple nucleoside analogue antiretroviral therapy
Therapeutic potential of blocking HIV entry into cells
New research in macrolides and ketolides since 1997
Mechanism of action of the oxazolidinones
Antisense agents and bacterial infections
Nitric oxide scavengers


Drug Evaluations

Daptomycin
Tacrolimus
AIT-082

Meeting Highlights

Meeting of the American Society for Microbiology

PETITIONERS
ISSN 1354-3784



PROPERTY OF THE
NATIONAL
LIBRARY OF
MEDICINE

Page 1 of 1

The material on this page was copied from the collection of the National Library of Medicine by a third party and may be protected by U.S. Copyright law.



<http://www.ashley-pub.com>

Drug Evaluation

1. Introduction
 2. Microbiology
 - 2.1 Mechanism of action
 - 2.2 *In vitro* antibacterial spectrum
 - 2.3 Bactericidal activity
 - 2.4 Post-antibiotic and sub-MIC effects
 - 2.5 Resistance
 - 2.6 Antimicrobial interactions
 3. *In vivo* models
 4. Animal toxicology
 5. Clinical studies
 - 5.1 Clinical pharmacokinetics
 - 5.2 Clinical efficacy
 - 5.3 Clinical safety
 6. Conclusions
- Bibliography

Monthly Focus: Anti-infectives

Daptomycin: a novel agent for Gram-positive infections

Francis P Tally, Michael Zeckel, Margaret M Wasilewski, Claudio Carini, Cindy L Berman, George L Drusano & Frederick B Oleson, Jr

The alarming increase in the incidence of Gram-positive infections, including those caused by resistant bacteria, has sparked renewed interest in novel antibiotics. One such agent is daptomycin, a novel lipopeptide antibiotic with proven bactericidal activity *in vitro* against all clinically relevant Gram-positive bacteria. These include resistant pathogens, such as vancomycin-resistant enterococci (VRE), methicillin-resistant *Staphylococcus aureus* (MRSA), glycopeptide intermediately susceptible *Staphylococcus aureus* (GISA), coagulase-negative staphylococci (CNS) and penicillin-resistant *Streptococcus pneumoniae* (PRSP), for which there are very few therapeutic alternatives. Daptomycin provides rapid, concentration-dependent killing and a relatively prolonged concentration-dependent post-antibiotic effect *in vitro*. Spontaneous acquisition of resistance to daptomycin occurs rarely. Daptomycin exhibits linear pharmacokinetics, minimal accumulation with once-daily dosing, and low plasma clearance and volume of distribution. Phase II clinical trials indicate that daptomycin at doses of 2 mg/kg q24 h and 3 mg/kg q12 h is efficacious against skin and soft tissue infections and bacteremia, respectively. In addition, results in endocarditis suggested potential efficacy with higher doses. On the basis of clinical trials to date, it appears that daptomycin has an excellent safety profile, with the incidence and nature of serious adverse events comparable to those observed with conventional therapy. Adverse events associated with other classes of antimicrobials (nephrotoxicity, local irritation, ototoxicity, hypersensitivity, and gastrointestinal effects) were uncommon with daptomycin. Minimal skeletal muscle toxicity was seen at only the highest dose tested (4 mg/kg q12 h), predicted by elevations in serum creatinine phosphokinase, and readily reversible upon discontinuation of treatment. There were no signs of toxicity in cardiac or smooth muscle. Phase II and III clinical trials are underway to evaluate daptomycin for the treatment of Gram-positive bacteremia and complicated skin and soft tissue infections, respectively. Daptomycin holds promise as a rapidly acting and highly effective antibiotic for Gram-positive infections.


Keywords: antibiotics, bacteremia, coagulase-negative staphylococci, daptomycin, endocarditis, glycopeptide intermediately susceptible *Staphylococcus aureus*, Gram-positive bacteria, methicillin-resistant *Staphylococcus aureus*, penicillin-resistant *Streptococcus pneumoniae*, resistance, skin infection, soft tissue infection, vancomycin-resistant enterococci

Exp. Opin. Invest. Drugs (1999) 8(8):1223-1238

1223

1999 © Ashley Publications Ltd. ISSN 1354-3784

Figure 1: Amino acid structure and location of decanoic acid side-chain of daptomycin.



1. Introduction

Daptomycin is a unique cyclic lipopeptide antibiotic derived from the fermentation of a strain of *Streptomyces roseosporus*. It is comprised of a decanoyl side-chain linked to the N-terminal tryptophan of a cyclic 13-amino acid peptide (**Figure 1**). Daptomycin is an antimicrobial agent with bactericidal activity against all clinically important Gram-positive bacteria, including resistant pathogens for which there are very limited therapeutic alternatives. These resistant pathogens include VRE, MRSA, GISA, CNS, and PRSP [1-6].

Daptomycin was discovered in the early 1980s by scientists at Eli Lilly and Company (Lilly) and was subsequently developed as an iv. drug for the treatment of serious Gram-positive infections. Nineteen Phase I and two Phase II clinical studies involving more than 370 subjects were conducted in the 1980s and early 1990s. The results with skin and soft tissue infections and bacteremia were highly encouraging. Results in the treatment of endocarditis suggest potential efficacy at higher doses. At a dose higher than those tested in Phase II trials, mild, reversible skeletal muscle toxicity was observed, prompting Lilly in 1991 to suspend voluntarily clinical investigation of daptomycin. This decision was made prior to the recent marked increase in the incidence of Gram-positive infections and bacterial antibiotic resistance to current therapies.

In 1997, Cubist Pharmaceuticals, Inc. (Cubist) licensed worldwide rights for daptomycin from Lilly. In today's environment of increasingly prevalent antimicrobial resistance, daptomycin may provide greater utility than existing therapies. Cubist is conducting Phase II and III trials evaluating daptomycin in bacteremia and complicated skin and soft tissue infections, respectively. Doses used in these trials are lower than those previously associated with transient skeletal muscle toxicity. In addition a Phase II study in hospitalised patients with Gram-positive urinary tract infections is planned. The results of these trials should be available at the end of the year 2000, and daptomycin may be approved for clinical use as early as 2001.

2. Microbiology

2.1 Mechanism of action

Although the precise mechanism of action of daptomycin is not completely understood, it is clearly distinct from that of other antibiotics, including β -lactams, aminoglycosides, glycopeptides and macrolides. Daptomycin kills Gram-positive bacteria by disrupting multiple aspects of bacterial plasma membrane function, while not penetrating into the cytoplasm. Potential antibacterial mechanisms include inhibition of peptidoglycan synthesis, inhibition of lipoteichoic acid synthesis [7-9] and dissipation of bacterial membrane potential [10-11].

2.2 *In vitro* antibacterial spectrum

The spectrum of activity for daptomycin includes most clinically significant Gram-positive bacteria; it does not act against Gram-negative bacteria. **Table 1** summarises the activity of daptomycin vs. vancomycin against recent (1996 - 1998) Gram-positive isolates chosen for their diverse susceptibilities to other antibiotics [6,12-15]. The clinical isolates were obtained from diverse locations, representing all geographic areas in the USA. Minimum inhibitory concentration (MIC)₉₀ values for daptomycin are below the proposed break point of 2 $\mu\text{g/ml}$ for staphylococci (including methicillin-susceptible *S. aureus*, MRSA, *Staphylococcus epidermidis* and CNS), streptococci (including penicillin-susceptible and penicillin-resistant *S. pneumoniae*, *Streptococcus pyogenes*, *Streptococcus agalactiae*, and viridans streptococci), enterococci (including vancomycin-susceptible and vancomycin-resistant


Table 1: *In vitro* activity against Gram-positive bacteria: daptomycin vs. vancomycin.

Micro-organism	N Total = 973	MIC range (µg/ml)	
		Daptomycin	Vancomycin
Classically 'susceptible' strains			
<i>S. aureus</i> (MSSA) [6,12,14]	134	0.03 - 2	0.25 - 2
<i>S. epidermidis</i> (MSSE) [12,14]	54	0.06 - 0.5	0.25 - 2
<i>Staphylococcus</i> spp. (Coag. Neg.) [6,14]	37	0.25 - 2	2 - 4
<i>S. pneumoniae</i> [6]	15	0.0075 - 0.06	0.015 - 5
<i>Staphylococcus haemolyticus</i> [14]	20	0.03 - 1	0.5 - 4
<i>S. pyogenes</i> [6,14]	61	0.015 - 0.5	0.25 - 1
<i>E. faecalis</i> [6,14]	56	0.03 - 2	0.5 - 4
<i>E. faecium</i> [14]	28	0.5 - 2	0.5 - 4
<i>S. agalactiae</i> [14]	31	0.12 - 0.25	0.25 - 0.5
Streptococci Group C [14]	9	0.03 - 0.25	0.25 - 0.5
Streptococci Group G [14]	10	0.015 - 0.06	0.25 - 0.5
Streptococci viridans [14]	37	0.12 - 2	0.5 - 1
<i>L. monocytogenes</i> [14]	25	2 - 8	1 - 1
<i>Corynebacterium jeikeium</i> [14]	21	0.12 - 0.5	0.12 - 1
<i>Clostridium clostridioforme</i> [15] [§]	11	1 - 8	0.5 - 1
<i>C. difficile</i> [15]	10	2 - 2	0.5 - 2
<i>Clostridium innocuum</i> [15]	11	8 - > 16	8 - 16
<i>Clostridium perfringens</i> [15]	10	0.5 - 2	1 - 1
<i>Clostridium ramosum</i> [15]	10	> 16	4 - 4
<i>Peptostreptococcus</i> spp. [15]	10	0.5 - 2	0.25 - 0.5
'Intermediately susceptible' strains			
<i>S. aureus</i> (GISA) [13]	8	0.25 - 2	2 - 8
<i>S. epidermidis</i> (GISE) [13]	3	0.5 - 1	8
<i>S. haemolyticus</i> [13]	1	1	8
Classically 'resistant' strains			
<i>S. aureus</i> (MRSA) [6,12,14]	155	0.06 - 2	0.25 - 2
<i>S. epidermidis</i> (MRSE) [12]	29	0.06 - 1	0.5 - 2
<i>Staphylococcus</i> spp. (MRS) [14]	56	0.004 - 1	0.25 - 4
<i>S. pneumoniae</i> (PRSP) [6]	40	< 0.03 - 1	0.06 - 2
<i>E. faecalis</i> (VRE) [14]	14	0.5 - 1	8 - > 64
<i>E. faecium</i> (VRE) [6,14]	67	0.5 - 2	> 64 - 1024

[§] Agar dilution method, 25 mg/l Ca²⁺.

Figure 2: *In vitro* kill rates: daptomycin vs. vancomycin [16].

In vitro kill rates for *S. aureus* and *E. faecalis*. Daptomycin minimum inhibitory concentration (MIC) = 1 µg/ml; vancomycin MIC = 4 µg/ml. (Adapted from HANBERGER H, NILSSON LE, MALLER R, ISAKSSON B. *Antimicrob. Agents Chemother.* (1991) 35:1710-1716. Reprinted with permission. © American Society for Microbiology.)



strains), *Clostridium difficile*, and *Propionibacterium acnes*. Antibacterial activity against classically 'resistant' strains is comparable to that against classically 'susceptible' strains. In addition, the MIC value for daptomycin against susceptible strains is typically 4-fold lower than that of vancomycin. These results are in agreement with previous studies, with the exception of the findings for *Listeria monocytogenes*, which in earlier studies appeared to be more susceptible [1,3].

2.3 Bactericidal activity

Unlike glycopeptide antibiotics, daptomycin exhibits rapid, concentration-dependent bactericidal activity *in vitro* against Gram-positive organisms, including enterococci. This has been demonstrated with both time-kill curves and broth dilution methodology. Using the standard definition of a 3-log reduction in viable organisms, daptomycin, but not vancomycin, is bactericidal against both *S. aureus* and *Enterococcus faecalis* (Figure 2) [16]. The *S. aureus* initial kill rate of daptomycin is exceptionally rapid: a greater than 3-log₁₀ reduction in viable organisms is typically achieved in less than 1 h at daptomycin concentrations 4 - 8 times the MIC value [6,16-17]. In contrast, vancomycin typically takes 6 - 24 h to achieve a 3-log₁₀ kill at equivalent concentrations. Faster bactericidal activity has been demonstrated with daptomycin vs. vancomycin in logarithmic as well as stationary phases of bacterial growth [18]. Daptomycin also is rapidly bactericidal *in vitro* against

enterococci (2 h for a 3-log₁₀ kill) whereas vancomycin generally does not exhibit bactericidal activity against enterococci [19-26].

2.4 Post-antibiotic and sub-MIC effects

Daptomycin produces a post-antibiotic effect (PAE), regrowth times and sub-MIC effects *in vitro* that are prolonged and concentration-dependent. PAEs lasting from 1 - 6 h were observed against *E. faecalis* and *S. aureus* following exposure to daptomycin concentrations ranging from 0.25 - 16 µg/ml (from 1- to 8-fold the MIC value). Using viable-cell counts, exposure to 15 µg/ml daptomycin for 2 h produced PAEs of 2.4 - 5.3 h and 3.5 - 3.9 h in four clinical *S. aureus* isolates and two clinical *E. faecalis* isolates, respectively [27]. A PAE of 6.3 - 6.7 h was observed for strains of both *S. aureus* and *E. faecalis* following use of bioluminescent techniques [16]. Effective regrowth time following the PAE appears to be concentration-dependent [28]. Daptomycin also exerts antimicrobial effects at sub-MIC concentrations. In one study, *S. aureus* exposed to daptomycin at concentrations one-quarter the MIC level were phagocytosed and killed in significantly greater numbers than were non-exposed bacteria [29].

2.5 Resistance

Spontaneous acquisition of resistance to daptomycin is rare. No spontaneously resistant mutants were obtained for any Gram-positive organism tested when challenged at eight times the MIC value (resistance

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.