

CPMP/ICH/380/95

ICH Topic Q 1 A Stability Testing Guidelines: Stability Testing of New Drug Substances and Products

Step 5

NOTE FOR GUIDANCE ON STABILITY TESTING: STABILITY TESTING OF NEW DRUG SUBSTANCES AND PRODUCTS (CPMP/ICH/380/95)

[This guideline replaces relevant section of previous guideline]

APPROVAL BY CPMP	December 1993
STUDIES CAN BE SUBMITTED ACCORDING TO THIS GUIDELINE	January 1994
DATE FOR COMING INTO OPERATION (DOSSIERS SUBMITTED AFTER)	1 January 1998

ICH - Technical Coordination - R. Bass
7 Westferry Circus, Canary Wharf, London E14 4HB, UK
Tel: (+44-171) 418 84 11 Fax: (+44-171) 418 85 51
E_Mail: mail@emea.eudra.org http://www.eudra.org/emea.html



STABILITY TESTING OF NEW DRUG SUBSTANCES AND PRODUCTS

ICH Harmonised Tripartite Guideline

[EMEA Status as of December 1993]

Preamble

The following guideline sets out the stability testing requirement for a Registration Application within the three areas of the EC, Japan and the USA. It does not seek necessarily to cover the testing that may be required for registration in or export to other areas of the world.

The guideline seeks to exemplify the core stability data package required for new drug substances and products. It is not always necessary to follow this when there are scientifically justifiable reasons for using alternative approaches.

The guideline provides a general indication on the requirements for stability testing, but leaves sufficient flexibility to encompass the variety of different practical situations required for specific scientific situations and characteristics of the materials being evaluated.

The principle that information on stability generated in any one of the three areas of the EC, Japan and the USA would be mutually acceptable in both of the other two areas has been established, provided it meets the appropriate requirements of this guideline and the labelling is in accord with national/regional requirements.

Details of the specific requirements for sampling, test requirements for particular dosage forms/packaging etc., are not covered in this guideline.

Objective

The purpose of stability testing is to provide evidence on how the quality of a drug substance or drug product varies with time under the influence of a variety of environmental factors such as temperature, humidity and light, and enables recommended storage conditions, re-test periods and shelf lives to be established.

Scope

The guideline primarily addresses the information required in Registration Applications for new molecular entities and associated drug products.

This guideline does not currently seek to cover the information required for abbreviated or abridged applications, variations, clinical trial applications, etc.

The choice of test conditions defined in this guideline is based on an analysis of the effects of climatic conditions in the three areas of the EC, Japan and the USA. The mean kinetic temperature in any region of the world can be derived from climatic data (Grimm, W. Drugs Made in Germany, 28, 196-202, 1985 and 29, 39-47, 1986).

CPMP/ICH/380/95



DRUG SUBSTANCE

General

Information on the stability of the drug substance is an integral part of the systematic approach to stability evaluation.

Stress Testing

Stress testing helps to determine the intrinsic stability of the molecule by establishing degradation pathways in order to identify the likely degradation products and to validate the stability indicating power of the analytical procedures used.

Formal Studies

Primary stability studies are intended to show that the drug substance will remain within specification during the re-test period if stored under recommended storage conditions.

Selection of Batches

Stability information from accelerated and long term testing is to be provided on at least three batches. The long term testing should cover a minimum of 12 months duration on at least three batches at the time of submission.

The batches manufactured to a minimum of pilot plant scale should be by the same synthetic route and use a method of manufacture and procedure that simulates the final process to be used on a manufacturing scale.

The overall quality of the batches of drug substance placed on stability should be representative of both the quality of the material used in pre-clinical and clinical studies and the quality of material to be made on a manufacturing scale.

Supporting information may be provided using stability data on batches of drug substance made on a laboratory scale.

The first three production batches of drug substance manufactured post approval, if not submitted in the original Registration Application, should be placed on long term stability studies using the same stability protocol as in the approved drug application.

Test Procedures and Test Criteria

The testing should cover those features susceptible to change during storage and likely to influence quality, safety and/or efficacy. Stability information should cover as necessary the physical, chemical and microbiological test characteristics. Validated stability-indicating testing methods must be applied. The need for the extent of replication will depend on the results of validation studies.

Specification

Limits of acceptability should be derived from the profile of the material as used in the preclinical and clinical batches. It will need to include individual and total upper limits for impurities and degradation products, the justification for which should be influenced by the levels observed in material used in pre-clinical studies and clinical trials.

CPMP/ICH/380/95



Storage Conditions

The length of the studies and the storage conditions should be sufficient to cover storage, shipment and subsequent use. Application of the same storage conditions as applied to the drug product will facilitate comparative review and assessment. Other storage conditions are allowable if justified. In particular, temperature sensitive drug substances should be stored under an alternative, lower temperature condition which will then become the designated long term testing storage temperature. The six months accelerated testing should then be carried out at a temperature at least 15°C above this designated long term storage temperature (together with the appropriate relative humidity conditions for that temperature). The designated long term testing conditions will be reflected in the labelling and re-test date.

Conditions Minimum Time

Period at Submission

Long term testing $25^{\circ}\text{C} \pm 2^{\circ}\text{C} / 60\% \text{ RH} \pm 5\%$ 12 Months

Accelerated Testing $40^{\circ}\text{C} \pm 2^{\circ}\text{C}/75\% \text{ RH} \pm 5\%$ 6 Months

Where 'significant change' occurs during six months storage under conditions of accelerated testing at $40^{\circ}\text{C} \pm 2^{\circ}\text{C}/75$ percent RH \pm 5 percent, additional testing at an intermediate condition (such as $30^{\circ}\text{C} \pm 2^{\circ}\text{C}/60$ percent RH \pm 5 percent) should be conducted for drug substances to be used in the manufacture of dosage forms tested long term at 25°C/60 percent RH and this information included in the Registration Application. The initial Registration Application should include a minimum of 6 months data from a 12 months study.

'Significant change' at 40°C/75 percent RH or 30°C/60 percent RH is defined as failure to meet the specification.

The long term testing will be continued for a sufficient period of time beyond 12 months to cover all appropriate re-test periods, and the further accumulated data can be submitted to the Authorities during the assessment period of the Registration Application.

The data (from accelerated testing or from testing at an intermediate condition) may be used to evaluate the impact of short term excursions outside the label storage conditions such as might occur during shipping.

Testing Frequency

Frequency of testing should be sufficient to establish the stability characteristics of the drug substance. Testing under the defined long term conditions will normally be every three months, over the first year, every six months over the second year and then annually.

Packaging/Containers

The containers to be used in the long term, real time stability evaluation should be the same as or simulate the actual packaging used for storage and distribution.

CPMP/ICH/380/95



Evaluation

The design of the stability study is to establish, based on testing a minimum of three batches of the drug substance and evaluating the stability information (covering as necessary the physical, chemical and microbiological test characteristics), a retest period applicable to all future batches of the bulk drug substance manufactured under similar circumstances. The degree of variability of individual batches affects the confidence that a future production batch will remain within specification until the retest date.

An acceptable approach for quantitative characteristics that are expected to decrease with time is to determine the time at which the 95% one-sided confidence limit for the mean degradation curve intersects the acceptable lower specification limit. If analysis shows that the batch to batch variability is small, it is advantageous to combine the data into one overall estimate and this can be done by first applying appropriate statistical tests (for example, p values for level of significance of rejection of more than 0.25) to the slopes of the regression lines and zero time intercepts for the individual batches. If it is inappropriate to combine data from several batches, the overall retest period may depend on the minimum time a batch may be expected to remain within acceptable and justified limits.

The nature of any degradation relationship will determine the need for transformation of the data for linear regression analysis. Usually the relationship can be represented by a linear, quadratic or cubic function on an arithmetic or logarithmic scale. Statistical methods should be employed to test the goodness of fit of the data on all batches and combined batches (where appropriate) to the assumed degradation line or curve.

The data may show so little degradation and so little variability that it is apparent from looking at the data that the requested retest period will be granted. Under the circumstances, it is normally unnecessary to go through the formal statistical analysis but merely to provide a full justification for the omission.

Limited extrapolation of the real time data beyond the observed range to extend expiration dating at approval time, particularly where the accelerated data supports this, may be undertaken. However, this assumes that the same degradation relationship will continue to apply beyond the observed data and hence the use of extrapolation must be justified in each application in terms of what is known about the mechanism of degradation, the goodness of fit of any mathematical model, batch size, existence of supportive data etc.

Any evaluation should cover not only the assay, but the levels of degradation products and other appropriate attributes.

Statements/ Labelling

A storage temperature range may be used in accordance with relevant national/regional requirements. The range should be based on the stability evaluation of the drug substance. Where applicable, specific requirements should be stated, particularly for drug substances that cannot tolerate freezing. The use of terms such as 'ambient conditions' or 'room temperature' is unacceptable.

A re-test period should be derived from the stability information.

CPMP/ICH/380/95



DOCKET

Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.

