

Guidance for Industry and for FDA Reviewers

**Guidance Document for Premarket
Notification Submissions for Nitric Oxide
Delivery Apparatus, Nitric Oxide Analyzer
and Nitrogen Dioxide Analyzer**

Document Issued on: January 24, 2000



U.S. Department of Health and Human Services
Food and Drug Administration
Center for Devices and Radiological Health

Anesthesiology, Respiratory, and Defibrillator Devices Group
Division of Cardiovascular, Respiratory, and Neurological Devices
Office of Device Evaluation

Preface

Public Comment

For 90 days following the date of publication in the Federal Register of the notice announcing the availability of this guidance, comments and suggestions regarding this document should be submitted to Docket No. 99D-5297, Dockets Management Branch, Division of Management Systems and Policy, Office of Human Resources and Management Services, Food and Drug Administration, 5630 Fishers Lane, (HFA-305), Room 1061, Rockville, MD 20852.. Such comments will be considered when determining whether to amend the current guidance.

After 90 days following the date of publication in the Federal Register of the notice announcing the availability of this guidance comments and suggestions may be submitted at any time for Agency consideration to: Michael Bazaral, M.D., Ph.D., Center for Devices and Radiological Health (HFZ-450), Food and Drug Administration, 9200 Corporate Blvd., Rockville, MD 20850. Comments may not be acted upon by the Agency until the document is next revised or updated. For questions regarding the use or interpretation of this guidance contact Michael Bazaral, M.D., Ph.D. at 301-443-8609.

Additional Copies

World Wide Web/CDRH/home page at <http://www.fda.gov/cdrh/ode/1157.pdf> or CDRH Facts on Demand at 1-800-899-0381 or 301-827-0111, specify number 1157 when prompted for the document shelf number.

Table of Contents

1. INTRODUCTION.....	1
1.1 PURPOSE.....	1
1.2 BACKGROUND	1
1.3 SCOPE	1
2. DEVICE DESCRIPTION	3
2.1 NITRIC OXIDE ADMINISTRATION APPARATUS.....	3
2.2 NITRIC OXIDE GAS ANALYZER.....	3
2.3 NITROGEN DIOXIDE GAS ANALYZER.....	4
3. SPECIFIC CRITERIA AND TESTING	5
3.1 NITRIC OXIDE DELIVERY APPARATUS	5
3.1.1 <i>Loss of nitric oxide therapy and incorrect nitric oxide concentration.....</i>	<i>5</i>
3.1.2 <i>Insufficient or excess ventilation or oxygenation.....</i>	<i>8</i>
3.1.3 <i>Excessive nitrogen dioxide administration.....</i>	<i>9</i>
3.1.4 <i>Potential for catastrophic release of nitric oxide.....</i>	<i>10</i>
3.1.5 <i>Adulteration of the nitric oxide.....</i>	<i>11</i>
3.1.6 <i>Electrical hazards</i>	<i>11</i>
3.1.7 <i>Adverse effects on other electronic devices.....</i>	<i>11</i>
3.1.8 <i>Release of nitric oxide and release and generation of nitrogen dioxide.....</i>	<i>11</i>
3.2 NITRIC OXIDE ANALYZER.....	12
3.2.1 <i>Nitric Oxide measurement error.....</i>	<i>12</i>
3.2.2 <i>Electrical hazards</i>	<i>14</i>
3.2.3 <i>Adverse effects on other electronic devices.....</i>	<i>14</i>
3.3 NITROGEN DIOXIDE ANALYZER.....	14
3.3.1 <i>Nitrogen Dioxide measurement error.....</i>	<i>15</i>
3.3.2 <i>Electrical hazards</i>	<i>17</i>
3.3.3 <i>Adverse effects on other electronic devices.....</i>	<i>17</i>
4. GENERAL CRITERIA AND TESTING.....	18
4.1 GENERAL CRITERIA	18
4.2 GENERAL TEST METHODS	18
5. ELECTRICAL SAFETY.....	20
5.1 PERFORMANCE CRITERIA	20
5.1.1 <i>Battery power.....</i>	<i>20</i>
5.1.2 <i>Electrical power indicators.....</i>	<i>20</i>
5.1.3 <i>Overcurrent protection.....</i>	<i>21</i>
5.1.4 <i>Dielectric withstand.....</i>	<i>21</i>
5.1.5 <i>AC power grounding and polarity.....</i>	<i>21</i>
5.1.6 <i>Leakage current.....</i>	<i>21</i>
5.1.7 <i>Auxiliary output.....</i>	<i>22</i>
5.2 TEST METHODS	22
5.2.1 <i>Battery power.....</i>	<i>22</i>
5.2.2 <i>Electrical power indicators.....</i>	<i>23</i>
5.2.3 <i>Overcurrent protection.....</i>	<i>23</i>

5.2.4	<i>Dielectric withstand</i>	24
5.2.5	<i>AC power grounding and polarity</i>	24
5.2.6	<i>Leakage current</i>	24
5.2.7	<i>Auxiliary output</i>	24
6.	ELECTROMAGNETIC COMPATIBILITY	25
6.1	PERFORMANCE CRITERIA	25
6.1.1	<i>Emissions</i>	25
6.1.1.1	Radiated and conducted electromagnetic energy.....	25
6.1.1.2	Magnetic fields.....	25
6.1.2	<i>Immunity</i>	25
6.1.2.1	Electrostatic discharge	26
6.1.2.2	Radiated electromagnetic fields	26
6.1.2.3	AC voltage fluctuations, transients, and surges	26
6.1.2.4	Conducted electromagnetic energy.....	27
6.1.2.5	Magnetic fields.....	28
6.1.2.6	Quasi-static electric fields	28
6.2	TEST METHODS	28
6.2.1	<i>Emissions</i>	28
6.2.1.1	Radiated and conducted electromagnetic energy.....	28
6.2.1.2	Magnetic fields.....	29
6.2.2	<i>Immunity</i>	29
6.2.2.1	Electrostatic discharge	29
6.2.2.2	Radiated electromagnetic fields	30
6.2.2.3	AC voltage fluctuations, transients, and surges	33
6.2.2.4	Conducted electromagnetic energy.....	35
6.2.2.5	Magnetic fields.....	36
6.2.2.6	Quasi-static electric fields	36
7.	PERFORMANCE SPECIFICATIONS, ENVIRONMENTAL AND MECHANICAL SAFETY	38
7.1	PERFORMANCE CRITERIA	38
7.1.1	<i>Controls protection</i>	38
7.1.2	<i>Connector protective incompatibility</i>	38
7.1.3	<i>Mechanical safety</i>	38
7.1.4	<i>Mechanical vibration and shock resistance</i>	39
7.1.5	<i>Fluid spill resistance</i>	39
7.1.6	<i>High and low temperature and humidity</i>	39
7.1.7	<i>Surface temperature</i>	39
7.1.8	<i>Toxic materials</i>	40
7.1.9	<i>Strangulation</i>	40
7.1.10	<i>Determination of Endurance</i>	40
7.1.11	<i>Material Compatibility</i>	40
7.1.12	<i>Medical Gas Cylinder Connections</i>	40
7.2	TEST METHODS	40
7.2.1	<i>Controls protection</i>	40
7.2.2	<i>Connector protective incompatibility</i>	40
7.2.3	<i>Mechanical safety</i>	41
7.2.4	<i>Mechanical vibration and shock resistance</i>	41
7.2.5	<i>Fluid spill resistance</i>	41

7.2.6	<i>High and low temperature and humidity</i>	41
7.2.7	<i>Surface temperature</i>	42
7.2.8	<i>Toxic materials</i>	42
7.2.9	<i>Strangulation</i>	42
7.2.10	<i>Determination of Endurance</i>	42
7.2.11	<i>Material Compatibility</i>	42
8.	HARDWARE DOCUMENTATION	43
9.	SOFTWARE DOCUMENTATION	44
10.	LABELING	45
10.1	IDENTIFICATION OF MEDICAL GAS CYLINDERS AND CONNECTIONS	45
10.2	INSTRUCTIONS FOR USE.....	45
10.2.1	<i>Intended Use</i>	45
10.2.2	<i>Validated Ventilators</i>	45
10.2.3	<i>Installation Instructions</i>	46

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.