

Technical Quality Dextran

Technical Quality Dextran is a high purity dextran fraction with selected average molecular weight and molecular weight distribution.

Easily soluble Pharmacosmos Technical Quality Dextran fractions are used as starting or intermediate reagents in several different processes in biotech, photographic and chemical manufacturing industries.

- [Dextran T1](#)
- [Dextran T1.5](#)
- [Dextran T3.5](#)
- [Dextran T5](#)
- [Dextran T6](#)
- [Dextran T10](#)
- [Dextran T20](#)
- [Dextran T25](#)
- [Dextran T40](#)
- [Dextran T60](#)
- [Dextran T70](#)
- [Dextran T110](#)
- [Dextran T150](#)
- [Dextran T250](#)
- [Dextran T500](#)
- [Dextran T750](#)
- [Dextran T2000](#)

Suggested reading

5/27/2015

Technical Quality Dextran | Pharmacosmos Dextran

;

Dextran T1 Technical Quality

Dextran T1

Technical Quality Dextran T1 is a high purity dextran fraction with selected average molecular weight and molecular weight distribution.



- [Documentation](#)
- [Item no. & Packaging](#)
- [Use of Dextran](#)

If you require a dextran product manufactured according to Good Manufacturing Practice (GMP), we recommend Pharmacosmos' [Pharmaceutical Quality Dextrans](#).

Download Section

[BRC Dextran T1](#)

[Batch Release Certificate | Dextran 1](#)

[GMP Active Substances Human | April 2014](#)

[Statement on Good Manufacturing Practices for Dextran API](#)

[MSDS on Dextran Powder](#)

Safety Data Sheet

Technical Quality Dextran from Pharmacosmos are packed and stored in safe and clean conditions.

Technical Quality Dextran T1 is available in the following pack sizes:

Item no	Weight	Packaging
5510 0001 9006	100 g.	Polyethylene container, FDA registered, DMF type III.
5510 0001 9007	500 g.	Polyethylene container, FDA registered, DMF type III.
5510 0001 9001	5 kg.	Polypropylene container. Dextran powder sealed in double polyethylene bags
5510 0001 9005	50 kg.	Fiber drum with the dextran powder sealed in double polyethylene bags.

Pharmacosmos supplies Technical Quality Dextran to all countries. The product packing meets all national and international regulations.

Technical Quality Dextran fractions are used as starting or intermediate reagents in several different processes in biotech, photographic and chemical manufacturing industries.

Suggested reading

Batch Release Certificate

Product name: Dextran T1

Specification No.: 40036

Batch No.: ABxxx

Manufacturing sites: Pharmacosmos A/S, Roervangsvej 30, DK-4300 Holbaek, Denmark
 DKMA* No.: 254629
 GMP certificate No.: N/A
 FDA establishment No.: N/A
 FDA facility classification: N/A
 EDQM* certificate No.: N/A

Method:	Parameter:	Results of analysis:	Limits:
–	Description:	Complies	White Powder
LI030-1	Average molecular mass, Mw:	xx,100	800 – 1,200
DF001	Loss on drying (105°C, 5h), % w/w:	x.x	≤7
DF030	Color of solution (Absorbance at 375 nm, 6% sol., 1 cm):	x.xx	≤0.12
DF005	Acidity or Alkalinity:	Complies	Complies
DF009	Specific rotation (+/-) °:	+x.x	+140 – +170
DF012	Heavy metals, ppm lead:	x	≤10
DF019	Nitrogen containing impurities, ppm N:	x	≤100

*) EDQM refers to 'European Directorate for the Quality of Medicines and healthcare'. DKMA refers to 'Danish Medicines Agency'.

We hereby confirm that no class 1, class 2 and class 3 solvent, cf. ICH Q3C and VICH GL 18, is used in the manufacturing of this product.

I hereby certify that the above information is authentic and accurate. This batch has been manufactured, including packaging and quality control in full compliance with the above mentioned specifications. NOTE: This product is not manufactured according to Good Manufacturing Practice (GMP).

Date (dd.mm.yyyy):

 Hans Juhl, M.Sc., Quality Control

 PHARMACOSMOS A/S
 ROERVANGSVEJ 30
 DK-4300 HOLBAEK
 DENMARK

 TELEPHONE: +45 59 48 59 59
 TELEFAX: +45 59 48 59 60
 E-MAIL: info@pharmacosmos.com
 WEB-SITE: www.pharmacosmos.com

 CVR NO.: DK15517085
 VAT NO. EXPORT: DK29127204

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.