

SETTING THE STANDARD WITH A COMPREHENSIVE FAMILY OF PORTS



DOCKET A L A R M



Find authenticated court documents without watermarks at <u>docketalarm.com</u>.

As a single-source provider, Bard Access Systems sets the standard for implanted ports and vascular access devices. We make it possible for you to address the unique clinical requirements of individual cases by offering you a comprehensive family of implanted ports, which include:


- More than 100 port configurations to address your access needs
- Thirteen distinct port designs
- Three catheter types: Groshong,[®] Hickman[®] and ChronoFlex[®] polyurethane
- Catheter sizes ranging from 6.0 12.0 French with single- or dual-lumen configurations
- Attachable or pre-attached catheters
- Three kit options on most ports: Basic, Intermediate and MicroIntroducer
- The convenience of BardPort's Cath-Lock compression lock for easy catheter assembly
- More BardPorts[™] are placed worldwide than any other port
- Ninety-eight percent of the top 50 cancer facilities in the United States use BardPorts

DOCKET

Δ

X-Port *isp*[™] Implanted Port

- · Compact, intermediate-sized port.
- · Large septum is easy to palpate and provides easy access during needle insertion.
- Silicone-filled suture holes reduce the risk of tissue in-growth and facilitate port removal.
- Tapered design eases port introduction into surgical pocket.
- Plastic material reduces MRI distortion.
- · Small footprint enables small port pocket.
- The firmness and rounded tip of Chronoflex polyurethane catheters facilitate advancement into vein.





X-Port isp M.R.I. Plastic, Single-Lumen **Implanted Ports with Attachable Groshong Catheter**

Port with MicroIntroducer Kit

DOCKE.

Δ

7757540	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable Groshong 8.0 French Single-Lumen Venous Catheter	
Port with	Intermediate Kit	
7707540	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable Groshong 8.0 French Single-Lumen Venous Catheter	
X-Port <i>isp</i> M.R.I. Plastic, Open-Ended, Single-Lumen Implanted Ports		
Ports with MicroIntroducer Kit		

0657500	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter	P 0
0657510	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 8.0 French Single-Lumen Venous Catheter	06
0657520	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Pre-Attached Open-Ended Silicone 9.6 French Single-Lumen Venous Catheter	06
0657525	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable Open-Ended Silicone 9.6 French Single-Lumen Venous Catheter	06

Large septum for easy needle insertion.

Tapered port design for ease of insertion.

Radiopaque ring ensures accurate visualization of

locking mechanism.

Cath-Lock mechanism easily secures catheter.

Silicone-filled suture holes secure and stabilize port.

Ports with Intermediate Kit X-Port isn M.R.I. Plastic Impla


0607520

0607530	Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter
0607540	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 8.0 French Single-Lumen Venous Catheter
0607550	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Pre-Attached Open-Ended Silicone 9.6 French Single-Lumen Venous Catheter
0607555	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable Open-Ended Silicone 9.6 French Single-Lumen Venous Catheter
Ports with	Basic Tray
Ports with 0607500	Basic Tray X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter
	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 6.0
0607500	X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter X-Port <i>isp</i> M.R.I. Plastic Implanted Port with Attachable ChronoFlex Polyurethane 8.0

X-Port isp Port

X-Port[™] Implanted Port

- Sized for arm placement.
- Ideal for chest placement in smaller patients.
- Tapered design eases port introduction into surgical pocket.
- Low-profile allows for a smaller port pocket.
- Rear suture hole placement facilitates easy suturing for secure placement.
- The firmness and rounded tip of the ChronoFlex polyurethane catheter facilitate advancement into the vein.
- The Cath-Lock mechanism is engineered for security and ease of attachment.
- Radiopaque ring ensures accurate visualization of locking mechanism via fluoroscopy or X-ray.



Radiopaque ring ensures accurate visualization of locking mechanism.

Rear suture holes secure and stabilize port.

ChronoFlex polyureth catheter provides sm outside diameter and maximum flow rates.

Cath-Lock mechanism easily secures cather

Silicone-filled suture holes

Stable base secures placement.

Titanium X-Port Single-Lumen Implanted Ports

Port with MicroIntroducer Kit

0655870 X-Port Titanium Implanted Port with Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter

Large, elevated septum guides needle insertion.

Tapered port design


for ease of insertion.

Port with Intermediate Kit

0605850 X-Port Titanium Implanted Port with Attachable ChronoFlex Polyurethane 6.0 French Single-Lumen Venous Catheter

Port with Basic Tray

DOCKE.



Find authenticated court documents without watermarks at docketalarm.com.

DOCKET A L A R M



Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.