

5/26/52

① will DDA stick to plastic.

1	2	3	4	5	6	7	8	9
	Ad 2 DDA - biot. w R7							
	denatured in				- diluted w/o carrier			
					0.3M NaOH			10ng → 8 serial
					Neut. to 1M Amm acetate			dilutions
								w/o
								carrier
	on polylysine				on plate			
	coated plate				straight			
					denatured in			
					0.3M NaOH			
					on plate			
					straight			

Goat Ab.
Rabbit Goat - APase
Tetrazolium

polylysine coated plate - buffer turns purple
red peaks come off smeared spots.

NaOH - good
NaOH + Amm. Acet. better } no bkgnd

Dot blots.
Bkgnd in buffer
= polylysine.

• • • 2.5

• • • 3 0.150

• • • 37, 0.018, 0.009

① Double paper may be better

↑ it has 2 where no paper found. [Signature]

5/26/82

EFFICIENCY LINE # 22-206

Ad 2. NaOH - AmmAcet
 in Elisa plates
 PNPP
 read.
fantastic.

	1	2	3	4	5	6	7	8	9
1									
2									
3									
4									
5									
6									
7									
8									
9									
10	1	A	BLANK						
11	1	B	0.006						
12	1	C	0.012						
13	1	D	0.009						
14	1	E	0.013						
15	1	F	0.006						
16	1	G	0.007						
17	1	H	0.005						
18			<i>NaOH</i>						
19	2	A	0.268						
20	2	B	0.194						
21	2	C	0.182						
22	2	D	0.284						
23	2	E	0.471						
24	2	F	0.301						
25	2	G	0.697						
26	2	H	0.838						
27			<i>NaOH + Amm Acet.</i>						
28	3	A	OVER						
29	3	B	OVER						
30	3	C	OVER						
31	3	D	OVER						
32	3	E	OVER						
33	3	F	OVER						
34	3	G	OVER						
35	3	H	OVER						
36									
37									
38									
39									
40									
41									
42									
43									
44									
45									
46									
47									
48									
49									
50									
51									
52									
53									
54									
55									
56									
57									
58									
59									
60									
61									
62									
63									
64									
65									
66									
67									
68									
69									
70									
71									
72									
73									
74									
75									
76									
77									
78									
79									
80									
81									
82									
83									
84									
85									
86									
87									
88									
89									
90									
91									
92									
93									
94									
95									
96									
97									
98									
99									
100									

8	A	<i>BLK</i>	0.001	
8	B	<i>✓1</i>	0.435	1.045
8	C	<i>✓4</i>	0.211	2.0383
8	D	<i>✓7</i>	0.051	30.504
8	E	<i>10</i>	0.017	40.211
8	F	<i>13</i>	0.026	5.0192
8	G		0.018	6.0071
8	H		0.015	7.0051
9	A	<i>BLK</i>	0.011	8.0036
9	B	<i>✓2</i>	0.383	
9	C	<i>✓5</i>	0.192	
9	D	<i>18</i>	0.036	
9	E	<i>11</i>	0.005	
9	F	<i>14</i>	0.016	
9	G		0.021	
9	H		0.019	
10	A	<i>BLK</i>	0.022	
10	B	<i>✓3</i>	0.504	
10	C	<i>✓6</i>	0.071	
10	D	<i>±9</i>	0.031	
10	E	<i>12</i>	0.024	
10	F	<i>15</i>	0.018	
10	G		0.040	
10	H		0.033	

Handwritten signature

5/27/82

Can we hybridize in microElioplates

22 ng

1	① Ad2 DNA	2	② Ad2 DNA	3	③ Ad2 DNA	4	④ Ad2	5	⑤ Ad2	6	⑥ Ad2
	② Ad2 DNA		AD ₂ +BSA								

③ Ad2 + SSC carrier

Denat 0.3M NaOH
to 14 Amm Acet
Dry 39°

Block


1, 2, 3
Blocking
buff
+SS

4, 5, 6
1% BSA

200 ul.
37° - 1HR

Blocking buff. | Block BSA

Ad →
AD BSA →
Ad carrier →



Hybridize with
40 ng DNA

[Handwritten signature]
A.H.J.

6/11/82
BY

— Burwood exp —

SC
10/1/82

5/28/82

6/2/82

32
P in mineral
+ Reaflucor

32pm plastic in mineral
+ Reaflucor

32P + POPOP

32P + plastic + POPOP

32p + toluene

32p + plastic + toluene

we can count "plastic" things

[Handwritten signature]
ATG

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.