

- 33 -

Example 28: 40-O-[2-(4',5'-dicarboethoxy-1',2',3'-triazol-1'-yl)-ethyl]-rapamycin

98 mg of 40-O-(2-azidoethyl)-rapamycin and 32 mg diethylacetylene dicarboxylate are suspended in 0.5 ml toluene and heated at 65 C for 5h. The reaction mixture is then cooled at room temperature, loaded on 10 g silica gel and eluted with hexane/ethyl acetate 1/1 to afford the title product: MS (FAB) m/z 1175 (20%, M+Na); 1121 (15%, M-MeOH); 1103 (60%, M-(MeOH+H₂O))

H-NMR (CDCl₃) δ: 0.62 (1H, q, J=12 Hz); 1.40 (3H, t, J=8 Hz); 1.42 (3H, t, J=8 Hz); 3.13 (3H, s); 3.25 (3H, s); 3.33 (3H, s)

MBA (rel. IC₅₀): 2.7

IL-6 dep. prol. (rel. IC₅₀): 12

The previous examples may also be made using as starting material instead of rapamycin, 9-deoxo-rapamycin, 26-dihydro rapamycin, or 9-deoxo-, 26-dihydro-rapamycin. Alternatively, and preferably, as described e.g., in example 20, the rapamycin compounds of the above examples may be hydrogenated or reduced, using suitable protecting groups where necessary. The following novel methods for reducing the keto at C9, or hydrogenating the keto at C26 are provided:

Example 29: Removal of keto at C9

A stream of hydrogen sulfide is passed at room temperature through a stirred solution of 3.2 g (3.5 mmol) of rapamycin in 50 ml pyridine and 2.5 ml DMF. The solution turns from colorless to yellow. After two hours, the introduction of hydrogen sulfide is stopped and stirring is continued for five days, during which time the solution turns gradually orange. TLC and HPLC analysis verifies complete consumption of the starting material and the presence of a single new compound. The solution is purged with nitrogen for one hour and concentrated under reduced pressure. The residue is taken up in ethyl acetate, washed with cold 1N HCl solution (3x), saturated sodium bicarbonate solution and saturated brine. The organic layer is dried over anhydrous sodium sulfate and filtered and concentrated under reduced pressure. The residue is taken up in ether and the precipitated

- 34 -

sulfur is filtered off. Concentration of the ethereal solution followed by column chromatography on silica gel (10:4:1 CH₂Cl₂/i-Pr₂O/MeOH) yields 9-deoxorapamycin as a colorless foam. The identity of the product is confirmed by nuclear magnetic resonance spectroscopy (NMR), mass spectrometry (MS), and/or infrared spectroscopy (IR). 9-deoxorapamycin is found to exhibit the following characteristic physical data: ¹H NMR (CDCl₃) δ 1.61 (3H,d,J = 1 Hz, C17-CH₃), 1.76 (3H,d,J = 1.2 Hz,C29-CH₃), 2.42 (1H,d,J = 14.5 Hz, H-9), 2.74 (1H,d,J = 14.5 Hz, H-9), 3.13 (3H,s,C16-OCH₃) 3.5 (3H,s,C27-OCH₃), 3.40 (3H,s,C39-OCH₃), 5.40 (1H,d,J = 10 Hz, H-30), 5.57 (1H,dd,J₁ = 8.6 Hz, J₂ = 15 Hz, H-22), 5.96 (1H,d,J = 9 Hz, H-18), 6.09 (1H,d,J = 1.7 Hz, 10-OH), 6.15 (1H,dd,J₁ = 10 Hz, J₂ = 15Hz, H-21), 6.37 (1H,dd,J₁ = 1.5 Hz, J₂ = 5 Hz, H-19), 6.38 (1H,J = 9.5 Hz, H-20). ¹³C NMR (CDCl₃) δ 38.5 (C-9), 98.0 (C-10), 170.7 (C-1), 173.0 (C-8), 208.8 (C-32), 216.9 (C-26).

MS(FAB) m/z 922 8[M+Na⁺], 899 (M⁺), 881 ([M-H₂O]⁺), 868 ([M-OCH₃]⁺), 850 ([M-(H₂O+OCH₃)]⁺).

IR (major peaks)(cm⁻¹) 987, 1086, 1193, 1453, 1616, 1717, 1739, 3443.

MBA (rel. IC₅₀): 1

MLR (rel. IC₅₀): 14

IL-6 dep. prol. (rel. IC₅₀): 9

Example 30: Dihydrogenation of keto at C26

To a stirred solution of 421 mg (1.6 mmol) of tetramethylammonium triacetoxymethylborohydride in 2 ml of acetonitrile is added 2 ml of acetic acid. The resulting mixture is stirred for 30 minutes at room temperature and cooled to -35°C. At this temperature a solution of 180 mg (0.2 mmol) of 9-deoxo-rapamycin in 1 ml of acetonitrile is added and the resulting mixture is allowed to stir for 24 hours. The mixture is quenched with a saturated sodium potassium tartrate solution and allowed to warm to room temperature. Stirring is continued until both layers are clear and ethyl acetate is added. The layers are separated and the aqueous layer is extracted twice with ethyl acetate. The resulting organic solution is washed once with a 10% sodium bicarbonate solution and twice with

- 35 -

saturated brine, then dried over anhydrous sodium sulfate, filtered and concentrated under reduced pressure. The residue is purified by column chromatography on silica gel (90:10 AcOEt-hexane). As the starting material in this case was 9-deoxorapamycin, the final compound is 9-deoxorapamycin, 26-dihydrorapamycin is produced as a colorless foam, having the following characteristic spectroscopic data: ^1H NMR (CDCl_3) (major isomer) δ : 0.93 (3H,d,J = 6.9 Hz, CHCH_3), 1.00 (3H,d,J = 6.9 Hz CHCH_3), 1.07 (3H,d,J = 6.9 Hz, CHCH_3), 1.17 (3H,d,J = 6.9 Hz, CHCH_3), 1.61 (3H,d,J = 1Hz, C17- CH_3), 1.73 (3H,d,J = 1.2 Hz, C29- CH_3), 2.43 (1H,dd,J = 4.1 and 16.0 Hz, H-33), 2.46 (1H,d,J = 13.8 Hz, H-9), 2.58 (1H,m,H-25), 2.77 (1H,d,J = 13.8 Hz, H-9), 2.82 (1H,dd,J = 8.3 and 16.0 Hz, H-33), 3.17 (1H,dd,J = 4.1 and 9.2 Hz, H-27), 3.61 (2H,m, H-14 and H28), 5.19 (1H,ddd,J = 4.1, 4.6 and 8.3 Hz, H-34), 5.49 (1H, broad d,J = 5.0 Hz, H-2), 5.56 (1H,d,J = 9.1 Hz, H-30), 5.75 (1H,dd,J = 6.9 and 14.7 Hz, H-22), 5.76 (1H,s,10-OH), 5.99 (1H,broad d,J = 9.2 Hz, H-18), 6.10 (1H,m,H-21), 6.36 (2H,m,H-19 and H-20); MS (FAB) m/z 924 ($[\text{M} + \text{Na}]$), 852 ($[\text{M} - (\text{H}_2\text{O} + \text{CH}_3\text{O})]^+$).
MBA (rel. IC_{50}): 47
MLR (rel. IC_{50}): 134
IL-6 dep. prol. (rel. IC_{50}): 78

26-dihydrorapamycin is prepared in the same manner, using rapamycin in place of 9-deoxorapamycin. This product has the following characteristic spectroscopic data: ^{13}C -NMR (CDCl_3) (major isomer) δ = 208.3 (C-32); 194.0 (C-9); 169.3 (C-1); 166.6 (C-8); 140.9 (C-22); 136.5 (C-29); 136.2 (C-17); 133.5 (C-20); 129.1 (C-21); 128.7 (C-18); 126.2 (C-30); 125.3 (C-19); 98.6 (C-10); 84.4 (C-39); 83.9 (C-16); 81.6 (C-27); 75.4 (C-34); 74.3 (C-28); 73.9 (C-40); 72.9 (C-26); 67.4 (C-14); 59.1 (27- OCH_3); 56.6 (39- OCH_3); 55.9 (16- OCH_3); 51.3 (C-2); 46.8 (C-31); 44.3 (C-6); 40.4 (C-33); 40.4 (C-25); 39.5 (C-24); 38.8 (C-15); 38.0 (C-36); 34.3 (C-23); 34.2 (C-38); 33.5 (C-11); 33.3 (C-37); 33.2 (C-35); 31.5 (C-42); 31.3 (C-41); 30.9 (C-13); 27.1 (C-12); 27.0 (C-3); 25.2 (C-5); 21.4 (23- CH_3); 20.7 (C-4); 17.3 (11- CH_3); 16.1 (31- CH_3); 15.9 (35- CH_3); 14.4 (25- CH_3); 14.2 (29- CH_3); 10.3 (17- CH_3).

- 36 -

MS (FAB) m/z : 884 (M-OCH₃, 35%); 866 (M-[OCH₃ + H₂O], 100%; 848 (M-[OCH₃ + 2 H₂O], 40%).

MBA (rel. IC₅₀): 1.7

MLR (rel. IC₅₀): 1

IL-6 dep. prol. (rel. IC₅₀): 7.5

MULTIPLE
FEE
(FOR

PENDING CLAIM
STATEMENT SHEET
(FORM PTO-375)

08/416623

APPLICANT

AS FILED		AFTER 1st AMENDMENT		AFTER 2nd AMENDMENT		CLAIMS					
IND.	DEP.	IND.	DEP.	IND.	DEP.	IND.	DEP.	IND.	DEP.	IND.	C
1								51			
2								52			
3								53			
4								54			
5								55			
6								56			
7								57			
8								58			
9								59			
10								60			
11								61			
12								62			
13								63			
14								64			
15								65			
16								66			
17								67			
18								68			
19								69			
20								70			
21								71			
22								72			
23								73			
24								74			
25								75			
26								76			
27								77			
28								78			
29								79			
30								80			
31								81			
32								82			
33								83			
34								84			
35								85			
36								86			
37								87			
38								88			
39								89			
40								90			
41								91			
42								92			
43								93			
44								94			
45								95			
46								96			
47								97			
48								98			
49								99			
50								100			
TOTAL IND.								TOTAL IND.			
TOTAL DEP.								TOTAL DEP.			
TOTAL CLAIMS								TOTAL CLAIMS			

PTO-1360 (3-78)

*MAY BE USED FOR ADDITIONAL CLAIMS OR AMENDMENTS

U.S. DEPARTMENT OF COMMERCE
Patent and Trademark Office

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.