NDA 208073

Page 5

HIGHLIGHTS OF PRESCRIBING INFORMATION These highlights do not include all the information needed to use XIIDRA safely and effectively. See full prescribing information for XIIDRA.

 $\mathbf{XIIDRA^{TM}}$ (liftegrast ophthalmic solution) 5%, for topical ophthalmic use

Initial U.S. Approval: 2016

-----INDICATIONS AND USAGE---

Xiidra (liftegrast ophthalmic solution) 5% is a lymphocyte function-associated antigen-1 (LFA-1) antagonist indicated for the treatment of the signs and symptoms of dry eye disease (DED). (1)

 -----DOSAGE FORMS AND STRENGTHS-------Ophthalmic solution containing liftegrast 5% (50 mg/mL). (3)

-----CONTRAINDICATIONS------

None (4)

To report SUSPECTED ADVERSE REACTIONS, contact Shire US Inc. at 1-800-828-2088 or FDA at 1-800-FDA-1088 or www.fda.gov./medwatch.

See 17 for PATIENT COUNSELING INFORMATION and FDA-approved patient labeling.

Revised: 06/2016

- FULL PRESCRIBING INFORMATION: CONTENTS*
- 1 INDICATIONS AND USAGE
- 2 DOSAGE AND ADMINISTRATION3 DOSAGE FORMS AND STRENGTHS
- 4 CONTRAINDICATIONS
- 6 ADVERSE REACTIONS
- 6.1 Clinical Studies Experience
- 8 USE IN SPECIFIC POPULATIONS
 - 8.1 Pregnancy
 - 8.2 Lactation
 - 8.4 Pediatric Use
 - 8.5 Geriatric Use

- 11 DESCRIPTION
- 12 CLINICAL PHARMACOLOGY 12.1 Mechanism of Action
 - 12.3 Pharmacokinetics
- 13 NONCLINICAL TOXICOLOGY
- 13.1 Carcinogenesis, Mutagenesis, Impairment of Fertility14 CLINICAL STUDIES
- 16 HOW SUPPLIED/STORAGE AND HANDLING

17 PATIENT COUNSELING INFORMATION

*Sections or subsections omitted from the full prescribing information are not listed.

FULL PRESCRIBING INFORMATION

1 INDICATIONS AND USAGE

Xiidra[™] (lifitegrast ophthalmic solution) 5% is indicated for the treatment of the signs and symptoms of dry eye disease (DED).

2 DOSAGE AND ADMINISTRATION

Instill one drop of Xiidra twice daily (approximately 12 hours apart) into each eye using a single-use container. Discard the single-use container immediately after using in each eye.

Contact lenses should be removed prior to the administration of Xiidra and may be reinserted 15 minutes following administration.

3 DOSAGE FORMS AND STRENGTHS

Ophthalmic solution containing lifitegrast 5% (50 mg/mL).

4 CONTRAINDICATIONS

None.

6 ADVERSE REACTIONS

6.1 Clinical Studies Experience

Because clinical studies are conducted under widely varying conditions, adverse reaction rates observed in clinical studies of a drug cannot be directly compared to rates in the clinical trials of another drug and may not reflect the rates observed in practice.

In five clinical studies of dry eye disease conducted with liftegrast ophthalmic solution, 1401 patients received at least 1 dose of liftegrast (1287 of which received liftegrast 5%). The majority of patients (84%) had \leq 3 months of treatment exposure. 170 patients were exposed to liftegrast for approximately 12 months. The majority of the treated patients were female (77%). The most common adverse reactions reported in 5-25 % of patients were instillation site irritation, dysgeusia and reduced visual acuity.

Other adverse reactions reported in 1% to 5% of the patients were blurred vision, conjunctival hyperemia, eye irritation, headache, increased lacrimation, eye discharge, eye discomfort, eye pruritus and sinusitis.

8 USE IN SPECIFIC POPULATIONS

8.1 Pregnancy

Risk Summary

There are no available data on Xiidra use in pregnant women to inform any drug associated risks. Intravenous (IV) administration of lifitegrast to pregnant rats, from pre-mating through gestation day 17, did not produce teratogenicity at clinically relevant systemic exposures. Intravenous administration of lifitegrast to pregnant rabbits during organogenesis produced an increased incidence of omphalocele at the lowest dose tested, 3 mg/kg/day (400-fold the human plasma exposure at the recommended human ophthalmic dose [RHOD], based on the area under the curve [AUC] level). Since human systemic exposure to lifitegrast following ocular administration of Xiidra at the RHOD is low, the applicability of animal findings to the risk of Xiidra use in humans during pregnancy is unclear [see Clinical Pharmacology (12.3)].

<u>Data</u>

Animal Data

Lifitegrast administered daily by intravenous (IV) injection to rats, from pre-mating through gestation day 17, caused an increase in mean preimplantation loss and an increased incidence of several minor skeletal anomalies at 30 mg/kg/day, representing 5,400-fold the human plasma exposure at the RHOD of Xiidra, based on AUC. No teratogenicity was observed in the rat at 10 mg/kg/day (460-fold the human plasma exposure at the RHOD, based on AUC). In the rabbit, an increased incidence of omphalocele was observed at the lowest dose tested, 3 mg/kg/day (400-fold the human plasma exposure at the RHOD, based on AUC), when administered by IV injection daily from gestation days 7 through 19. A fetal No Observed Adverse Effect Level (NOAEL) was not identified in the rabbit.

8.2 Lactation

Risk Summary

There are no data on the presence of liftegrast in human milk, the effects on the breastfed infant, or the effects on milk production. However, systemic exposure to liftegrast from ocular administration is low *[see Clinical Pharmacology (12.3)]*. The developmental and health benefits of breastfeeding should be considered, along with the mother's clinical need for Xiidra and any potential adverse effects on the breastfed child from Xiidra.

8.4 Pediatric Use


Safety and efficacy in pediatric patients below the age of 17 years have not been established.

8.5 Geriatric Use

No overall differences in safety or effectiveness have been observed between elderly and younger adult patients.

11 DESCRIPTION

The chemical name for liftegrast is (S)-2-(2-(benzofuran-6-carbonyl)-5,7-dichloro-1,2,3,4-tetrahydroisoquinoline-6-carboxamido)-3-(3-(methylsulfonyl)phenyl)propanoic acid. The molecular formula of liftegrast is $C_{29}H_{24}Cl_2N_2O_7S$ and its molecular weight is 615.5. The structural formula of liftegrast is:



Lifitegrast is a white to off-white powder which is soluble in water.

Xiidra (lifitegrast ophthalmic solution) 5% is a lymphocyte function-associated antigen-1 (LFA-1) antagonist supplied as a sterile, clear, colorless to slightly brownish-yellow colored, isotonic solution of lifitegrast with a pH of 7.0–8.0 and an osmolality range of 200–330 mOsmol/kg.

Xiidra contains Active: lifitegrast 50 mg/mL: Inactives: sodium chloride_sodium phosphate dibasic anhvdrous

Find authenticated court documents without watermarks at docketalarm.com.

12 CLINICAL PHARMACOLOGY

12.1 Mechanism of Action

Lifitegrast binds to the integrin lymphocyte function-associated antigen-1 (LFA-1), a cell surface protein found on leukocytes and blocks the interaction of LFA-1 with its cognate ligand intercellular adhesion molecule-1 (ICAM-1). ICAM-1 may be overexpressed in corneal and conjunctival tissues in dry eye disease. LFA-1/ICAM-1 interaction can contribute to the formation of an immunological synapse resulting in T-cell activation and migration to target tissues. *In vitro* studies demonstrated that lifitegrast may inhibit T-cell adhesion to ICAM-1 in a human T-cell line and may inhibit secretion of inflammatory cytokines in human peripheral blood mononuclear cells. The exact mechanism of action of lifitegrast in dry eye disease is not known.

12.3 Pharmacokinetics

In a subset of dry eye disease patients (n=47) enrolled in a Phase 3 trial, the pre-dose (trough) plasma concentrations of liftegrast were measured after 180 and 360 days of topical ocular dosing (1 drop twice daily) with Xiidra (liftegrast ophthalmic solution) 5%. A total of nine (9) of the 47 patients (19%) had plasma liftegrast trough concentrations above 0.5 ng/mL (the lower limit of assay quantitation). Trough plasma concentrations that could be quantitated ranged from 0.55 ng/mL to 3.74 ng/mL.

13 NONCLINICAL TOXICOLOGY

13.1 Carcinogenesis, Mutagenesis, Impairment of Fertility

Carcinogenesis

Animal studies have not been conducted to determine the carcinogenic potential of lifitegrast.

<u>Mutagenesis</u>

Lifitegrast was not mutagenic in the *in vitro* Ames assay. Lifitegrast was not clastogenic in the *in vivo* mouse micronucleus assay. In an *in vitro* chromosomal aberration assay using mammalian cells (Chinese hamster ovary cells), lifitegrast was positive at the highest concentration tested, without metabolic activation.

Impairment of fertility

DOCKE

Lifitegrast administered at intravenous (IV) doses of up to 30 mg/kg/day (5400-fold the human plasma exposure at the recommended human ophthalmic dose (RHOD) of lifitegrast ophthalmic solution, 5%) had no effect on fertility and reproductive performance in male and female treated rats.

14 CLINICAL STUDIES

The safety and efficacy of lifitegrast for the treatment of dry eye disease were assessed in a total of 1181 patients (1067 of which received lifitegrast 5%) in four 12-week, randomized, multi-center, double-masked, vehicle-controlled studies. Patients were randomized to Xiidra or vehicle (placebo) in a 1:1 ratio and dosed twice a day. Use of artificial tears was not allowed during the studies. The mean age was 59 years (range, 19–97 years). The majority of patients were female (76%). Enrollment criteria included, minimal signs (i.e., Corneal Fluorescein Staining (CFS) and non-anesthetized Schirmer Tear Test (STT)) and symptoms (i.e., Eye Dryness Score (EDS) and Ocular Discomfort Score (ODS)) severity scores at baseline.

Effects on Symptoms of Dry Eye Disease

Eye dryness Score (EDS) was rated by patients using a visual analogue scale (VAS) (0 = no discomfort, 100 = maximal discomfort) at each study visit. The average baseline EDS was between 40 and 70. A larger reduction in EDS favoring Xiidra was observed in all studies at Day 42 and Day 84 (see Figure 1).

Find authenticated court documents without watermarks at docketalarm.com.

Study 1 Visit	Vehicle (N = 58)	Xiidra (N = 58)	Difference ^[1] (95% Cl)	← Favors Xiidra	Study 2 Visit	Vehicle (N = 295)	Xiidra (N = 293)	[1] Difference (95% Cl)	← Favors Xiidra
Baseline	51.8 (23.55)	51.6 (24.69)			Baseline	41.6 (29.69)	40.2 (28.64)		
Day 14	-3.9 (25.46)	-8.9 (21.72)	-5.1 (-13.1, 3.0)	· • •	Day 14	-7.5 (29.01)	-6.7 (27.36)	0.1 (-3.9, 4.1)	
Day 42	-7.9 (19.60)	-17.3 (24.96)	-9.4 (-17.0, -1.9)	·	Day 42	-9.1 (30.03)	-12.6 (30.71)	-4.2 (-8.5, 0.0)	
Day 84	-7.2 (25.29)	-14.4 (25.36)	-7.3 (-16.1, 1.4)		Day 84	-11.2 (28.78)	-15.2 (31.48)	-4.7 (-8.9, -0.4)	
				-20 -10 0 5					-20 -10 0 5
Study 3 Visit	Vehicle (N = 360)	Xiidra (N = 358)	[1] Difference (95% CI)	- Fauera Vildea	Study 4 Visit	Vehicle (N = 356)	Xiidra (N = 355)	[1] Difference (95% CI)	Son Fastara Wildon
Baseline	69.2 (16.76)	69.7 (16.95)		← Favors Xiidra	Baseline	69.0 (17.08)	68.3 (16.88)		← Favors Xiidra
Day 14	-13.1 (24.04)	-19.7 (26.49)	-6.4 (-10.0, -2.8)		Day 14	-14.9 (22.35)	-22.7 (25.41)	-8.0 (-11.4, -4.5)	
Day 42	-18.2 (26.51)	-28.3 (27.69)	-10.0 (-13.8, -6.1)		Day 42	-23.7 (25.98)	-33.0 (27.46)	-9.6 (-13.4, -5.8)	
Day 84	-22.8 (28.60)	-35.3 (28.40)	-12.3 (-16.4, -8.3)		Day 84	-30.5 (28.03)	-37.7 (28.91)	-7.5 (-11.6, -3.5)	· • ·
				-20 -10 0 5					-20 -10 0 5

Figure 1: Mean Change (SD) from Baseline and Treatment Difference (Xiidra – Vehicle) in Eye Dryness Score in 12-Week Studies in Patients with Dry Eye Disease

[1] Based on ANCOVA model adjusted for baseline value in Study 1, and ANCOVA model adjusted for baseline value and randomization stratification factors in Studies 2-4. All randomized and treated patients were included in the analysis and missing data were imputed using last-available data. In Study 1, one Xiidra treated subject who did not have a baseline value was excluded from analysis.

Effects on Signs of Dry Eye Disease

DOCKE

Δ

RM

Inferior fluorescein corneal staining score (ICSS) (0 = no staining, 1 = few/rare punctate lesions, 2 = discrete and countable lesions, 3 = lesions too numerous to count but not coalescent, 4 = coalescent) was recorded at each study visit. The average baseline ICSS was approximately 1.8 in Studies 1 and 2, and 2.4 in Studies 3 and 4. At Day 84, a larger reduction in ICSS favoring Xiidra was observed in three of the four studies (see Figure 2).

Figure 2: Mean Change (SD) from Baseline and Treatment Difference (Xiidra – Vehicle) in Inferior Corneal Staining Score in 12-Week Studies in Patients with Dry Eye Disease.

Study 1 Visit	Vehicle (N = 58)	Xiidra (N = 58)	Difference ^[1] (95% CI)		Study 2	Vehicle (N = 295)	Xiidra (N = 293)	Difference ^[1] (95% CI)	
Baseline	1.65 (0.513)	1.77 (0.515)		← Favors Xiidra	Baseline	1.81 (0.599)	1.84 (0.597)		← Favors Xiidra
Day 14	0.24 (0.709)	0.06 (0.522)	-0.14 (-0.36, 0.08)		Day 14	0.08 (0.771)	0.04 (0.734)	-0.03 (-0.14, 0.08)	
Day 42	0.19 (0.694)	0.08 (0.591)	-0.05 (-0.28, 0.17)	·•	Day 42	-0.02 (0.893)	-0.14 (0.861)	-0.10 (-0.23, 0.02)	∴ ⊷ +
Day 84	0.38 (0.785)	0.04 (0.745)	-0.25 (-0.50, -0.00)	· · · · ·	Day 84	0.17 (0.819)	-0.07 (0.868)	-0.23 (-0.36, -0.10)	<u></u>
				-0.50 0.00 0.25					-0.50 0.00 0.2
Study 3 Visit	Vehicle (N = 360)	Xiidra (N = 358)	[1] Difference (95% CI)		Study 4 Visit	Vehicle (N = 356)	Xiidra (N = 355)	[1] Difference (95% Cl)	
Baseline	2.40 (0.722)	2.39 (0.763)		+ Favors Xiidra	Baseline	2.46 (0.746)	2.46 (0.681)		← Favors Xiidra
Day 14	-0.48 (0.798)	-0.48 (0.802)	-0.00 (-0.11, 0.11)	· • • •	Day 14	-0.44 (0.775)	-0.49 (0.914)	-0.05 (-0.17, 0.07)	
Day 42	-0.60 (0.899)	-0.69 (0.918)	-0.09 (-0.22, 0.04)		Day 42	-0.66 (0.927)	-0.69 (0.941)	-0.03 (-0.16, 0.10)	·•
Day 84	-0.71 (0.943)	-0.73 (0.926)	-0.03 (-0.16, 0.10)	-0.50 0.00 0.25	Day 84	-0.63 (0.911)	-0.80 (0.939)	-0.17 (-0.30, -0.03)	-0.50 0.00 0.2

[1] Based on ANCOVA model adjusted for baseline value in Study 1, and ANCOVA model adjusted for baseline value and randomization stratification factors in Studies 2-4. All randomized and treated patients were included in the analysis and missing data were imputed using last-available data. In Study 2, one Vehicle treated subject who did not have a study eye designated was excluded from analysis.

Find authenticated court documents without watermarks at docketalarm.com.

DOCKET A L A R M



Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.