



Affi

Chromatography

O]

TE CHROMATOGRAPHY IC ACIDS AND PROTEINS

Herbert Schott

Universität Tübingen stitut für Organische Chemie en, Federal Republic of Germany

Exhibit 2004 Page 2

DOCKET A L A R M

MARCEL DEKKER, INC.

New York and Basel

Library of Congress Cataloging in Publication Data

Schott, Herbert, [date]
Affinity chromatography.

(Chromatographic science; v. 27) Bibliography: p. Includes index.

1. Affinity chromatography. 2. Nucleic acids--Separation. 3. Proteins--Separation. I. Title. II. Series. [DNLM: 1. Chromatography, Affinity. 2. Nucleic Acids--isolation & purification. 3. Proteins-isolation & purification. W1 CH943 v.27 / QU 25 S375a] QP519.9.A35S34 1984 574.19'245 84-17461 ISBN 0-8247-7111-7

COPYRIGHT © 1984 by MARCEL DEKKER, INC. ALL RIGHTS RESERVED

Neither this book nor any part may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, microfilming, and recording, or by any information storage and retrieval system, without permission in writing from the publisher.

MARCEL DEKKER, INC. 270 Madison Avenue, New York, New York 10016

Current printing (last digit): 10 9 8 7 6 5 4 3 2 1

PRINTED IN THE UNITED STATES OF AMERICA

Exhibit 2004 Page 3



Contents

0	reface	ii
4	bbreviations	ix
i	INTRODUCTION	1
?	IMMOBILIZATION OF NUCLEOSIDE PHOSPHATES, NUCLEOSIDES, AND NUCLEOBASES	3
	I. Immobilization of Nucleoside PhosphatesII. Immobilization of Nucleosides and Nucleobases	3 12
3	IMMOBILIZATION OF OLIGONUCLEOTIDES	15
	 I. Covalent Binding of the Oligonucleotide Terminals II. Secondary Valence Binding of Polymer-Bound Oligonucleotides 	15 19
ļ	IMMOBILIZATION OF POLYNUCLEOTIDES AND NUCLEIC ACIDS	23
	I. Immobilization via Physical and Mechanical Forces or Irradiation	23
	II. Immobilization Using Chemical Activation III. Special Methods for the Immobilization of RNA	28
	and Polyribonucleotides	41
	Exhibit 2004 Page 4	



Special Methods 41

III. SPECIAL METHODS FOR THE IMMOBILIZATION OF RNA AND POLYRIBONUCLEOTIDES

The efficiency of immobilization of high-molecular-weight yeast RNA on aminated silica carriers has been studied [449]. More than 90% of the RNA added was bound to the carrier in undegraded form. The procedure allows us to link as much as 2.0 mg of RNA per gram of carrier. During nucleic acid hybridization (65°C, 2 × SSC) in the presence of trace amounts of pronase, chipping off of RNA was observed only during the first 6 hr of incubation. More than half of the immobilized poly(A) RNA remained tightly bound to the carrier. The immobilized messenger and ribosomal yeast RNAs were hybridized with total yeast [1251] DNA. Several methods for the immobilization of polyribonucleotides take advantage of the vicinal diol groups present or potentially present in polyribonucleotides.

A. Immobilization of Oxidized Polyribonucleotides

The free terminal 3'-hydroxyl groups of polyribonucleotides can be oxidized in aqueous solution by periodate as reported in Sec. I.B of Chap. 2 for the oxidation of mononucleotides. The resulting aldehydropolyribonucleotides form Schiff bases with the primary amino groups of AE-cellulose (Scheme 16). The unstable Schiff bases can be converted to stable structures by reduction using NaBH₄ [8,103-105]. The binding in these products is thought to arise from the formation of a substituted morpholine ring structure that includes the nitrogen of the AE-cellulose and the five atoms that originally constituted the ribose ring of the terminal nucleoside. The assignment of this linkage is based on the structure of the product obtained from the borohydride reduction of the complex formed between periodate-oxidized AMP and methylamide [106].

The aldehyde-polyribonucleotides may also be added to a hydrazine or a hydrazine-containing support to form a stable hydrazone bond (Scheme 16). For example, a mixture of agar and polyacrylic acid hydrazide has been used in this way to immobilize tRNA [107]. Periodate-oxidized mRNA has been coupled to a linear polymer of acrylic hydrazide entrapped in agar [108]. The reduction with NaBH₄ has also been used to stabilize the linkage in a similar complex obtained from periodate-oxidized tRNA and hydrazinyl-Sepharose [109].

The attachment of RNA by its 3' terminus to agarose has been accomplished [110] using alkaline activation of agarose with CNBr,

Exhibit 2004 Page 5



DOCKET

Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.

