

5/26/82

① will DNA stick to plastic.

1	2	3	4	5	6	7	8	9
Ad 2 DDA - biot. w R7 denatured in				- diluted w/o carrier 0.3M NaOH Neur. to 1M Amm acetate				10ng → 8 serial dilutions w/o carrier
↓ on poly lysine coated plate				↓ on plate straight				
denatured in				0.3M NaOH				
↓				on plate straight				

Goat Ab.
Rabbit Goat - APase
Tetra zolium

poly lysine coated plate - buffer turns
purple
red beads come off
smear spots.

NaOH - good
NaOH + Amm ref. better } no bk ground

Dot blots.
BK ground
in buffer
= poly lysine.

• • • 1.5

• • • 3 0.150

• • • 0.018, 0.004

① unlabeled
paper
may v.
better

↑ w/ the rug
is where
no paper
binned.
AHJ per

5/26/82

Ad 2. NaOH - AmmAcet
in Elisa plates
PNPP
read.
fantastic.

AGENCY LINE 22 206

1 A BLANK
1 B 0.006
1 C 0.012
1 D 0.009
1 E 0.013
1 F 0.006
1 G 0.007
1 H 0.005

NaOH

2 A 0.268
2 B 0.194
2 C 0.182
2 D 0.284
2 E 0.471
2 F 0.301
2 G ② 0.697
2 H X 0.838

NaOH + Amm Acet.

3 A OVER
3 B OVER
3 C OVER
3 D OVER
3 E OVER
3 F OVER

OVER

8 A BLK 0.001
8 B ✓ 1 0.435
8 C ✓ 4 0.211
8 D ✓ 7 0.051
8 E 10 0.017
8 F 13 0.026
8 G 0.018
8 H 0.015

9 A BLK 0.011
9 B ✓ 2 0.383
9 C ✓ 5 0.192
9 D ± 8 0.036
9 E 11 0.005
9 F 14 0.016
9 G 0.021
9 H 0.019

10 A BLK 0.022
10 B ✓ 3 0.504
10 C ✓ 6 0.071
10 D ± 9 0.031
10 E 12 0.024
10 F 15 0.018
10 G 0.040
10 H 0.033

1.047
2.0383
30.504
46.211
5.0192
6.0071
7.0051
8.0036

Handwritten signature/initials

5/27/82


Can we hybridize in microELISA plates

	1	2	3	4	5	6	7	8	9	10
	①	②	③	④	⑤	⑥				
<u>22 ng</u>	① Ad2 DNA	Ad2 DNA	Ad2 DNA	Ad2	Ad2	Ad2				
	② Ad2 DNA	AD ₂ +BSA								
		to 1% BSA								

③ Ad2 + ss carrier
 Denat 0.3M NaOH
 to 1M Amm Acet
 Dye 37°

Block
 1, 2, 3
 Blocking buff
 +SS
 4, 5, 6
 1% BSA

200 ul.
 37° - 1HR



Hybridize with
40 ng DNA

[Handwritten signature]

6/17/82

— Burmese exp. —

5/28/82

6/2/82

	1	2	3	4	5	6	7	8	9
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									
21									
22									
23									
24									
25									
26									
27									
28									
29									
30									
31									

32 p in mineral + Reafluc

32 p in plastic in mineral + Reafluc

32p + POPOP

32p + plastic + POPOP

32p + toluene

32p + plastic + toluene

we can count "plastic" things

Handwritten signature/initials

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.