

Patent Number:

US006125523A

United States Patent [19]

Brown et al.

[11] **Date of Patent:** Oct. 3, 2000 [45]

[54]	STENT C USE	RIMPING TOOL AND METHOD OF
[75]	Inventors:	Daniel G. Brown, Temecula; Stephen A. Morales, Santa Clara, both of Calif.
[73]	Assignee:	Advanced Cardiovascular Systems, Inc., Santa Clara, Calif.
[21]	Appl. No.:	09/196,534
[22]	Filed:	Nov. 20, 1998
[51] [52]		B21D 39/04 ; B23P 11/00 29/516; 29/282; 72/90; 606/1
[58]	Field of So	earch

[56] References Cited

U.S. PATENT DOCUMENTS

696,289	3/1902	Williams .
2,964,088	12/1960	Erath 72/90
4,455,854	6/1984	Ermolovich 72/90
4,468,224	8/1984	Enzmann et al
4,576,142	3/1986	Schiff.
4,644,936	2/1987	Schiff.
4,681,092	7/1987	Cho et al
4,697,573	10/1987	Schiff.
4,901,707	2/1990	Schiff.
4,907,336	3/1990	Gianturco 29/515
5,132,066	7/1992	Charlesworth et al
5,133,732	7/1992	Wiktor 606/1
5,183,085	2/1993	Timmermans 140/89
5,189,786	3/1993	Ishikawa et al 29/283
5,437,083	8/1995	Williams et al 606/1
5,546,646	8/1996	Williams et al 29/516
5,626,604	5/1997	Cottone, Jr 606/108
5,630,830	5/1997	Verbeek 606/108
5,653,691	8/1997	Rupp et al 606/194
5,672,169	9/1997	Verbeek 606/1
5,725,519	3/1998	Penner et al
5,738,674	4/1998	Williams et al 606/1
5,746,764	5/1998	Green et al 606/108
5,759,474	6/1998	Rupp et al
		= =

5,783,227	7/1998	Dunham .	
5,785,715	7/1998	Schatz	. 606/1
5,787,572	8/1998	Toms	29/282
5,810,873	9/1998	Morales	. 606/1
5,836,952	11/1998	Davis et al 6	506/108
5,920,975	7/1999	Morales	29/282

6,125,523

FOREIGN PATENT DOCUMENTS

630 623 A2	12/1994	European Pat. Off	
826 346 A1	3/1998	European Pat. Off	
0 873 731 A1	10/1998	European Pat. Off	
0 938 877 A2	9/1999	European Pat. Off	
464004	8/1928	Germany	72/88
159065	1/1921	United Kingdom .	
WO 98/14120	4/1998	WIPO .	
WO 98/19633	5/1998	WIPO .	

OTHER PUBLICATIONS


The eXTraordinary Stent, C.R. Bard Brochure (undated).

Primary Examiner—Daniel C. Crane Attorney, Agent, or Firm-Fulwider Patton Lee & Utecht, LLP


[57] **ABSTRACT**

A hand held tool for crimping a stent onto a balloon of a catheter is disclosed. The stent crimping tool is operated in one hand by squeezing two plates together while simultaneously displacing the plates linearly to crimp and roll the stent held between the plates. Specifically, the crimping tool includes a base plate with two bosses through which respective pins pass linking the base plate to a compression plate. The compression plate has elongated diameter openings at opposite sides thereof to receive the pins. With the elongated diameter openings, the compression plate can pivot at the pins and translate linearly relative to those pins thus enabling the rolling action during the crimping process. A compression profile pad and a tapered profile pad are attached to the crimping areas of the plates to grip and together apply pressure to the uncrimped stent held therebetween. The pads may include specific contours in order to impart a desired profile to the crimped stent.


23 Claims, 4 Drawing Sheets



















DOCKET A L A R M

Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.

