

Biology Contribution

Preclinical Evaluation of Genexol-PM, a Nanoparticle Formulation of Paclitaxel, as a Novel Radiosensitizer for the Treatment of Non-Small Cell Lung Cancer

Michael E. Werner, PhD,^{*,†} Natalie D. Cummings, BS,^{*,†} Manish Sethi, PhD,^{*,†}
Edina C. Wang, BS,^{*,†} Rohit Sukumar, BS,^{*,†} Dominic T. Moore, MPH,[‡]
and Andrew Z. Wang, MD^{*,†}

**Laboratory of Nano- and Translational Medicine, Department of Radiation Oncology, Lineberger Comprehensive Cancer Center, †Carolina Center for Cancer Nanotechnology Excellence, and ‡Division of Biostatistics and Data Management, Lineberger Comprehensive Cancer Center, University of North Carolina-Chapel Hill, Chapel Hill, North Carolina*

Received Sep 17, 2012, and in revised form Feb 1, 2013. Accepted for publication Feb 6, 2013

Summary

Genexol-PM was evaluated preclinically as a radiosensitizer for chemoradiation therapy of non-small cell lung cancer (NSCLC). Using NSCLC cell lines and mouse xenograft models of NSCLC, we demonstrated that Genexol-PM is an effective radiosensitizer and is more effective than Taxol. We also found that Genexol-PM leads to lower paclitaxel dose in normal lung when compared with Taxol. Our findings support the clinical evaluation

Purpose: A key research objective in radiation oncology is to identify agents that can improve chemoradiation therapy. Nanoparticle (NP) chemotherapeutics possess several properties, such as preferential accumulation in tumors, that are uniquely suited for chemoradiation therapy. To facilitate the clinical translation of NP chemotherapeutics in chemoradiation therapy, we conducted preclinical evaluation of Genexol-PM, the only clinically approved NP chemotherapeutic with a controlled drug release profile, as a radiosensitizer using non-small cell lung cancer (NSCLC) as a model disease.

Methods and Materials: The physical characteristics and drug release profile of Genexol-PM were characterized. Genexol-PM's efficacy as a radiosensitizer was evaluated in vitro using NSCLC cell lines and in vivo using mouse xenograft models of NSCLC. Paclitaxel dose to normal lung and liver after Genexol-PM administration were quantified and compared with that after Taxol administration.

Results: Genexol-PM has a size of 23.91 ± 0.41 nm and surface charge of -8.1 ± 3.1 mV. It releases paclitaxel in a controlled release profile. In vitro evaluation of Genexol-PM as a radiosensitizer showed it is an effective radiosensitizer and is more effective than Taxol, its small molecule counterpart, at the half maximal inhibitory concentration. In vivo study of Genexol-PM as a radiosensitizer demonstrated that it is more effective as a radiosensitizer than Taxol. We also found that Genexol-PM leads to lower paclitaxel exposure to normal lung tissue than Taxol at 6 hours postadministration.

Reprint requests to: Andrew Z. Wang, MD, 101 Manning Dr, Campus Box 7512, Chapel Hill, NC 27599. Tel: (919) 445-5208; E-mail: zawang@med.unc.edu

Authors M.E.W. and N.D.C. contributed equally to this work.

This work was supported by the University Cancer Research Fund from the University of North Carolina, a North Carolina Translational and Clinical Sciences Institute (NC TraCS) pilot grant, and a Lung Cancer Research Foundation Grant. A.Z.W. was supported by National Institutes of Health/National Cancer Institute K12 Career Development Award 5-K12-CA120780-01-05 and National Institutes of Health Center for

Nanotechnology Excellence Grant 1-U54-CA151652-01.

Conflict of interest: A.Z.W. is a consultant for Samyang Biopharmaceuticals, Inc.

Supplementary material for this article can be found at www.redjournal.org.

Acknowledgments—We thank the Microscopy Services Laboratory, Animal Studies Core, and the Analytical Chemistry Core (School of Pharmacy) at the University of North Carolina for their assistance with procedures. We also thank E. Claire Dees for advice on this manuscript.

of Genexol-PM in chemo-radiation therapy for NSCLC.

Conclusions: We have demonstrated that Genexol-PM is more effective than Taxol as a radiosensitizer in the preclinical setting and holds high potential for clinical translation. Our data support the clinical evaluation of Genexol-PM in chemoradiation therapy for NSCLC. © 2013 Elsevier Inc.

Introduction

The concurrent administration of chemotherapy and radiation therapy, also called chemoradiation therapy, is an important treatment paradigm in the curative management of many cancers (1). Chemoradiation therapy has not only consistently shown improved local tumor control but also improves survival when compared with either sequential treatment or sole administration of chemotherapy or radiation therapy in some malignancies (1, 2). However, chemoradiation therapy is not without limitations. The concurrent use of both chemotherapy and radiation therapy has significantly higher toxicities compared with either treatment alone or sequential use (3). Thus, the development of agents and approaches to further improve the therapeutic index of chemoradiation therapy has been a major research objective. A key challenge in this effort has been to selectively deliver chemotherapeutics to tumors while minimizing drug dose to normal tissue. Although traditional drug delivery techniques have failed to address this challenge (4), the development of nanoparticle (NP) formulations of chemotherapeutics offers an unprecedented opportunity. NP therapeutics possess important characteristics that are uniquely suited for chemoradiotherapy. NPs, as macromolecules, preferentially accumulate in tumors through the enhanced permeability and retention (EPR) effect, leading to high intratumoral drug concentrations (5). NPs are also unable to penetrate normal vasculatures and capillaries, thus leading to lower drug doses to normal tissues compared with their small molecule counterparts. Moreover, NP therapeutics have been found to have lower systemic toxicity than small molecules (6). Lastly, NP formulations of chemotherapeutics can release their cargo in a controlled fashion (5). Such prolonged release can increase the synergistic effects of chemotherapy and radiation therapy. Indeed, several preclinical studies have demonstrated that NP chemotherapeutics can improve the therapeutic efficacy of chemoradiotherapy (7, 8). However, there have been no reports evaluating NP therapeutics that are clinically approved or undergoing clinical investigation for use in chemoradiation therapy. Such studies are necessary for clinical translation of NP therapeutics in chemoradiation therapy. To fill this knowledge gap, we conducted a preclinical evaluation of Genexol-PM, the only clinically approved second-generation NP chemotherapeutic with a controlled drug release profile, in chemoradiation therapy using non-small cell lung cancer (NSCLC) as a model disease.

Genexol-PM is a polymeric NP micelle formulation of paclitaxel that has been approved in South Korea for the treatment of breast cancer and NSCLC (9, 10). It is composed of low-molecular-weight amphiphilic diblock copolymer, monomethoxy poly(ethylene glycol)-block-poly(D,L-lactide) (mPEG-PDLLA) and paclitaxel (11). Genexol-PM has shown lower toxicity than Taxol with its maximum tolerated dose identified as 2 to 3 times that of Taxol (9-11). Given that chemoradiation therapy is often used for locally advanced and unresectable NSCLC and paclitaxel has been shown to be an excellent radiosensitizer in NSCLC, we chose to study Genexol-PM in chemoradiation therapy using NSCLC as a model disease. Moreover, chemoradiation therapy for

of patients dying from treatment-related pulmonary toxicities (3, 12). Thus, there is a strong clinical need to improve chemoradiation therapy for NSCLC. In this study, Genexol-PM was compared with Taxol in chemoradiation therapy *in vitro* using NSCLC cell lines and *in vivo* using mouse flank xenograft models of NSCLC. We also quantified the paclitaxel dose delivered to normal lung and liver tissues after Genexol-PM and Taxol administration.

Methods and Materials

Materials

Genexol-PM was provided by Samyang Biopharmaceuticals Corporation (Seoul, Korea) as a gift. Each gram of Genexol-PM contains 115 mg of paclitaxel. Genexol-PM was resuspended in phosphate buffer saline based on clinical administration protocol before use. Taxol 6 mg/mL injectable solution (Teva, Sellersville, PA) was purchased from the University of North Carolina Hospital pharmacy.

Characterization of Genexol-PM

Genexol-PM size (diameter, nm) and surface charge (ζ -potential, mV) were characterized using a ZetaPALS dynamic light scattering detector (Brookhaven Instruments, Holtsville, NY). Transmission electron microscopy (TEM) images of Genexol-PM were obtained at the University of North Carolina (UNC) Microscopy Services Laboratory Core Facility. NPs were stained using a phosphotungstate stain before TEM imaging.

Paclitaxel release characterization

To measure the release profile of paclitaxel from Genexol-PM, 0.1 mL of Genexol-PM solution at a concentration of 1 mg/mL was aliquot equally into Slide-A-Lyzer MINI dialysis microtubes with a molecular weight cut-off of 10 kDa (Pierce, Rockford, IL) and subjected to dialysis as described previously (7). Paclitaxel content was quantitatively analyzed using an Agilent 1100 HPLC (Paolo Alto, CA) equipped with a C18 chromolith flash column (Merck, Darmstadt, Germany). Paclitaxel absorbance was measured by a UV-VIS detector at 227 nm and a retention time of 8.5 min in 0.25 mL/min gradient (from 0:100 to 100:0) of acetonitrile/water.

Cell culture

NSCLC cell lines were acquired from the Lineberger Comprehensive Cancer Center Tissue Culture Facility. A549 cells were cultured in Dulbecco's Modified Eagle Medium (DMEM) (Gibco, Invitrogen, Carlsbad, CA) supplemented with 10% fetal bovine serum (Mediatech, Manassas, VA) and penicillin/streptomycin

fetal bovine serum (FBS), 2 mM L-glutamine, 1.5 g/L Na bicarbonate, 10 mM HEPES (4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid) buffer, 1 mM Na pyruvate, penicillin/streptomycin (Mediatech), 4.5 g/L glucose.

In vitro x-ray irradiation

Radiation therapy was given using a Precision X-RAD 320 (Precision X-Ray, North Branford, CT) operating at 320 kVp and 12.5 mA as described previously (7).

Clonogenic survival assay

Cells were seeded at densities ranging from 100 to 200,000 cells in 4 mL of culture medium in 25-mL flasks 1 day before treatment. Cells were treated with IC_{50} doses of Taxol (1.77 μ g/mL for A549 and 17.7 μ g/mL for H460) or Genexol-PM (15 μ g/mL for A549 and 150 μ g/mL for H460) for 1 hour and washed 3 times with fresh media after incubation. Cells were irradiated with 2 fractions of 0, 1, 2, or 3 Gy separated by 12 hours. Cells were incubated and counted as described previously (7). Data were analyzed using Origin Pro 8.6 software. Polynomial curve fitting was performed using second polynomial order.

In vivo tumor assay

Tumors were established in the upper dorsal region of Nu/Nu mice by injecting 1×10^6 of H460 or A549 cells in a 1:1 RPMI:Matrigel solution. Tumors were incubated for 10 days (H460) or 14 days (A549) to reach approximate tumor volume of 100 mm³ before treatment. Mice (7 per experimental group) were treated intravenously with Taxol (2.5 mg/kg) or Genexol-PM (21.5 mg/kg) and subsequently irradiated 6 hours postinjection. The mice were randomized into experimental groups before treatment. The dose rate at a source-subject distance of 70 cm was 50 cGy/min. Mice were irradiated with 5 daily fractions of 3 Gy. The head and abdomen regions of mice were shielded using 0.5 cm of lead. Only the truncal region containing tumor xenograft was irradiated. Tumor volumes were calculated by measuring 2 perpendicular diameters with a caliper and using the formula of $V = 0.5 \times a \times b^2$ where, a and b are the larger and smaller diameters, respectively. The tumors were measured every 2 days, and the relative percent change in tumor volume was calculated using the relation $100 * (V_i - V_0)/V_0$, where V_i is the volume calculated and V_0 is the initial volume on day 1.

Tumor volume was measured every other day until the tumor reached 3 times the initial volume or 2 cm in the maximum dimension, at which point the animal was euthanized. All animal work was approved and monitored by the UNC Animal Care and Use Committee.

Statistical methods

We used response features analysis to analyze our serial measures of tumor growth (13). This method allows for the comparison of individual growth delay profiles by using summary measures. The particular summary measure used was the AUC (or area under the curve), which was particularly appropriate for the growth-delay

a significant difference in AUCs between chemoradiation therapy with Genexol-PM compared with chemoradiation therapy with Taxol. Statistical analyses were performed using SAS statistical software, version 9.2, from the SAS Institute (Cary, NC).

Quantification of paclitaxel dose in tissue

Mice bearing H460 tumors (5 animals per group) were treated intravenously with equivalent paclitaxel dose of Taxol (2.5 mg/kg) or Genexol-PM (21.5 mg/kg). Mice were euthanized after 6 h or 24 h, and organs were collected and weighed. Organs were homogenized 4:1 (v/w) in chloroform to extract paclitaxel. Homogenate was centrifuged, and the organic layer was collected, evaporated, and reconstituted in 0.1% acetic acid in methanol for analysis via liquid chromatography-mass spectrometry using an Applied Biosystems API 4000 triple quadrupole mass spectrometer with an atmospheric pressure chemical ionization interface.

Results

Characterization of Genexol-PM

Genexol-PM's physical properties were characterized as these properties have not been described in previous publications. Genexol-PM has a size of 23.91 ± 0.41 nm by dynamic light scattering. The NPs are also monodisperse with a polydispersity index of 0.08 ± 0.02 . The NP size and dispersity were confirmed by TEM (Fig. 1A). We found that Genexol-PM has a surface charge of -8.1




Fig. 1. Characterization of Genexol-PM. (A) Transmission electron microscopy image of Genexol-PM depicting a monodisperse population of particles with a narrow size distribution of 23.0 ± 4.5 . (B) Drug-release curve of Genexol-PM. Genexol-PM

± 3.1 mV. We also studied the drug release profile of Genexol-PM. Paclitaxel release from the Genexol-PM demonstrated first-order controlled release kinetics with $\sim 65\%$ drug released at 24 hours and 95% drug release at 48 hours (Fig. 1B).

In vitro radiosensitization with Genexol-PM

The therapeutic efficacy of Genexol-PM as a radiosensitizer in vitro was evaluated in 2 commonly studied NSCLC cell lines H460 (large cell) and A549 (adenocarcinoma). To determine the sensitivity of each cell line to Genexol-PM, a dose-response curve was performed with no radiation therapy for each cell line (Fig. S1). Radiosensitization experiments were conducted using IC_{50} concentrations as well as equivalent paclitaxel concentrations of Genexol-PM and Taxol. Clonogenic survival curves of cells after combined Genexol-PM and fractionated radiation therapy were generated for each cell line. They were compared to that of fractionated radiation therapy alone. We demonstrated that despite the lack of EPR effect in vitro, Genexol-PM is an effective radiosensitizer and is more effective than Taxol in both cell lines at IC_{50} doses (Fig. 2). The sensitizer enhancement ratio (SER) of Taxol at 10% survival in H460 cells is 1.03, whereas the SER of Genexol-PM is 1.12. The SER is greater in A549 cells at 10%; with Taxol at 1.12 and Genexol-PM at 1.23. Also, Genexol-PM was as effective a radiosensitizer as Taxol at equivalent paclitaxel doses (Fig. S2). These results confirmed that Genexol-PM is an effective radiosensitizer.

In vivo radiosensitization with Genexol-PM

To validate the in vitro result and to evaluate the efficacy of Genexol-PM as a radiosensitizer in vivo, we used a murine flank xenograft model of NSCLC. Mice bearing either H460 or A549 cell xenograft tumors were treated with Genexol-PM, Taxol, or saline followed by radiation therapy. After treatment with 5 daily fractions of 3 Gy, the tumor growth delay curves demonstrated that chemoradiation therapy with Genexol-PM leads to significantly longer tumor growth delay compared with chemoradiation therapy with Taxol in H460 cells ($P=.008$; Fig. 3A). Although the difference between the 2 treatments does not appear as significant in A549 cells ($P=.18$; Fig. 3B), the impression remains that chemoradiation therapy with Genexol-PM is more effective. Since tumors that were treated with chemotherapy only (Genexol-PM or Taxol) did not show significant difference in growth delay, the difference between Genexol-PM+XRT and Taxol+XRT is likely due to more effective radiosensitization by Genexol-PM (Fig. S3). Our results showed the high potential of Genexol-PM in improving chemoradiation therapy.

Evaluation of normal tissue drug exposure

Pulmonary toxicity is the main cause of mortality and morbidity from chemoradiation therapy for NSCLC. It is caused by normal lung tissue receiving both chemotherapy and radiation therapy. To determine whether Genexol-PM can lead to lower normal lung paclitaxel exposure, we compared the paclitaxel concentrations in mouse lungs after Genexol-PM and Taxol administration. We found that Genexol-PM leads to significantly lower lung paclitaxel concentration ($P=.01$) acutely after injection (6 hours). This




Fig. 2. Efficacy of Genexol-PM as a radiosensitizer in non-small cell lung cancer cell lines in vitro. Radiation survival curves for (A) H460 and (B) A549 non-small cell lung cancer cells treated with saline, Taxol, or Genexol-PM. Cells were irradiated with 2 fractions of indicated doses 12 hours apart. Significant differences between Taxol and Genexol-PM at a given dose are indicated (* $P<.01$, † $P<.05$). Error bars correspond to standard error of the mean (3 samples per time point). XRT = radiation therapy.

and the concern over hepatotoxicity, we also compared the paclitaxel concentration in mouse liver after Genexol-PM and Taxol administration. Although not statistically significant, Genexol-PM lead to lower paclitaxel exposure to mouse liver than paclitaxel at 6 hours ($P=.30$). Similar to the drug exposure in the lung, there was no difference in liver paclitaxel concentration between Genexol-PM and Taxol at 24 hours.

Discussion




Fig. 3. Efficacy of Genexol-PM as a radiosensitizer in chemo-radiation therapy for non-small cell lung cancer cell lines in vivo. Mice bearing flank tumor xenografts were administered (via tail-vein intravenous injection) saline, Taxol, or Genexol-PM followed by radiation therapy (XRT). Tumors were irradiated in 5 daily fractions of 3 Gy beginning at 6 hours after chemotherapy injection. Changes in tumor volume were measured and the tumor growth delay curves for mice bearing (A) H460 and (B) A549 tumors were generated.

widely used as chemotherapeutics, there has been limited clinical evaluation of these drugs in chemoradiation therapy because of toxicity concerns (15). Recent nanomedicine clinical translation efforts have focused on the development of NP formulation of taxanes, a class of chemotherapeutics that is also used in chemoradiation. Today, there are 2 clinically approved NP taxane formulations, nab-paclitaxel (Abraxane) and Genexol-PM, with several additional therapeutics under clinical investigation (paclitaxel poliglumex [Opaxio], BIND-014 [biologically targeted polymeric NP docetaxel]) (5). Among these, Genexol-PM is the only NP therapeutic that possesses controlled drug release and is approved for clinical use. Furthermore, the only other approved NP taxane, nab-paclitaxel (nanoparticle albumin-bound paclitaxel), dissociates upon administration and cannot take full advantage of the EPR effect (16). Therefore, Genexol-PM is the best NP therapeutic candidate for clinical translation in chemoradiation therapy.

Preclinical studies have shown that the optimal NP physical characteristics for tumor targeting are sizes <100 nm with negative surface charges (17). Genexol-PM's properties (23.0 ± 4.5 nm and -8.1 ± 3.1 mV) suggest that it is excellent for tumor targeting. Another advantage of using Genexol-PM in chemoradiotherapy is its controlled drug release. Although previous clinical investigations have demonstrated that prolonged drug exposure can increase the therapeutic efficacy of chemoradiation therapy (18), most small molecule chemotherapeutics, including taxanes, are too toxic to be administered continuously. NP therapeutics, conversely, can increase drug exposure time in tumors without increasing toxicity, both to local healthy tissue and systemically. In this study, we found Genexol-PM releases paclitaxel in a controlled and slow fashion, with only 40% drug release after 16 hours and >90% of drug release after 48 hours. Such slow drug release should lead to increased synergistic effects between paclitaxel and radiation therapy, which in turn should result in improved therapeutic efficacy.

The efficacy of Taxol as a radiosensitizer has been evaluated previously (19). Our study compared the efficacy of Genexol-PM as a radiosensitizer with that of Taxol. At the IC_{50} dose,




Fig. 4. Quantification of paclitaxel dose in mouse lungs after Genexol-PM administration. Mice were administered equivalent paclitaxel

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.