

BIOPHARMACEUTICS  
AND  
RELEVANT  
PHARMACOKINETICS

# BIOPHARMACEUTICS AND RELEVANT PHARMACOKINETICS

by JOHN G. WAGNER, Phm.B., B.S.P., B.A., Ph.D.

*PROFESSOR OF PHARMACY, COLLEGE OF PHARMACY, AND  
ASSISTANT DIRECTOR FOR RESEARCH AND DEVELOPMENT,  
PHARMACY SERVICE, UNIVERSITY HOSPITAL,  
THE UNIVERSITY OF MICHIGAN, ANN ARBOR, MICHIGAN*

*WITH A SPECIAL CHAPTER ON 'QUALITY CONTROL'*

by M. Pernarowski, B.S.P., M.S., Ph.D.

*PROFESSOR OF PHARMACEUTICAL CHEMISTRY,  
FACULTY OF PHARMACY, UNIVERSITY OF BRITISH COLUMBIA,  
VANCOUVER, BRITISH COLUMBIA, CANADA*

*FIRST EDITION 1971*

*DRUG INTELLIGENCE PUBLICATIONS, HAMILTON, ILLINOIS 62341*

COPYRIGHT © BY DRUG INTELLIGENCE PUBLICATIONS 1971

*All rights reserved.*

*No part of this book may be reproduced  
in any form without written  
permission of the copyright holder*

Library of Congress Catalog Number 75-160736

*Printed in the United States of America by The Hamilton Press, Hamilton, Illinois 62341*

# Enteric Coatings

## *Definition*

AN ENTERIC COATING IS ONE THAT RESISTS THE ACTION OF stomach fluids and disintegrates in the intestines.

## *Reasons for Enteric Coating*

These are: (1) to prevent gastric digestion of a susceptible compound such as protein, or hydrogen ion catalyzed decomposition of a susceptible drug such as erythromycin; (2) to prevent nausea and vomiting caused by a drug, *e.g.*, emetine, atabrine, iron salts and diethylstilbestrol; (3) to prevent dilution of a drug before reaching its site of action in the intestine, *e.g.*, intestinal antiseptics and anthelmintics; (4) to prevent hindrance of gastric digestion of normal foodstuff by a drug, *e.g.*, alkaline medicaments; (5) to provide delayed action of a drug, *e.g.*, barbiturates, amphetamines and aspirin; (6) to deliver medication to the intestinal tract for optimum absorption in the duodenum and jejunum.

## *Requirements of a Good Enteric Coating*

These are: (1) the coating must be nontoxic; (2) the components of the coating and their degradation products, if any, must be physiologically inactive; (3) optimally, the coating should not disintegrate or dissolve in

the stomach during the time that the enteric coated dosage form remains in the stomach; (4) the enteric coated dosage form should disintegrate and release its contained medication as rapidly as possible once the dosage form empties from the stomach into the intestines.

## *Composition of Enteric Coatings*

The patent and non-patent literature on enteric coatings is very extensive and much of it has been reviewed by Wagner (1956 and 1966); hence the reader is referred to the two chapters in *Remington's Practice of Pharmacy* for the review of early literature and references.

In general, there are two types of enteric coatings. First, those which are intended to be digestible by enzymes in the intestinal tract. These include such substances as keratin, formalized gelatin, glycerides, waxes, etc. This type of coating usually does not meet the requirements of a good enteric coating since the surface area of the "digestible substance" exposed to the enzyme is so small that the digestion is much too slow. Second, most current enteric coatings contain an ionizable polyacid and, most frequently, a long-chain polymer with ionizable carboxyl groups which is also a film-forming agent. Usually the polymeric agent is applied to the




Figure 23.1. Coating pans in a Coating Development Laboratory (courtesy of The Upjohn Company)

dosage units dissolved in an organic solvent. The coating solution may also contain other substances, such as a plasticizer to aid in film formation, and a surface active agent. A dusting powder, such as talcum, is added during the coating of the dosage units to prevent the dosage units from sticking together; sometimes powder such as talcum is also suspended in the solution of the polymeric agent in the organic solvent.

#### *Methods of Application of Enteric Coatings*

There are four principal methods of applying enteric coatings to tablets. These are: (1) manual pan coating; (2) programmed pan coating; (3) fluidized-bed coating; and (4) compression coating. The oldest method, and probably the one still most widely used, is manual pan coating. Figure 23.1 shows various sizes of coating pans in a pharmaceutical development laboratory. Figure 23.2 shows a row of production-size coating pans which have a diameter of 36 inches. Larger pans are also used including some giant doughnut-shaped coating pans which handle the volume of a dozen regular commercial-size pans. In programmed pan coating most of the operations normally performed manually are done by machines but

the revolving coating pan is similar to that used in manual pan coating. In fluidized-bed coating the tablets are suspended in a column by rapidly moving air introduced at the bottom of the column and the coating solution is sprayed onto the tablets while they are suspended in the air. In compression coating the components of the enteric coating are mixed in the dry state and applied by compression around the core tablet thus producing essentially a tablet-within-a-tablet.

#### *Theory of Enteric Coatings and Their Properties*

In man the stomach contents are usually in the pH range 1 to 3.5 and pH 1 to 2.5 is the most common range. During the night the human stomach usually has the lowest acidity. During the day the taking of meals is the most frequent cause of changes in acidity of the stomach contents. The high acidities, near pH 1, are usually recovered within a short time after the meals are eaten. In man the duodenal contents are usually

# Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

## Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

## Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

## API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

## LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

## FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

## E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.