

# Elucidation of Phosphatidylcholine Composition in Krill Oil Extracted from *Euphausia superba*

Bjørn Winther · Nils Hoem · Kjetil Berge ·  
Léon Reubsaet

Received: 9 June 2010 / Accepted: 30 August 2010 / Published online: 17 September 2010  
© The Author(s) 2010. This article is published with open access at Springerlink.com

**Abstract** High performance liquid chromatography-electrospray tandem mass spectrometry was used to elucidate the phospholipids in krill oil extracted from *Euphausia superba*, an emerging source for human nutritional supplements. The study was carried out in order to map the species of the choline-containing phospholipid classes: phosphatidylcholine and lyso-phosphatidylcholine. In addition, the prevalent phosphatidylcholine class was quantified and the results compared with prior analysis. The qualification was performed with separation on a reverse phase chromatography column, while the quantification was obtained with class separation on a normal phase chromatography column. An Orbitrap system was used for the detection, and pulsed-Q dissociation fragmentation was utilized for the identification of the species. An asymmetrical exclusion list was applied for detection of phospholipid species of lower concentration, significantly improving the number of species observed. A total of 69 choline-containing phospholipids were detected, whereof 60 phosphatidylcholine substances, among others seven with probable omega-3 fatty acids in both *sn*-1 and *sn*-2. The phosphatidylcholine concentration was estimated to be  $34 \pm 5$  g/100 g oil ( $n = 5$ ). These results confirm the complexity of the phospholipid composition of krill oil, and the presence of long chained, heavily unsaturated fatty acids.

**Keywords** Fish oil · Krill oil · Mass spectrometry · Omega-3 · Phosphatidylcholine · Phospholipid

## Abbreviations

EPA	Eicosapentaenoic acid
DHA	Docosahexaenoic acid
lyso-PtdCho	Lyso-phosphatidylcholine
NPLC	Normal phase liquid chromatography
PtdCho	Phosphatidylcholine
PtdEtn	Phosphatidylethanolamine
PtdIns	Phosphatidylinositol
PtdSer	Phosphatidylserine
PL	Phospholipid
RPLC	Reverse phase liquid chromatography

## Introduction

Krill oil has emerged as an important source of omega-3 fatty acids for human consumption during the last decade, and the amount sold on the world market is rapidly increasing. In contrast to traditional omega-3 supplements on today's market, which are based on omega-3 fatty acids bound to triglycerides (such as cod liver oil and fish oil) or bound as ethyl esters (Omacor/Lovaza), krill oil contains a high proportion of omega-3 fatty acids bound to phospholipids.

Krill oil has been investigated in several preclinical and clinical studies [1–4], and there is growing evidence that the molecular form of the omega-3 fatty acids (i.e. triglycerides, ethyl-esters, phospholipids) might be of importance for their biological effect as well as distribution of the omega-3 fatty acids in the body. In one animal study,

B. Winther · L. Reubsaet (✉)  
Department of Pharmaceutical Chemistry,  
School of Pharmacy, University of Oslo, Oslo, Norway  
e-mail: j.l.reubsaet@farmasi.uio.no

N. Hoem · K. Berge  
Aker BioMarine ASA, Fjordalléen 16, Vika,  
P.O. Box 1423, 0115 Oslo, Norway

it was demonstrated that when krill oil and fish oil were administered to Zucker rats with an equimolar dose eicosapentaenoic acid (EPA) + docosahexaenoic acid (DHA), krill oil had stronger and in some instances different effects than fish oil on specific parameters related to the metabolic syndrome [1]. The lipid level in both heart and liver was significantly lower in rats treated with krill oil, when compared to rats fed the fish oil diet. The authors suggest that this difference may be linked to differences in the incorporation of omega-3 fatty acids into membranes, and consequently a reduction of inflammatory molecules and endocannabinoids, which might be relevant for the differences observed between fish oil and krill oil. Further, in the same study, it was demonstrated that the level of DHA in the brain increased significantly after krill oil administration, but not after fish oil administration, when compared to control animals [2]. Thus, omega-3 fatty acids linked to phospholipids may be differently distributed in the body compared to omega-3 fatty acids in other molecular forms. Moreover, in a clinical safety study, the presence of EPA and DHA in the blood plasma was determined after daily administration of 2 g krill oil or 2 g menhaden oil for 4 weeks [3]. The authors concluded that EPA and DHA from krill oil are absorbed at least as well as that from menhaden oil.

The aim of the current study was to characterize the phospholipids in krill oil in more detail to evaluate the composition of the fatty acids present in the phospholipids. The composition was determined using LC/ESI-MS(/MS), a technique which has lately played an important role in characterization of the lipidome in tissues and organisms [5]. An inherent limitation in the use of ESI for the ionization of long chained fatty acids has been described by Koivusalo et al. [6]. The study showed that the instrument response is affected by the acyl chain length. This is a consideration which is important particularly in the quantification of the lipids.

The elucidation of the phospholipid species is often performed either by doing a precursor ion scan or a neutral loss scan with triple quadrupole instrumentation [7–9], or with MS<sup>n</sup> fragmentation with systems based on ion traps [10–13]. Normal phase liquid chromatography (NPLC) and reverse phase liquid chromatography (RPLC) are both frequently used for the separation of the components [10, 14–17]. Of these two separation techniques, RPLC has been shown to be more suitable for species separation and characterization [8].

Different ionization and fragmentation techniques can be used for the evaluation of phospholipids. Ionization of the phospholipids may be performed in negative- and positive-ionization mode. In general, fragmentation of phospholipids in the positive mode provides information about the phospholipid head group, while fragmentation in

the negative mode is the source of structural information. For phospholipids containing choline-headgroups, the choline-specific fragment  $m/z$  184 has been used in precursor ion scanning operating in the positive ionization mode for class determination [7, 18]. Also in the negative mode, class-specific fragments may be used in the characterization. All phospholipid classes, except those containing choline, yield molecular ions  $[M-H]^-$  when a formate-based mobile phase is used. On the other hand, the choline-containing classes form stable adducts with formic acid in the mobile phase, yielding  $[M + FA-H]^-$  ions ( $m/z = M + 45$ ) [19, 20]. With fragmentation, this adduct dissociates with the loss of  $(HCOO + CH_3)$  into the fragment ion  $[M-CH_3]^-$ . This is in particularly useful in methods utilizing RPLC for separation. Although the chromatographic class information is lost in such setups, the class-specific fragments may be used in the characterization of the species [11].

Two ion activation techniques may be used for MS analysis utilizing ion traps: collision-induced dissociation (CID) and pulsed-Q dissociation (PQD) techniques. While CID has a low mass cut off below 28% of the  $m/z$  for the precursor ion, the novel PQD technique eliminates the potential loss of low mass fragments [21, 22]. This difference could be crucial in the fragmentation of larger molecules into low mass, specific fragments, as shown with detection of iTRAQ fragments with a linear ion trap [23].

The fatty acid composition of phosphatidylcholine (PtdCho) from krill oil has previously been investigated by Le Grandois et al. [24]. This study was performed with a method based on the ESI operated in the positive mode with triple quadrupole detection of lithium adduct ions, and showed the presence of a higher number of PtdCho species with long chained unsaturated fatty acids, than seen in egg yolk, ox liver and soy.

We believe the current study verifies previously presented findings and offer new insights into the composition of krill oil. In addition; it shows the advantage of performing an additional fragmentation using an exclusion list in the identification of low prevalent species.

## Experimental Procedures

### Chemicals

Phospholipid standards of lyso-phosphatidylcholine (lyso-PtdCho), PtdCho, phosphatidylethanolamine (PtdEtn), phosphatidylinositol (PtdIns) and phosphatidylserine (PtdSer) were purchased from Sigma-Aldrich (St. Louis, MO, USA). Lyso-PtdCho, PtdCho and PtdEtn were lyophilized powders obtained from egg yolk, whereas the PtdIns source was *glycine max* and the PtdSer source was

bovine brain. EPAX 6000 TG<sup>®</sup> fish oil was donated by EPAX (Ålesund, Norway), and Superba<sup>™</sup> krill oil was obtained from Aker BioMarine (Oslo, Norway). All other chemicals were of MS grade.

#### Instrumentation

The chromatography was carried out on a Dionex system consisting of an Ultimate 3000 pump, an Ultimate 3000 RS autosampler, and an Ultimate 3000 flow manager. Detection was obtained using a linear ion trap LTQ XL coupled to an Orbitrap Discovery, LC-operation, data acquisition and processing were carried out using Chromelion SDK 6.80 SP2 Build 2327 and Xcalibur version 2.0.7 coupled with DCMS<sup>Link</sup> 2.5 (all Instrument-Teknikk AS, Østerås, Norway).

#### Mass Spectrometry

The LTQ Orbitrap system was operated with a spray voltage of 5.00 kV, nitrogen as the sheath gas with flow rate set to 30 arbitrary units, and helium as the collision gas. The quantification of the PtdCho class was performed with a scan from  $m/z = 400$  to  $m/z = 1,000$  operated in negative ionization mode. MS<sup>n</sup> experiments for identification of the choline-containing phospholipids were performed using data dependent PQD for the first fragmentation step. The molecular ion selected for each fragmentation in this step, was the most intense ion detected by the Orbitrap analyzer with target mass resolution of 30,000 and a scan window from  $m/z = 400$  to  $m/z = 1,000$ . The normalized collision energy was 200 and the isolation width 2.00 Da. Subsequently, the most intense fragment ion detected was further fragmented using CID, with the normalized collision energy at 35 and an isolation width of 2.00 Da (MS<sup>3</sup>). The LTQ was utilized for the detection of the fragments and the  $m/z$  range was relative to the  $m/z$  of the molecular ion. An alternative method was used in order to be able to observe species that were not selected for fragmentation in this way. The overall setup of this method was as described above, with the distinction of adding an asymmetric exclusion list. The exclusion list was generated with the purpose of the LTQ to ignore already identified substances. The list was based on the  $m/z$  of the molecular ions, with an exclusion window from this mass-to-charge ratio, up to  $m/z + 1$ . The width of the exclusion window was selected in order to diminish the detection of isotopes of the molecular ions.

#### Chromatographic Conditions RPLC

Chromatographic separation was performed on a ZORBAX Eclipse Plus C18 column with particle diameter of 5  $\mu\text{m}$  and the column dimensions were 150  $\times$  2.1 mm i.d.

The mobile phase A consisted of 90 parts 1% TEA and 0.2% formic acid in water, and 10 parts mobile phase B (v/v). Mobile phase B consisted of 1% TEA and 0.2% formic acid in 60 parts methanol and 40 parts acetonitrile (v/v).

A linear gradient was used for the separation. The system was first kept isocratic at 65% mobile phase B for 5 min after injection of sample. The gradient was then run from 65 to 100% mobile phase B in 5 min and was kept isocratic at 100% mobile phase B for 20 min, before it was returned to the initial condition in 0.1 min. The column was regenerated with 65% mobile phase B for 16 min. The mobile phase flow was set to 0.2 mL/min and the injection volume was 20  $\mu\text{L}$  throughout the study.

#### Chromatographic Conditions NPLC

NPLC was performed on a HiCHROM LiChrospher 100 DIOL column with a particle diameter of 5  $\mu\text{m}$  and column dimensions of 250  $\times$  2 mm i.d. Mobile phase C was 100% chloroform, and the mobile phase D consisted of 0.05% TEA, 0.05% ammonia and 0.1% formic acid in methanol (v/v). For the class separation of the phospholipids, a linear gradient was used. The gradient was run from 5 to 27.5% mobile phase D in 15 min, followed by a rise to 80% in 2 min to flush the column. This concentration was kept isocratic for 4 min, before it was returned to the initial condition in 2 min. The column was regenerated with 5% mobile phase D for 12 min. The mobile phase flow was set to 0.3 mL/min, and the injection volume was 20  $\mu\text{L}$  throughout the study.

#### Sample Preparation

Samples of krill oil and stock solutions of standards were prepared by dissolving the lipids in a mixture of chloroform and methanol at a ratio of 2:1. These solutions were stored at  $-32^\circ\text{C}$  and excessive heating cycles were avoided. Samples were prepared by further dilution with solvents compatible with the mobile phases used. For NPLC, this was achieved with chloroform:MeOH 95:5, while it was attained by dilution in mobile phase A for RPLC.

#### Calibration Curve

For the quantification, a calibration curve was established with samples of PL free fish oil (EPAX<sup>®</sup>) spiked with a PtdCho standard purified from egg yolk to concentrations of 100  $\mu\text{g}/\text{mL}$ . The spiking of PL free fish oil was performed in order to produce comparable matrixes in the standards and the krill oil samples. Stock solutions were made by dissolving PtdCho standard, PL free fish oil and krill oil separately in mixtures of chloroform:MeOH 2:1. The concentration of krill oil and PtdCho standard was

1 mg/mL and for PL free fish oil 10 mg/mL for these solutions. Respectively, 100  $\mu\text{L}$  of PtdCho standards and of krill oil was added to 900  $\mu\text{L}$  of the PL free fish oil, producing samples with concentrations of 100  $\mu\text{g}/\text{mL}$ . For the calibration curve, the PtdCho standard samples were consecutively diluted to the desired concentrations of 10.0, 5.00, 2.50, 1.00, 0.50, 0.25, and 0.10  $\mu\text{g}/\text{mL}$  ( $n = 5$ ) with a mixture of chloroform:MeOH 95:5. The krill oil samples were diluted in the same way to a concentration of 1.00  $\mu\text{g}/\text{mL}$  in order to measure the PtdCho content within the linear area of response of the calibration curve.

## Results

### Selection of MS-Mode for PtdCho-Classification


Initial experiments with standards of PtdEtn, PtdIns, PtdSer, PtdCho and lyso-PtdCho, were performed in both

positive- and negative- ion mode using the LTQ Orbitrap as a detector. Separation of these substances on a NPLC column yielded class-specific baseline separation (data not shown). The results of these tests indicated a minor difference in the signal intensities obtained between the two settings, with a slightly higher response in negative mode.

### Identification of Choline-Containing Phospholipids in Krill Oil

Since krill oil, as established in Fig. 1, is predominantly composed of choline-containing phospholipids, the emphasis of the work was focused on the elucidation of the species in the PtdCho and lyso-PtdCho classes. In NPLC, class separation of the phospholipids is achieved. A clear tendency of the elution order from the column was seen from high  $m/z$  to lower  $m/z$ . PtdCho class species eluted from 6.5 to 8 min, and the elution of the lyso-PtdCho class occurred between 10.5 and 12 min. Some species separation

**Fig. 1** Base peak chromatogram and three dimensional (3D) map of phospholipid class separation of krill oil, performed with a linear gradient NPLC/ESI-MS on a HiCHROM LiChrospher 100 DIOL column (250  $\times$  2 mm i.d., 5  $\mu\text{m}$ ). MS was operated in the negative ionization mode and set to scan  $m/z$  400–1,000. Only the relevant part of the chromatogram is shown (5–15 min). For the major part of the species, the map show adduct ions in the form of  $[M + \text{FA}-\text{H}]^-$



was seen within the PtdCho class, however, this was not adequate for identification of the diverse species within the classes.

As the identification of the components is performed with data dependent fragmentation, chromatographic separations of the substances are critical for detection of the less prominent species. Hence, the separation for characterization of the species was performed utilizing a RPLC C18 column as described under “[Experimental Procedures](#)”. This improved the chromatographic performance for species separation compared to NPLC (Fig. 2). As lyso-PtdCho only carry one fatty acyl group, these components elute earlier in the chromatogram than the PtdCho species. Lyso-PtdCho dominate the region between 12.0 min and 14.5 min in the chromatogram, while the heavier PtdCho components dominate the chromatogram after 14.5 min.

The identification of the species was performed utilizing a MS<sup>3</sup> data dependent fragmentation method, with an initial PQD fragmentation resulting in the loss of methyl formate, followed by a CID fragmentation. Analysis of the fragmentation spectra obtained typically revealed the identity of the substances without ambiguity. However, co-elution of isobaric compounds could potentially complicate the interpretation of the spectra. This challenge is minimized by applying a set of criteria for reliable identification. The following criteria were applied for reliable interpretation of a choline-containing phospholipid: based on the mass-to-charge ratio of the molecular ion, it is likely to be a choline-containing phospholipid (i.e.  $m/z$  being an even number). Following the first CID with PQD, a daughter ion should be produced by the loss of (HCOO + CH<sub>3</sub>) as 60 Da. Further fragmentation with CID of the resulting product should produce specific fragments revealing the nature of the fatty acyl groups in both the *sn*-1 and the *sn*-2 position, either by the occurrence of the fragment for the fatty acyl group itself, or by the presence of the corresponding fragments of the lyso-compound. The sum of the fatty acyl groups elucidated in this matter should yield a mass matching the initial molecular mass. This is illustrated in Fig. 3, showing the elucidated fragment identity for the fragmentation of the 20:5–22:6 diacyl PtdCho.

The spectra were generally dominated by fatty acyl fragments from both the *sn*-1 and the *sn*-2 positions in addition to their corresponding fragments of the lyso-compound, ensuring identification of the species. The lyso-PtdCho and PtdCho substances identified by applying this method are presented in Tables 1 and 2. The relative intensity of the molecular ions is also presented.

As described earlier, the use of signal intensities in MS, for concentration comparison of the different substances, is only semi-quantitative. However, it provides a valuable indication of the composition of the PtdCho and lyso-PtdChoclasses. Chromatograms and fragmentation patterns

are presented in Figs. 4 and 5 for the 10 foremost substances characterized from the PtdCho class.

As data dependent fragmentation methods are, by nature, biased in the selection of the most prevalent substances, the experiments were repeated with the use of an asymmetric exclusion list added to the MS-method. The exclusion list was based on the data sets obtained with the initial settings (i.e. Tables 1, 2). This method allowed the detection and identification of the additional substances presented in Table 3.

#### Quantification of PtdCho-Class in Krill Oil

Krill oil predominantly contains phospholipids from the PtdCho class (Fig. 1). It was therefore attempted to quantify the absolute concentration of this class by use of class separation with NPLC. Quantification of the PtdCho class was performed with a method developed “in-house” with a LTQ Orbitrap mass spectrometer for the detection.

In the construction of the calibration curve, PtdCho concentrations above 1.00 µg/mL resulted in a relative decrease in the MS signal response, yielding a quadratic polynomial curve ( $y = ax^2 + bx + c$ ) where  $a = -9,766$ ,  $b = 330,360$  and  $c = 4763.8$  with  $r^2 = 0.9995$ . From 0.10 µg/mL to 1.00 µg/mL, the calibration curve showed a high degree of linearity ( $r^2 = 0.9995$ ) with a linear regression curve ( $y = bx + c$ ) where  $b = 368,737$  and  $c = -16,547$ .

The latter area was chosen for quantification purpose. From this, the PtdCho content of the undiluted krill oil was determined to be  $34 \pm 5\%$  ( $w/w$ ) ( $n = 5$ ). Comparisons of the mass spectrum of the sample with the spectrum of the PtdCho standard indicated a difference in PtdCho class composition (Fig. 6). From the results, the average acyl chain lengths appear to be higher in krill oil than in egg yolk. This has also previously been shown by others [24]. As mentioned above, instrument response is affected by the acyl chain length of the PtdCho. These differences in chain lengths could therefore influence the quantification of the PtdCho class as discussed later. The quantitative results were compared with an earlier analysis of the krill oil, performed by the accredited analytical company Nofima (Bergen, Norway). They reported the PtdCho concentration in the krill oil sample to be 35 g/100 g oil.

#### Discussion

The fact that fatty acyl chain lengths of the PtdCho species are relatively long, affects both the choice of fragmentation technique and the effect of standards used for quantification purposes. Due to the low mass cut-off limit at 28% of the molecular ion mass with CID fragmentation in ion

# Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

## Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

## Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

## API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

## LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

## FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

## E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.