


applications.

January 24, 2017 in Surface

Treatment (/primary-  
keywords/surface-treatment-0)

[f \(/#facebook\)](#) [t \(/#twitter\)](#) [in \(/#linkedin\)](#) [e \(/#email\)](#)



**An expert offers a primer on popular medical device coatings and their applications.**

*Kevin Guenther*

*A catheter with a lubricious coating.*

In the world of low friction, lubricious coatings for medical devices are diverse and complex. There are many different types of coatings, each with their own advantages, disadvantages, and applications. This is a basic overview of the most popular medical device coatings and what can be expected from each type.

Lubricious coatings are most commonly distinguished by their degree of lubricity, or the amount of reduction in friction they provide. To understand what these measurements mean, we must first learn how these measurements are obtained. In the medical device industry, friction is measured by what is known as a pinch test. In a pinch test, a coated device is secured on an instrument that pulls the device between the jaws of a clamp which produces a load on the part. While clamped, the coated

the higher the CoF and the lower the lubricity of the surface. Moreover, the higher the pull force, the higher the CoF and the higher the lubricity.

When measuring friction between two surfaces, another characteristic is that of static friction pull force vs. dynamic friction pull force. Static friction is the force needed to break the initial friction force so that the two surfaces can move relative to each other. Dynamic friction pull force is that force needed to maintain the motion once it has been initiated. Static friction (sometimes referred to as "stiction") is almost always significantly higher than dynamic friction.

Now that we have an understanding of how the CoF is obtained and the significance of the different amounts of grams of pull force required, we can take a look at some commonly used coatings and compare them.

Before we look at the different types of coatings, let's begin by looking at the properties of an uncoated sample to provide a reference point. A very common material used in medical devices is PEBAX® (a registered trademark of Arkema, Inc.), the brand name of polyether block amide. When uncoated, PEBAX® (hardness of 55 ShoreD) generally produces pinch test results of about 500 to 600 grams of pull force. Assuming a 500-gram clamp force was used, the CoF would be (1.0-1.2), a very low level of lubricity. Other uncoated catheter materials of construction such as nylon 12, polyethylene, and others exhibit similar CoF. Such a high level of friction would make many of today's catheter-delivered minimally invasive procedures nearly infeasible.

Device manufacturers have turned to a variety of coatings to provide the surface lubricity not offered by the bare catheter material. The first class of coating is **silicone oil (siloxane)**. This type of coating has long been used to reduce surface friction on medical catheters, introducers, etc. It is generally applied by the manufacturer, and is relatively inexpensive. It offers a much reduced CoF, at about 0.13, and approximately 125-175 grams of pull force. However, devices coated with silicone oil can be difficult to apply as they are very mobile and can spread to other manufacturing areas where they are undesirable. For example, the presence of silicone oil on a surface can prevent other surface-enhancing coatings from adhering to the device surface. Once silicone oil is present, it can be very difficult to clean from a device or work surface. As a result, these coatings can be a major contaminant risk to the manufacturer.

The second class of coating is **PTFE (Polytetrafluoroethylene)**, known best by the brand name Teflon®, another popularly used material to reduce friction. It has been employed as a material of construction for catheter-type products or as an added coating for guidewires and other devices made from metal. PTFE typically exhibits between 150 and 200 grams of pull force, or a CoF of about 0.3 - 0.4. On bare metals such as stainless steel and nitinol, PTFE offers performance similar to that of silicone oil. The main advantage of PTFE coatings is that they provide some level of lubricity while dry and do not require any source of moisture to perform.



*A friction test system to measure lubricity.*

characteristics as long as the coating remains wet. Compared to PTFE coatings, which are hydrophobic and repel water, hydrophilic coatings provide superior lubricity (lower CoF) when wet. Devices coated with these coatings demonstrate pull forces as low as 3 grams and CoF values as low as 0.01. In clinical use, devices coated with hydrophilic material demonstrate superior control, exhibit lower stiction (static friction) behavior, and require far less force to perform any given movement. Ranging from 3-25 grams of pull force, hydrophilic coatings offer the best CoF, averaging between .01 and .05. This high level of lubricity enables devices to navigate tortuous anatomical pathways while reducing tissue irritation and adding to patient comfort.

Some types of hydrophilic coatings may generate particulates from their surface in use. For these applications where particle generation from a coating is critical, such as cardiovascular and neurovascular, the selected formulation must be one designed specifically for low particle generation. The only real limitation of this category of coatings is that they must be hydrated to perform.

Commercially available hydrophilic coatings average between 8 and 12 grams of pull force that reduces friction over an uncoated surface by as much as 98%. This type of low-friction performance increases the device's ability to navigate through tortuous anatomical pathways, improves device control, reduces tissue damage, and adds to patient comfort.


Advances in surface enhancements, especially low-friction hydrophilic coatings, have accelerated minimally invasive surgery's growth. These coatings now impact a variety of surgical specialties, including cardiovascular, orthopedic, urological, peripheral, and neurological procedures.

Kevin Guenther is the methods director at [Harland Medical Systems](http://directory.qmed.com/harland-medical-systems-inc-comp237833.html) (<http://directory.qmed.com/harland-medical-systems-inc-comp237833.html>), a medical device coating technologies company based in Eden Prairie, MN. He can be reached at [info@harlandmedical.com](mailto:info@harlandmedical.com) (<mailto:info@harlandmedical.com>).

[Images courtesy of [HARLAND MEDICAL SYSTEMS](http://www.harlandmedical.com/) (<http://www.harlandmedical.com/>)]

Filed Under:  
Surface Treatment (/primary-keywords/surface-treatment-0)

## RECOMMENDED FOR YOU



(/supplier-stories-week-july-7)

[\(/how-dry-lubricants-improve-medical-device-assembly\)](#)

**How Dry Lubricants Improve Medical Device Assembly** [\(/how-dry-lubricants-improve-medical-device-assembly\)](#)

Recent medical device design incorporates single-use devices that make...

**Nu-Med Files Patent for a Reactive Tube Coating** [\(/nu-med-files-patent-reactive-tube-coating\)](#)

Nu-Med Plus, a Salt Lake City, UT-based medical device company, recently filed...

**(/could-protect-implantable-devices-dangerous-bacteria) Could This Protect Implantable Devices from Dangerous Bacteria?** [\(/could-protect-implantable-devices-dangerous-bacteria\)](#)

Bacteria may need to find a new home inside the body, as researchers from...

**Using Laser Marking to Meet UDI Requirements** [\(/using-laser-marking-meet-udi-requirements\)](#)

Laser marking is an effective and safe solution for marking unique device...

[See All in Surface Treatment » \(/taxonomy/term/109\)](#)

Comment \*

Post

500 characters remaining

Get Medtech News in Your



Inbox Daily

(subscribe)

Subscribe Now

**TRENDING**

FDA strikes hard at coronavirus outbreak new policy)

Thermo Fisher Scientific Set to Acquire Qiagen for \$11.5B (/thermo-fisher-scientific-set-acquire-qiagen-115b)

The Anatomy of a Medical Device Startup (/anatomy-medical-device-startup)

SweynTooth Cybersecurity Vulnerabilities Put Dozens of Medical Devices at Risk (/sweyntooth-cybersecurity-vulnerabilities-put-dozens-medical-devices-risk)

Livongo is Still Blowing Minds with Its Success (/livongo-still-blowing-minds-its-success)

Why NuVasive's Pulse Delay Is Actually a Good Thing (/why-nuvasives-pulse-delay-actually-good-thing)

Search for products & manufacturers:

Search MDDI/Qmed



(<https://directory.qmed.com/index.html>)

### OUR NEXT EVENT



Boston Convention & Exhibition Center,

Boston

May 06, 2020, to May 07, 2020

([https://biomedboston.com/?\\_mc=arti\\_x\\_mdir\\_edt\\_aud\\_allen\\_bmner\\_med\\_902\\_x-eventspage](https://biomedboston.com/?_mc=arti_x_mdir_edt_aud_allen_bmner_med_902_x-eventspage))



# Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

## Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

## Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

## API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

## LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

## FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

## E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.