



**Table of Contents**

	<i>page</i>	<i>date</i>
<b>601: General Information</b>		
601 A: Principles	601-1	1/94
601 B: Modes of Operation	601-2	1/94
Liquid-Solid Chromatography	601-2	1/94
Liquid-Liquid Chromatography	601-3	1/94
Bonded Phase Chromatography	601-3	1/94
Ion Exchange Chromatography	601-4	1/94
Size Exclusion Chromatography	601-4	1/94
601 C: Instrumentation and Apparatus	601-5	1/94
Basic Components	601-5	1/94
HPLC System Plumbing	601-7	1/94
601 D: Solvents and Reagents	601-10	1/94
Potential Problems	601-10	1/94
Specific Solvents	601-12	1/94

	<i>page</i>	<i>date</i>
Reagent Blanks	601-14	1/94
Safety Precautions	601-14	1/94
601 E: Sample Preparation	601-14	1/94
Sample Cleanup	601-14	1/94
Sample Filtration	601-14	1/94
Sample Solvent Degassing	601-15	1/94
Choice of Sample Solvent	601-15	1/94
601 F: Reference Standards	601-15	1/94
Stock Solutions	601-15	1/94
Working Standard Solutions	601-16	1/94
Storage	601-16	1/94
References	601-16	1/94
<b>602: Columns</b>		
602 A: Column Selection	602-1	1/94
602 B: Analytical Columns	602-1	1/94
Liquid-Solid Chromatography	602-2	1/94
Bonded Phases	602-2	1/94
Ion Exchange	602-3	1/94
Ion Pair	602-4	1/94
Size Exclusion	602-4	1/94
602 C: Column Evaluation	602-5	1/94
602 D: Column Specifications	602-6	1/94
602 E: Analytical Column Protection	602-8	1/94
Filters	602-8	1/94
Precolumns	602-9	1/94
Guard Columns	602-9	1/94
602 F: Column Maintenance and Troubleshooting	602-9	1/94
Column Care	602-9	1/94
Column Evaluation by Injection of Test Mixtures	602-10	1/94
Column Storage	602-11	1/94
Column Regeneration	602-11	1/94
References	602-12	1/94
<b>603: Mobile Phase Selection, Preparation, and Delivery</b>		
603 A: Mobile Phase Selection	603-1	1/94
Normal Phase Chromatography	603-1	1/94
Reverse Phase Chromatography	603-3	1/94
Ion Exchange Chromatography	603-4	1/94
Ion Pair Chromatography	603-5	1/94
Size Exclusion Chromatography	603-5	1/94
Gradient Elution in HPLC	603-5	1/94
603 B: Mobile Phase Preparation	603-6	1/94
Filtering Solvents	603-6	1/94
Degassing Solvents	603-6	1/94

	<i>page</i>	<i>date</i>
Preparation of Multisolvent Mobile Phases	603-7	1/94
Solvent Reservoirs	603-7	1/94
603 C: Mobile Phase Delivery Systems	603-8	1/94
Pumps	603-8	1/94
Gradient Programming Systems	603-10	1/94
603 D: Maintenance and Troubleshooting	603-10	1/94
Problems with Pumps	603-10	1/94
References	603-13	1/94
<b>604: Injection Systems</b>		
604 A: Injection Valves	604-1	1/94
604 B: Automatic Injectors	604-3	1/94
604 C: Operation, Maintenance, Troubleshooting, and Repair of Injection Valves	604-3	1/94
References	604-4	1/94
<b>605: Detectors</b>		
605 A: UV/VIS Absorbance Detectors	605-2	1/94
Fixed Wavelength UV Detectors	605-3	1/94
Variable Wavelength UV Detectors	605-4	1/94
Solvents	605-4	1/94
Performance Characteristics	605-4	1/94
Multichannel or Photodiode Array Detectors	605-5	1/94
Applications	605-5	1/94
Problems, Maintenance, and Troubleshooting	605-5	1/94
605 B: Fluorescence Detectors	605-6	1/94
Detector Design	605-6	1/94
Solvents	605-7	1/94
Performance Characteristics	605-7	1/94
Parameter Adjustments	605-8	1/94
Applications	605-8	1/94
Detector Maintenance	605-8	1/94
605 C: Electrochemical Detectors	605-8	1/94
Conductivity Detectors	605-9	1/94
Amperometric and Coulometric Detectors	605-9	1/94
Performance Characteristics	605-10	1/94
Applications	605-11	1/94
605 D: Photoconductivity Detectors	605-11	1/94
Apparatus	605-11	1/94
Performance Characteristics	605-11	1/94
Applications	605-13	1/94
605 E: Mass Spectrometric Detectors	605-13	1/94
605 F: Derivatization for Detection Enhancement	605-14	1/94
Comparison of Pre- and Post-Column Derivatization	605-14	1/94

	<i>page</i>	<i>date</i>
Post-Column Reactor Design	605-14	1/94
References	605-16	1/94
<b>606: Residue Identification and Quantitation</b>		
606 A: Residue identification	606-1	1/94
Co-chromatography	606-1	1/94
Use of Alternative Columns	606-2	1/94
Spectrometric Confirmation	606-2	1/94
606 B: Quantitation	606-3	1/94
Reference	606-3	1/94
<b>607: Quality Assurance and Troubleshooting</b>		
607 A: Liquid Chromatograph Monitoring and Performance Testing	607-1	1/94
607 B: Troubleshooting from Chromatograms	607-2	1/94
<b>608: Bibliography</b>		
General Texts	608-1	1/94
Columns	608-1	1/94
Detectors	608-2	1/94
Troubleshooting	608-2	1/94
Application to Pesticides	608-2	1/94
<b>Figures</b>		
601-a Chromatographic Separation Techniques	601-1	1/94
601-b HPLC Modes of Operation	601-2	1/94
601-c Guide to Selection of HPLC Mode	601-5	1/94
601-d Block Diagram of HPLC System	601-6	1/94
601-e Column Outlet Fittings	601-7	1/94
601-f Low Dead Volume Fitting	601-7	1/94
601-g Standard Internal Fitting	601-8	1/94
601-h Internal Thread Low Dead Volume Fitting	601-8	1/94
602-a Calculation of Column Performance Parameters	602-5	1/94
603-a Reciprocating Pump	603-8	1/94
603-b Gradient System I	603-9	1/94
603-c Gradient System II	603-9	1/94
603-d Gradient System III	603-10	1/94
604-a External Loop Injector: Six-Port Injection Valve	604-1	1/94
604-b Internal Loop Injector	604-2	1/94
605-a UV/VIS Detector	605-3	1/94
605-b Variable Wavelength UV/VIS Detector	605-4	1/94
605-c Fluorescence HPLC Detector	605-6	9/96
605-d Three Electrode Electrochemical Detector	605-10	1/94
605-e Post-Column Reactors	605-15	1/94



# Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

## Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

## Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

## API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

## LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

## FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

## E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.