

Product Catalog





GUIDANT CORONARY STENT SYSTEMS	
CORONARY GUIDE WIRES	
CODONADY CHIDING CATHETERS	
CORONARY GUIDING CATHETERS	
CORONARY DILATATION CATHETERS	
ACCESSORIES	
DIRECTIONAL CORONARY ATHERECTOMY	



ORDERING INFORMATION

Return Policy

- 1. Guidant Corporation's Vascular Intervention division will accept return of any undamaged and unused standard catalog item within 90 days of invoice date. Defective product can be returned at any time, with an Authorized Return Goods Number (RGA).
- 2. Returned product is subject to examination for acceptable condition or alleged defects.
- 3. Acceptable returned product will be credited to customer account subject to a restocking fee.
- 4. A Returned Goods Authorization Number is required for all returns and replacements. This number can be obtained by calling your account representative, or Customer Service at 800 227-9902.
- 5. The RGA Number must be clearly marked on the outside of all boxes returned. A copy of the packing slip and/or Purchase Order Number should be included with the merchandise
- 6. Send non-complaint returns to:

42301 Zevo Drive

Suite E

Temecula, CA 92590

A biohazard box is required for all complaint returns. This box can be obtained by calling your account representative, or Customer Service.

Send complaint returns to: 26531 Ynez Road Temecula, CA 92591


Shipping Policy

Guidant Corporation's Vascular Intervention division offers a variety of shipping options, including overnight delivery. All shipment charges are paid by Guidant. Please specify your shipping preference when placing your order with Guidant Customer Service Representatives.

WARNING: Manufactured with CFC-113, a substance which harms public health and environment by destroying ozone in the upper atmosphere.

Guidant Corporation's Vascular Intervention division is diligently working to eliminate the use of CFC-113 in the manufacturing process.





Small Vessel



*ML PIXEL $^{\rm IM}$ is indicated for patients presenting with abrupt or threatened abrupt closure.

Small vessel stents combine low profiles with a flexible stent design specifically for treatment of vessels that are 2.0 mm to 2.5 mm in size.

ATTRIBUTES

- Thinner struts for increased flexibility and larger lumens
- Low system profiles for excellent deliverability and crossability
- Noncompliant balloon for controlled growth

Workhorse



The flexibility, conformability and deliverability of workhorse stents make them ideally suited for the majority of cases.

ATTRIBUTES

- Low profiles for easy lesion access, enhanced track and stent recross
- Clinically proven stent pattern designed to maintain vessel shape in tortuous anatomies while providing optimal scaffolding
- Minimal balloon overhang to reduce balloon interaction outside of the stent upon expansion

Large Vessel



The scaffolding and radial strength provided by large vessel stents make them a superb fit for large native coronary arteries and saphenous vein grafts.

ATTRIBUTES

- Low unsupported surface area and high metal-to-artery ratio for excellent scaffolding in large native vessels and saphenous vein grafts
- High radial strength to minimize recoil



DOCKET

Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.

