



US006415166B1

(12) **United States Patent**
Van Hoy et al.

(10) **Patent No.:** **US 6,415,166 B1**
(45) **Date of Patent:** **Jul. 2, 2002**

(54) **PHOTOPLETHYSMOGRAPHIC DEVICE WITH REMOTE FACSIMILE**

(75) Inventors: **Gilbert W. Van Hoy**, Broomfield;
Charles A. Gonzales, Westminster;
David L. Newcomb, Louisville;
Michael K. Brashears, Denver; **Tricia A. Dessel**, Boulder, all of CO (US)

5,348,004 A 9/1994 Hollub
5,515,176 A * 5/1996 Galen et al. 128/904
5,544,649 A * 8/1996 David et al. 128/904
5,581,369 A * 12/1996 Righter et al. 128/904
5,701,894 A 12/1997 Cherry et al.

* cited by examiner

(73) Assignee: **Datex-Ohmeda, Inc.**, Louisville, CO (US)

Primary Examiner—Eric F. Winakur
(74) *Attorney, Agent, or Firm*—Marsh Fischmann & Breyfogle LLP

(*) Notice: Subject to any disclaimer, the term of this patent is extended or adjusted under 35 U.S.C. 154(b) by 0 days.

(57) **ABSTRACT**

(21) Appl. No.: **09/386,691**

(22) Filed: **Aug. 30, 1999**

(Under 37 CFR 1.47)

Related U.S. Application Data

(63) Continuation-in-part of application No. 08/938,224, filed on Sep. 26, 1997, now abandoned.

(51) **Int. Cl.⁷** **A61B 5/00**

(52) **U.S. Cl.** **600/323; 600/333; 128/904**

(58) **Field of Search** 600/300, 301, 600/310, 322, 323, 324, 333; 128/903, 904


A photoplethysmographic monitoring system such as a pulse oximeter collects data regarding the blood analyte concentration and pulse rate of a patient through the analysis of light transmitted from an emitter through tissue to a photo detector. It is often necessary to review the collected data, such as oxygen saturation, pulse rate and pulsatility value at a location remote to the patient being monitored. The photoplethysmographic system formats the selected data for transmittal to a remote facsimile machine via standard telephone communication systems using an internal or external modem. The formatter is able to function within the processing and memory constraints of pulse oximeters by formatting and transmitting the data in sections. Formatted data may be sent via hard-wired telephone, cellular phone, PCS digital telephones or through satellite communication systems enabling transmittal of data from a portable photoplethysmographic system.

(56) **References Cited**

U.S. PATENT DOCUMENTS

5,078,136 A 1/1992 Stone et al.

19 Claims, 20 Drawing Sheets






FIG.1




FIG. 2




FIG.3




FIG. 4

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.