

Biotechnology and Applied Biochemistry

HEALTH SCIENCES LIBRARY
University of Wisconsin

JUN 20 1994

Editor in Chief

P. N. Campbell

Department of Biochemistry and Molecular Biology, University College, Gower Street, London WC1E 6BT, U.K.

Regional Editors

J. R. DeLoach

USDA, Agriculture Research Service, Route, 5, Box 810, College Station, TX 77845, U.S.A.

Editors

A. Aitken

National Institute for Medical Research, Mill Hill, London, U.K.

C. Bucke

University of Westminster, London, U.K.

A. De Flora

Università Degli Studi de Genova, Genova, Italy

G. W. Gould

Unilever Limited, Sharnbrook, Bedford, U.K.

M.-R. Kula

Institut für Enzymtechnologie der Universität Düsseldorf, Jülich, Germany

B. Mattiasson

University of Lund, Lund, Sweden

J. D. Newman

Cranfield Institute of Technology, Cranfield, Bedford, U.K.

D. O'Shannessy

SmithKline Beecham, King of Prussia, U.S.A.

T. Oshima

Tokyo Institute of Technology, Yokohama, Japan

W. Ostrowski

Instytut Biochemii Lekarskiej, Kraków, Poland

Honorary Editors **E. Katchalski-Katzir** • **J. Lederberg** • **J. Porath**

Staff Editor **Stuart Hobday**

Assistant Editor

P. M. Brickell

Medical Molecular Biology Unit, Department of Molecular Pathology, UCL Medical School, London W1P 6DB, U.K.

K. Soda

Laboratory of Microbial Biochemistry, Institute for Chemical Research, Kyoto University, Gokasho, Uji, Kyoto-fu 611, Japan

N. A. Rao

Indian Institute of Science, Bangalore, India

E. S. Severin

Research Center of Molecular Diagnostics, Ministry of Health, Moscow, Russia

B. A. Stone

La Trobe University, Bundoora, Australia 3083

B. Svensson

Carlsberg Laboratory, Valby, Denmark

Y. H. Tan

National University of Singapore, Singapore

M. Uhlen

Royal Institute of Technology, Stockholm, Sweden

S. D. Varfolomeyev

M. V. Lomonosov Moscow University, Moscow, Russia

R. Virden

The Medical School, Newcastle upon Tyne, U.K.

J. R. Wild

Texas A&M University, College Station, TX, U.S.A.

1305 Linden Drive
Madison, WI 53706

-
- *Biotechnology and Applied Biochemistry* (ISSN 0885-4513) is published bimonthly in February, April, June, August, October and December by Portland Press Ltd, 59 Portland Place, London W1N3AJ, U.K., on behalf of the International Union of Biochemistry and Molecular Biology.
 - Subscription orders should be sent to Portland Press Ltd, P.O. Box 32, Commerce Way, Colchester CO2 8HP, U.K., or to your usual agent. Volumes 19 and 20 (three issues each) will appear in 1994. The institutional subscription rate for 1994 is £100, or US\$175 for subscribers in North America. Subscription rates for North America and Japan include airfreight and Accelerated Surface Post delivery respectively. Air-mail delivery is available at additional cost; details are available on application. Special personal subscription rates for members of the Biochemical Society are also available.
 - All rights reserved. Apart from any fair dealing for the purposes of research or private study, or criticism or review, as permitted under the Copyright, Designs and Patents Act, 1988, this publication may be reproduced, stored or transmitted, in any forms or by any means, only with the prior permission in writing of the publishers, or in the case of reprographic reproduction in accordance with the terms of permissions granted by copyright clearing houses such as the Copyright Licensing Agency (U.K.) or the Copyright Clearance Center (U.S.A.). The CCC item-fee code for this publication is 0885-4513/94/\$5.00 + 0.00. Inquiries concerning reproduction outside those terms should be sent to the publishers at the London address.
 - Although, at the time of going to press, the information contained in this publication is believed to be correct, the authors, the editors, the publisher and the International Union of Biochemistry and Molecular Biology assume no responsibility for any errors or omissions herein contained. All material published herein represents the opinions of the author(s) and not necessarily those of the editors, the publishers or the International Union of Biochemistry and Molecular Biology.
 - Typeset by Unicus Graphics Ltd., Horsham, and Printed in Great Britain by Whitstable Litho Printers Ltd., Kent.
 - Second-class postage paid at Rahway, NJ, U.S.A.
 - U.S. Postmaster: send address changes to *Biotechnology and Applied Biochemistry*, c/o Mercury Airfreight International Ltd, 2323 Randolph Avenue, Avenel, NJ 07001.
-

Purification and properties of a β -D-galactosidase from *Bifidobacterium bifidum* exhibiting a transgalactosylation reaction¹

Vincent Dumortier,* Colette Brassart,† Stéphanie Bouquelet²

*Département Agro-Alimentaire et Biotechnologies, Institut Supérieur d'Agriculture de Lille, Lille, France, and †Laboratoire de Chimie Biologique, Unité mixte du CNRS n° 111 Université des Sciences et Technologies de Lille, 59655 Villeneuve d'Ascq Cédex, France

β -D-Galactosidase was extracted by ultrasonic treatment in the presence of Triton X-100. The transgalactosylase (TGase) was separated from other β -D-galactosidases and purified 1720-fold, as measured by the TGase assay. The enzyme was shown to be homogeneous on native PAGE. The complex nature of the enzyme was demonstrated by SDS/PAGE, with four major protein bands of 163, 170, 178 and 190 kDa. The molecular mass was estimated to be 362 kDa by gel chromatography, and the pI was 5.25. The purified enzyme was stable at temperatures below 45°C and over the pH range from 6.5–8. Lactose hydrolysis by the purified enzyme and by *Bifidobacterium bifidum* cells took place at the same pH and temperature, at pH 6.5 and 37°C respectively. For the TGase, however, *B. bifidum* cells exhibited specific galactosyltransferase activity at pH 4.5 and at 40°C, even though the transgalactosylation (TG) reaction was preferentially established by the purified enzyme at pH 4.8 and 45°C. The amount of galacto-oligosaccharides (GOS) could be as high as 29% of the total sugar, and maximum production was achieved when 60% of the initial lactose was transformed by *B. bifidum* cells. The K_m values of the purified β -D-galactosidase were 2 mM for p-nitrophenyl galactosidase and 13 mM for lactose. In the TG reaction, the Michaelis constant was 800 mM. Calcium and EDTA affected the TGase/hydrolase ratio. The purified enzyme exhibited a specificity towards glycosidic acceptors, and all tested β -D-galacto-oligosides can be used as galactosyl residue donor to produce GOS.

β -D-Galactosidase (EC 3.2.1.23) is an enzyme which is widespread in Nature [1] and is able to catalyse not only the hydrolysis of β -D-galactosidic linkages, but also the transgalactosylation (TG) reaction to produce oligosaccharides [2–14].

Transfer reaction products do not seem to have a predominant physiological role, even though allolactose formed by *Escherichia coli* galactosidase is a natural inducer of the *lac* operon [15]. Galacto-oligo-

¹This work was supported in part by the Centre National de la Recherche Scientifique UMR n° 111 (Relations structure-fonction des Constituants membranaires; Director: Professeur André Verbert) and by the Université des Sciences et Technologies de Lille.

²To whom correspondence should be sent.

³Abbreviations used: pNP-Gal, p-nitrophenyl β -D-galactopyranoside; GOS, galacto-oligosaccharide; TG, transgalactosylation; TGase, transgalactosylase.

saccharides (GOS) might be useful as a bifidobacterial growth factor [16] or for attenuating the intestinal disorders of elderly patients [17]. Nevertheless, production of GOS by the TG reaction has recently become of interest. Glycoside moieties of glycoconjugates participate in recognition phenomena, and the synthesis of related structures can be performed using β -galactosidase [18]. The structure elucidation of TG products showed generally the prevailing formation of $\beta(1,6)$ linkages. Nevertheless, the synthesis reaction also can form other linkages, such as $\beta(1,4)$ galactosidic linkage [19].

In previous papers it was reported that *B. bifidum*, a bacterium from human colon, is able to exert a beneficial influence on the health of breast-fed infants [20–21]. We showed that a β -D-galactosidase from *B. bifidum* exhibited a transgalactosylase (TGase) activity [22]. The structure determination of ten GOS synthesized from lactose showed the predominance of $\beta(1,3)$ linkages; such GOS had never previously been shown to be products of TG. Hence, it was of interest to know whether *B. bifidum* was able to synthesize *in vivo* such GOS structures. The GOS production was optimized, and the properties of the enzyme are described and discussed. Previous experiments indicated the presence of different β -D-galactosidases in *B. bifidum* cells [23]. The purification of these enzymes allowed us to determine whether all or only one β -D-galactosidase was responsible for the specific TG reaction.

Experimental

Materials

BSA, *p*-nitrophenyl glycosides and chemicals were purchased from Sigma. Lactose, glucose and galactose were obtained from Merck. Hydroxyapatite Ultrogel, and Q Fast Flow, Mono P, Superose 6 and Polybuffer 74 were products of Sepracor and Pharmacia respectively. Glucose 6-phosphate dehydrogenase, hexokinase, NADP, ATP, and $MgSO_4$ were purchased from Boehringer Mannheim. Aminex HPX 42C columns (300 mm \times 7.8 mm), Bio-Rad protein assay reagent and molecular-mass standards for SDS/PAGE were obtained from Bio-Rad, and molecular-mass standards for gel chromatography were from Pharmacia.

Organisms

Cultivation conditions for *B. bifidum* (DSM 20082) were optimized by Boutry [24]. Cultures were grown in anaerobic medium in a Setric fermenter with a final working volume of 200 litres at constant pH 6 and at 37 °C. Cells were harvested by centrifugation at 7 000 *g*, 10 h after the end of the exponential phase of growth.

Protein assay

Protein was measured by the method of Bradford [25], with BSA as standard.

Enzyme assays

Hydrolytic activity was determined with 1.25 mM *p*-nitrophenyl galactoside (pNP-Gal) as substrate in 0.1 M sodium phosphate buffer, pH 6.5. When lactose was used as the substrate, the glucose was assayed by the method of Finch [26]. One unit of β -D-galactosidase was defined as the amount of enzyme producing 1 μ mol of *p*-nitrophenol or 1 μ mol of glucose/min at 37 °C and at pH 6.5.

TGase activity was measured by an h.p.l.c. method. The incubation mixture was composed of 100 μ l of diluted enzyme, 300 μ l of 0.05 M sodium citrate/0.1 M disodium phosphate buffer, pH 4.25, and 100 μ l of 0.5 M lactose. After 4 h at 45 °C, the enzymic reaction was stopped by heating for 2 min on a boiling-water bath, then the incubation mixture was centrifuged for 30 min at 35 000 *g*. The supernatant was diluted four times with water before injection (10 μ l) on to an h.p.l.c. column (Aminex HPX 42C) thermostatically controlled at 85 °C in a SP 8100 apparatus (Thermo Separation Products). Oligosaccharides were eluted by de-aerated water at a flow rate of 0.75 ml/min and detected by refractometry. Oligosaccharide concentrations were calculated by using 10 μ l of external standard (20 mM lactose) as reference. The results were expressed as trisaccharide equivalents. A unit of TGase was defined as the amount of enzyme catalysing the formation of 1 μ mol of trisaccharide/min at 45 °C and at pH 4.25.

Determination of optimal conditions

The optimal temperature was determined according to the standard procedure at temperatures ranging from 4 to 60 °C. For the TG reaction, the enzyme was incubated with 0.24 M lactose in 0.05 M sodium citrate/0.1 M disodium phosphate, pH 4.8, for 20 h in the same range of temperature. Saccharides were separated by t.l.c. on Kieselgel 60 (Merck) in the solvent system butanol/water/acetic acid (2:1:1, v/v). Staining was done by spraying orcinol (0.2%) in 20% (v/v) H₂SO₄ and heating for 5 min at 120 °C. Saccharides were estimated by densitometry (Bio-Rad 620 densitometer).

The optimal pH for hydrolytic activity was determined using pNP-Gal as substrate in 0.1 M sodium citrate/phosphate/borate buffer, which had a pH range from 4 to 10. This buffer was chosen because it covers a wide range of pH without changing the nature of the salt. For the TG reaction, 50 μ l of β -D-galactosidase were incubated with 0.24 M lactose at different pH values. After 20 h at 45 °C, the incubation was stopped by heating and the reaction products were estimated with TLC, as described above.

Influence of lactose concentration on GOS production

Cells were used for kinetic studies with various lactose concentrations (100, 250 or 400 mM) in 0.05 M sodium citrate/0.1 M disodium phosphate buffer, pH 4.25. Samples (0.4 ml) were withdrawn at various times from 30 min up to 48 h and the reactions stopped by heating. After centrifugation and dilution with water, the reaction products were analysed by h.p.l.c.

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.