WO 2011/109433

PCT/US2011/026741

38. The method of Claim 33, wherein the section of tissue comprises an ultimate tensile strength of greater than about 25 MegaPascals.

39. The method of Claim 33, wherein the section of tissue comprises a treated pericardium tissue.

5

40. A method of preparing a section of tissue for medical use, comprising:

(a) contacting the section of tissue with distilled water;

(b) contacting the section of tissue with isopropyl alcohol for a pre-fixation period of time of not less than about 3 days; and

(c) contacting the section of tissue with one of

10

(i) a formalin solution, or

(ii) a glutaraldehyde solution

for a fixation period of time of not less than about 3 days; and

(d) contacting the section of tissue with isopropyl alcohol for a post-fixation period of time of not less than about 3 days;

wherein step (b) occurs sometime after step (a), wherein step (c) occurs sometime after step (b), and wherein step (d) occurs sometime after step (c).

41. The method of Claim 40, wherein for step (c):

if the formalin solution is used, then the formalin solution comprises a concentration of

about 1 - 37.5% formalin; and

20

15

if the glutaraldehyde solution is used, then the glutaraldehyde solution comprises a concentration of about 0.1 - 25% glutaraldehyde.

42. The method of Claim 40, wherein for step (c):

if the formalin solution is used, then the formalin solution comprises a concentration of about 8-12% formalin; and

25

OCKE.

if the glutaraldehyde solution is used, then the glutaraldehyde solution comprises a concentration of about 0.1-0.5% glutaraldehyde.

43. The method of Claim 40, wherein the section of tissue comprises a treated pericardium tissue.








FIG. 1

DOCKET A L A R M Find authenticated court documents without watermarks at <u>docketalarm.com</u>. 2/5




DOCKET A L A R M Find authenticated court documents without watermarks at <u>docketalarm.com</u>.


DOCKET



DOCKET

4/5





A L A R M Find authenticated court documents without watermarks at <u>docketalarm.com</u>.

DOCKET A L A R M



Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.