

# Guanylin peptides: cyclic GMP signaling mechanisms

L.R. Forte<sup>1,2</sup>,  
R.H. Freeman<sup>3</sup>,  
W.J. Krause<sup>4</sup> and  
R.M. London<sup>1,2</sup>

<sup>1</sup>Harry S. Truman Veterans' Hospital, Departments of  
<sup>2</sup>Pharmacology, <sup>3</sup>Physiology and <sup>4</sup>Pathology and Anatomical Sciences,  
School of Medicine, Missouri University, Columbia, MO, USA

## Abstract

Guanylate cyclases (GC) serve in two different signaling pathways involving cytosolic and membrane enzymes. Membrane GCs are receptors for guanylin and atriopeptin peptides, two families of cGMP-regulating peptides. Three subclasses of guanylin peptides contain one intramolecular disulfide (lymphoguanilin), two disulfides (guanylin and uroguanylin) and three disulfides (*E. coli* stable toxin, ST). The peptides activate membrane receptor-GCs and regulate intestinal Cl<sup>-</sup> and HCO<sub>3</sub><sup>-</sup> secretion via cGMP in target enterocytes. Uroguanylin and ST also elicit diuretic and natriuretic responses in the kidney. GC-C is an intestinal receptor-GC for guanylin and uroguanylin, but GC-C may not be involved in renal cGMP pathways. A novel receptor-GC expressed in the opossum kidney (OK-GC) has been identified by molecular cloning. OK-GC cDNAs encode receptor-GCs in renal tubules that are activated by guanylins. Lymphoguanilin is highly expressed in the kidney and heart where it may influence cGMP pathways. Guanylin and uroguanylin are highly expressed in intestinal mucosa to regulate intestinal salt and water transport via paracrine actions on GC-C. Uroguanylin and guanylin are also secreted from intestinal mucosa into plasma where uroguanylin serves as an intestinal natriuretic hormone to influence body Na<sup>+</sup> homeostasis by endocrine mechanisms. Thus, guanylin peptides control salt and water transport in the kidney and intestine mediated by cGMP via membrane receptors with intrinsic guanylate cyclase activity.

## Key words

- Kidney
- Intestine
- Guanylate cyclase
- Chloride secretion
- Sodium excretion

## Correspondence

L.R. Forte  
Department of Pharmacology  
School of Medicine  
Missouri University  
M-515 Medical Sciences Building  
Columbia, MO 65212  
USA  
Fax: + 1-573-884-4558  
E-mail: lrf@missouri.edu

Presented at the Meeting  
"NO Brazil, Basic and Clinical  
Aspects of Nitric Oxide",  
Foz do Iguaçu, PR, Brazil,  
March 10-13, 1999.

Received May 28, 1999  
Accepted June 22, 1999

## Introduction

Guanylin, uroguanylin and lymphoguanilin are heat-stable peptides that regulate the enzymatic activity of cell-surface guanylate cyclase signaling molecules. Guanylin was the first endogenous peptide identified and was isolated as a 15 amino acid peptide from rat intestine that stimulates cGMP production in T84 intestinal cells (1). Uroguanylin was next isolated from opossum urine

as biologically active, 13, 14 and 15 amino acid peptides using the T84 cell cGMP bioassay (2). Lymphoguanilin was recently identified using a PCR-based homology cloning strategy to isolate cDNAs from opossum spleen and other lymphoid tissues, which encode a 109 amino acid precursor that is most similar in its primary structure to preprouroguanylin (3). Guanylin and uroguanylin have two intramolecular disulfide bonds, but lymphoguanilin has only one disulfide

bond. A synthetic form of lymphoguanilin also activates guanylate cyclase (GC)-receptors in human T84 intestinal and opossum kidney (OK) cells. The family of guanylin regulatory peptides share similarities both in primary structures and biological activities with the heat-stable enterotoxin (stable toxin, ST) peptides secreted by strains of enteric microorganisms that cause a watery form of diarrhea similar to that seen in cholera (4). Thus, bacteria-derived ST peptides are molecular mimics of guanylin peptides acting to stimulate fluid and electrolyte secretion into the lumen of the intestine by activation of native guanylin receptors located on the apical surfaces of enterocytes lining the intestine. Biologically active peptides in the guanylin family of proteins are found at the COOH-termini of longer precursor polypeptides that are secreted as biologically inactive prohormones or protoxins (5-8). Proteolytic enzymes serve as converting enzymes to activate the proguanylins or proSTs, but these important enzymes have not been identified thus far.

### **Physiological actions of guanylin peptides**

Physiological actions of the endogenous guanylin peptides include the regulation of intestinal fluid secretion during digestion, and neutralization of HCl in the duodenum and of organic acids derived from enteric bacteria in the large intestine (9,10). Guanylins and STs stimulate the electrogenic secretion of both chloride and bicarbonate anions, which provides the physiological driving force to accomplish fluid secretion into the intestinal lumen. Control of these intestinal functions by guanylin peptides is mediated via intracellular cGMP through activation of cGMP-dependent protein kinase II (cG-kinase II) and/or cAMP-dependent protein kinase II (cA-kinase II) with subsequent phosphorylation of the cystic fibrosis transmembrane conductance regulator (CFTR)

protein (11,12). CFTR and GC-receptors are localized together in apical plasma membranes of target enterocytes. The first GC-receptor for guanylin family peptides was identified at the molecular level (i.e., GC-C) by molecular cloning of cDNAs encoding an intestinal membrane protein that binds  $^{125}\text{I}$ -ST with high affinity and is activated by ST (13). Transgenic suckling mice with disabled GC-C genes have a marked reduction in the intestinal fluid secretion response to *E. coli* ST, indicating that GC-C is responsible for a major component of the intestinal actions of guanylin peptides to enhance fluid secretion (14,15). However, about 10% of specific  $^{125}\text{I}$ -ST binding to receptor sites on intestinal membranes still remain functional in GC-C-knock out (GC-C-KO) animals (15). This suggests that an additional gene or multiple genes encoding GC-receptors for guanylin agonists exist in the mouse genome.

Bacterial STs were the first peptides shown to activate GC-receptor signaling molecules in the intestine, thus causing secretory diarrhea (4). Moreover, the intestinal GC targets for ST peptides in the intestine were considered unique and not existing outside the intestinal epithelium until we discovered that *E. coli* ST also stimulates GC-receptors found on the surface of OK and potaroo kidney (PtK-2) cell lines (Figure 1; Ref. 16). It can be seen from one of our original experiments that these two kidney cell lines have remarkable cGMP responses to *E. coli* ST, but only small cGMP responses to atriopeptin-A and no detectable responses to the nitric oxide donor, sodium nitroprusside. The additional demonstration that kidney cortex has specific GC-receptors that are activated by ST implies that this class of cell surface GC-receptors plays a much broader physiological role in the body than was previously recognized (17-19). It should be emphasized that these findings were made well before the first guanylin peptides were isolated (1,2). Putative "ST

receptors” were also found in other epithelia of the opossum in addition to the receptors localized to brush border membranes (BBM) of epithelial cells lining the intestinal tract and within renal tubules (16-19). The existence of renal, hepatic, airway and testicular “ST receptors” predicted that endogenous ST-like peptides exist to regulate the activity of the extra-intestinal as well as intestinal GC-receptors. These seminal experiments led directly to the subsequent isolation of guanylin and uroguanylin peptides from intestine and urine, respectively (1,2). Both guanylin and uroguanylin are produced in the intestine, but uroguanylin is the major bioactive peptide found in urine, which contains either no guanylin or very small amounts of the peptide (2,20-22). Uroguanylin and ST stimulate the enzymatic activity of renal tubular GC-receptors and increase the urinary excretion of sodium, potassium and water in both the perfused rat kidney *ex vivo* and the mouse *in vivo* (16-19,23-25). Guanylin is less potent in the stimulation of urinary  $\text{Na}^+$  and water excretion compared to either uroguanylin or ST, but guanylin does have marked kaliuretic activity in the perfused rat kidney (24). Thus, uroguanylin has biological activity consistent with a peptide hormone that influences renal function by regulating the urinary excretion of sodium chloride as a physiological mechanism that contributes to the maintenance of  $\text{Na}^+$  balance in the body. A natriuretic peptide such as uroguanylin was predicted to exist in the digestive system for release into the bloodstream for the purpose of stimulating the urinary excretion of  $\text{NaCl}$  following a salty meal (26,27). Uroguanylin is a prime candidate for this “intestinal natriuretic hormone” because it is produced at extraordinarily high concentrations in the upper small intestine and is released following a high salt meal (28-32). Uroguanylin mRNAs are most abundant in the small intestine compared to guanylin mRNA levels, which peak in the large intestine. Uroguanylin and prourogua-

nylin are also found in the circulation and it is likely that the gastrointestinal (GI) tract is a main source of the plasma peptides (7,33, 34). Secretion of uroguanylin from GI mucosa into the plasma in response to oral  $\text{NaCl}$  may explain the prolonged increase in urinary sodium excretion that occurs following a high salt meal (26,27).

### Identification of a kidney GC-receptor for uroguanylin

We sought to elucidate the primary structure of a membrane GC expressed in cultured OK cells and in the opossum kidney because the prior discovery of this renal GC-receptor was a stimulus that ultimately led to the isolation of guanylin and uroguanylin (1,2,16-19). A PCR-based cloning strategy was used to isolate 3762-bp cDNAs from RNA/cDNAs expressed in OK cells and opossum kidney cortex (35,36). Transfection and expression of the OK-guanylate cyclase (OK-GC) cDNA into COS and HEK293 cells produces a cell surface GC-receptor of ~160 kDa size that is activated by uroguanylin, guanylin and *E. coli* ST peptides. OK-GC cDNA contains an open reading frame encoding a 1049 residue mature




Figure 1 - Activation of surface receptors on OK and PtK-2 cells by Escherichia coli stable toxin (ST) and atriopeptin-A. Confluent cells were treated for 15 min at 37°C in DMEM-HEPES-MIX with 0.5  $\mu\text{M}$  ST, 0.1  $\mu\text{M}$  atriopeptin-A (ANP-A) or 1 mM sodium nitroprusside (Na-NP). The data are reported as the mean of 4 experiments with opossum kidney (OK) cells and 2 experiments with potoroo kidney (PtK-2) cells.

protein belonging to the family of membrane GC-receptor signaling molecules. OK-GC is similar to other membrane GC proteins containing NH<sub>2</sub>-terminal agonist-binding domains, a single membrane span and intracellular kinase-like and GC catalytic domains. OK-GC is 72, 76 and 75% identical in its overall structure compared to the intestinal BBM-localized GC-C receptors for guanylin peptides found in rats, humans and pigs, respectively (13,37,38). The catalytic domains of OK-GC and GC-C receptors of rat, human and porcine intestine share 92, 94 and 95% identity, respectively. The most highly variable region of membrane GC-receptors occurs within the NH<sub>2</sub>-terminal ligand-binding domains of these proteins. OK-GC shares only 55-59% identity in this domain when compared to GC-C intestinal receptors for guanylin peptides. The GC-C receptors of the rat are more closely related to human and pig GC-C in the ligand-binding domains with these proteins sharing ~70% identity in this region. Thus, OK-GC is a distinctive renal GC-receptor that, along with the intestinal GC-C receptor, provides two distinctive molecular subtypes of GC-receptors for a growing family of membrane receptors for the guanylin peptides.

OK-GC mRNA levels were measured in total RNA prepared from tissues of opossums by Northern and RT-PCR. A 3.8-kb mRNA was detected in kidney, OK cells, urinary bladder, adrenal gland, heart and intestine. Lower levels of OK-GC mRNA were detected in renal medulla compared to cortex. Tissues with lower levels of OK-GC mRNA are urinary bladder, adrenal gland and both the ventricles and atria of the heart. In the GI tract, high levels of OK-GC mRNA were measured in both the small and large intestine. Thus, OK-GC mRNA is highly expressed in the kidney cortex and intestinal mucosa with mRNA transcripts occurring at lower levels in a number of other organs.

OK-GC is a candidate for the renal tubular receptor that is activated by uroguanylin

and/or guanylin peptides that signals via cGMP to regulate the urinary excretion of sodium, potassium and water (23-25). OK-GC is the first kidney receptor for guanylin family peptides to be fully defined at the molecular level. Prior to the identification of GC-receptor signaling molecules for *E. coli* ST in the OK and PtK-2 cells and in the opossum kidney, it was thought that ST-stimulated GCs were restricted to intestinal mucosa (16-19). Thus, identification of OK-GC in OK cells and opossum kidney opened up a new field of inquiry that culminated recently with the discoveries of guanylin, uroguanylin and lymphoguanylin peptides (1-3). These peptides activate OK-GC and may serve as endogenous agonists for this membrane GC-receptor. The OK cell line has the differentiated properties of renal proximal tubules and receptor autoradiography experiments with <sup>125</sup>I-ST have clearly shown that specific binding sites for this uroguanylin-like radioligand are found in cells of both convoluted and straight portions of proximal tubules (16-19). High levels of <sup>125</sup>I-ST-labeled receptors are also found in BBMs isolated from the kidney and intestine, indicating that this GC-receptor is preferentially localized to apical plasma membranes of both kidney and intestinal target cells. It is likely that OK-GC serves as a physiological receptor for uroguanylin, which is the major bioactive guanylin family peptide in urine (2,20-22). However, guanylin was also isolated from opossum urine indicating that this peptide may influence renal function *in vivo* via cGMP (2). Active uroguanylin and inactive prouroguanylin peptides circulate in plasma, thus providing a source for the urinary forms of bioactive uroguanylin in opossum, human and rat urine (2,20-22). The kidney also expresses uroguanylin and guanylin mRNAs, which offers the possibility that an intra-renal signaling pathway exists for the guanylin peptides as well as providing another potential source of uroguanylin in the urine (35,39). Physi-

ologically, uroguanylin may regulate kidney function via an endocrine axis linking the intestine to the kidney and/or through an intra-renal paracrine mechanism. Both possibilities involve the activation of OK-GC as a key signaling molecule in renal tubular target cells that possess the guanylate cyclase-cGMP signaling machinery.

### Identification of lymphoguanylin

A third member of the guanylin family of peptides was sought since biologically active uroguanylin peptides were isolated from opossum urine (2). This inquiry culminated in the isolation of cDNAs encoding opossum prepolymphoguanylin, a third and unique member of the guanylin family of endogenous regulatory peptides (3). The identification of lymphoguanylin stems from recent experiments showing that mRNA transcripts for guanylin and/or uroguanylin and their GC-receptors are expressed broadly in tissues of the digestive, renal-urinary, cardiovascular, reproductive, lymphoid-immune and central nervous organ systems (7,35). Moreover, hybridization assays revealed the existence of a novel mRNA transcript that hybridizes with both uroguanylin and guanylin cDNA probes, but is substantially longer than either the 1.2-kb uroguanylin or 0.8-kb guanylin mRNAs (3). We identified a third guanylin-related mRNA transcript by molecular cloning using a homology cloning strategy based on the PCR with cDNAs isolated from opossum tissues. The deduced polypeptide was named prepolymphoguanylin because the first cDNAs isolated were derived from several different lymphoid organs and the encoded polypeptide is similar in primary structure to preprouroguanylin and preproguanylin (Figure 2). Sequence analyses of multiple independent cDNA clones reveals that the lymphoguanylin cDNA is 92.7% identical to the corresponding nucleotide sequences for preprouroguanylin reported previously (7). An open read-

ing frame within the lymphoguanylin cDNAs encode a 109 amino acid polypeptide that is 84% identical to preprouroguanylin and 40% identical to preproguanylin (7,35). At the COOH-terminus of the 109 amino acid prepolymphoguanylin is a 15 amino acid peptide that is 80% identical to uroguanylin, but shares only 40% identity with guanylin (Figure 3). Illustrated for comparison are the structures of opossum lymphoguanylin, uroguanylin, guanylin and an *E. coli* ST peptide, which form three different subclasses of guanylin peptides based on the number of intramolecular disulfides found within the active peptides. A major difference within lymphoguanylin is the tyrosine<sup>109</sup> residue because guanylins and uroguanylins have



Figure 2 - Prepeptide structures for guanylin, uroguanylin and lymphoguanylin. The amino acid sequences of the active peptides are shown using the single letter abbreviation of the amino acids within the COOH-terminal regions of each opossum polypeptide.




Figure 3 - Comparison of the primary structures of subclasses within the guanylin family of biologically active peptides. The peptide subclasses are provided from top to bottom in the order of their discovery as: Class I peptides containing 3 disulfides, class II peptides with 2 disulfides and class III peptides that have only 1 intramolecular disulfide bond. The single letter abbreviation is used to designate amino acid residues in each peptide. *E. coli* ST, *E. coli* stable toxin.

# Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

## Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

## Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

## Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

## API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

## LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

## FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

## E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.