



US005259381A

United States Patent [19]

[11] Patent Number: **5,259,381**

Cheung et al.

[45] Date of Patent: * **Nov. 9, 1993**

[54] **APPARATUS FOR THE AUTOMATIC CALIBRATION OF SIGNALS EMPLOYED IN OXIMETRY**

FOREIGN PATENT DOCUMENTS

0104772 4/1984 European Pat. Off. 128/633

[75] Inventors: **Peter W. Cheung**, Mercer Island; **Karl F. Gauglitz**, Kirkland; **Scott W. Hunsaker**, Seattle; **Stephen J. Prosser**, Lynnwood; **Darrell O. Wagner**, Monroe; **Robert E. Smith**, Edmonds, all of Wash.

Primary Examiner—Lee S. Cohen
Assistant Examiner—Robert L. Nasser, Jr.
Attorney, Agent, or Firm—Christensen, O'Connor, Johnson & Kindness

[73] Assignee: **Physio-Control Corporation**, Redmond, Wash.

[57] ABSTRACT

[*] Notice: The portion of the term of this patent subsequent to Apr. 3, 2007 has been disclaimed.

Under the present invention, a method and apparatus are provided for compensating for the effect temperature variations have on the wavelength of light emitted by the oximeter sensor light source (40, 42). In pulse oximetry, LEDs (40, 42) are typically employed to expose tissue to light at two different wavelengths. The light illuminating the tissue is received by a detector (38) where signals proportional to the intensity of light are produced. These signals are then processed by the oximeter circuitry to produce an indication of oxygen saturation. Because current oximetry techniques are dependent upon the wavelengths of light emitted by the LEDs (40, 42), the wavelengths must be known. Even when predetermined combinations of LEDs (40, 42) having relatively precise wavelengths are employed, variations in the wavelength of light emitted may result. Because the sensor (12) may be exposed to a significant range of temperatures while in use, the effect of temperature on the wavelengths may be significant. To compensate for this effect, a temperature sensor (50) is included in the sensor (12) to produce a signal indicative of sensor temperature. This signal is interpreted by the oximeter circuitry including, for example, a microcomputer (16), where the effect of temperature on wavelength is compensated for. In a preferred arrangement, this compensation takes the form of a computation of an alternative calibration curve from which the oxygen saturation is indicated, given a particular processing of signals from the detector (38).

[21] Appl. No.: 377,722

[22] Filed: Jul. 10, 1989

Related U.S. Application Data

[63] Continuation of Ser. No. 897,663, Aug. 18, 1986, Pat. No. 4,913,150.

[51] Int. Cl.⁵ **A61B 5/00**

[52] U.S. Cl. **128/633; 128/666; 356/41**


[58] Field of Search 128/633, 634, 664-667; 330/59, 308; 250/214 A

[56] References Cited

U.S. PATENT DOCUMENTS

3,709.612	1/1973	Clemens	128/633
4,586.513	5/1986	Hamaguri	128/633
4,685,464	8/1987	Goldberger et al.	128/633
4,710.631	12/1987	Aotsuka et al.	250/354.1
4,723.554	2/1988	Oman et al.	128/664

3 Claims, 9 Drawing Sheets






Fig. 1.




Fig. 2.




Fig. 3.




Fig. 6.




Fig. 4.




Fig. 5.




Fig. 9.




Fig. 10.




Fig. 7.




Fig. 8.

Explore Litigation Insights

Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

Real-Time Litigation Alerts



Keep your litigation team up-to-date with **real-time alerts** and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

Advanced Docket Research



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

Analytics At Your Fingertips



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

LAW FIRMS

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

FINANCIAL INSTITUTIONS

Litigation and bankruptcy checks for companies and debtors.

E-DISCOVERY AND LEGAL VENDORS

Sync your system to PACER to automate legal marketing.