# US Endodontics, LLC's Hearing Presentation

**US Endodontics, LLC** 

Petitioner,

V.

**Gold Standard Instruments, LLC,** 

Patent Owner.

Case PGR2015-00019

Patent 8,876,991 B2

October 19, 2016



### **Endodontic File**





### Grounds Upon Which Trial Was Instit

ORDERED that pursuant to 35 U.S.C. § 324(a), a post-grant review is hereby instituted as to claims 12–16 of the '991 patent on the following grounds of unpatentability:

- 1. Claims 12–16 under 35 U.S.C. § 112(a) as being unpatentable for lack of enablement;
- 2. Claims 12–16 under 35 U.S.C. § 112(a) as being unpatentable for lack of written description;
- 3. Claims 12–16 under 35 U.S.C. § 102 as being anticipated by Luebke 2008;
- 4. Claim 15 under 35 U.S.C. § 103 as being unpatentable over Luebke 2008 alone or in view of Heath or ISO 3630-1;
- 5. Claims 12–14 and 16 under 35 U.S.C. § 102 as being anticipated by Kuhn; and
- 6. Claim 15 under 35 U.S.C. § 103 as being unpatentable over Kuhn alone or in view of Heath or ISO 3630-1.



# US Pat. No. 8,876,991 – Claim 12

- 12. A method for manufacturing or modifying an endodo tic instrument for use in performing root canal therapy on tooth, the method comprising:
  - (a) providing an elongate shank having a cutting edextending from a distal end of the shank along an axilength of the shank, the shank comprising a superelast nickel titanium alloy, and
  - (b) after step (a), heat-treating the entire shank at a terperature above 25° C. up to but not equal to the melting point of the superelastic nickel titanium alloy,
  - wherein the heat treated shank has an angle greater than degrees of permanent deformation after torque at degrees of flexion when tested in accordance with IS Standard 3630-1.



## US Pat. No. 8,876,991 – Claims 13 ar

- 13. The method of claim 12 wherein:
- step (b) is performed in an atmosphere that is unread ambient or any other acceptable heat treatment pro-
- 14. The method of claim 12 wherein:
- the temperature is from 300° C. up to but not equal to melting point of the superelastic nickel titanium all



# DOCKET

# Explore Litigation Insights



Docket Alarm provides insights to develop a more informed litigation strategy and the peace of mind of knowing you're on top of things.

### **Real-Time Litigation Alerts**



Keep your litigation team up-to-date with **real-time** alerts and advanced team management tools built for the enterprise, all while greatly reducing PACER spend.

Our comprehensive service means we can handle Federal, State, and Administrative courts across the country.

### **Advanced Docket Research**



With over 230 million records, Docket Alarm's cloud-native docket research platform finds what other services can't. Coverage includes Federal, State, plus PTAB, TTAB, ITC and NLRB decisions, all in one place.

Identify arguments that have been successful in the past with full text, pinpoint searching. Link to case law cited within any court document via Fastcase.

### **Analytics At Your Fingertips**



Learn what happened the last time a particular judge, opposing counsel or company faced cases similar to yours.

Advanced out-of-the-box PTAB and TTAB analytics are always at your fingertips.

#### API

Docket Alarm offers a powerful API (application programming interface) to developers that want to integrate case filings into their apps.

#### **LAW FIRMS**

Build custom dashboards for your attorneys and clients with live data direct from the court.

Automate many repetitive legal tasks like conflict checks, document management, and marketing.

#### **FINANCIAL INSTITUTIONS**

Litigation and bankruptcy checks for companies and debtors.

### **E-DISCOVERY AND LEGAL VENDORS**

Sync your system to PACER to automate legal marketing.

